


## Software User Guide



Including beta (complete, reviewed) chapters.

Version 18, Revision: A-2436

Dated 09/04/2021

**Asylum Research**  
an Oxford Instruments company

# Contents

<b>I Software Overview</b>	<b>2</b>
1 Software Environment	4
2 Typical Use	9
3 All Things Video	14
4 Navigating your Data	17
5 Software Customization	22
6 Miscellaneous Procedures	23
<b>II Data and Image Analysis</b>	<b>28</b>
7 Introductory Image Analysis	31
8 Force Curve Analysis	92
9 Particle Analysis	130
10 Video Overlay	131
11 ARgyle 3D Data Visualization	132
12 Managing Image Data	152
13 Exporting Raw Data	173
<b>III Programming</b>	<b>176</b>
14 Basic Programming Tutorial	178
15 Using the Command Line	193
16 Help! Bug Reporting	196
<b>IV Bibliography, Glossary, and Index</b>	<b>203</b>

# Introduction

## Volumes of the AR SPM Software User Guide

The Asylum Research Software manual comes in volumes. To date these volumes are:

**Part I Overview.** This part of the user guide covers a basic overview of how the software is organized and what the basic controls are for.

**Part II Data and Image Analysis.** Once the data are collected, this part of the user guide describes many of the powerful analysis capabilities of the AR software. There are subsections on images and force plots.

**Part III Programming.** The AR software can be customized for your needs on many levels, from custom control panels to simple macros to completely open-source code for all the Igor Pro based code, which defines nearly all the AR software functionality.

**AR Software version** It is assumed that AR Software version 14 or later is installed on your system.

### Getting Help

For additional help with your Asylum Research instrument, including software support, refer to: <https://afm.oxinst.com/support/>

**Updates to the Manual** Bundled with the software, updates can be found at <http://support.asylumresearch.com>.

**Send Feedback** Send e-mail to [sba.manuals@oxinst.com](mailto:sba.manuals@oxinst.com) (← clickable link) and mention which version of the user guide you are using and what chapter and section you're commenting on.

## Prerequisites

We recommend that you have a running AFM, or at least a functioning copy of the AR software installed on your computer. For an overview of a properly set up MFP-3D AFM, please refer to *MFP-3D User Guide, Chapter: Installation*. Likewise, for the Cypher AFM a properly operating AFM system includes a PC with the AR software installed. In case this software requires an upgrade, or you want to install a copy on another computer (handy for image analysis), please see <https://support.asylumresearch.com/forum/content.php?150-Copies-of-Igor-instructions>.







## Part I

# Software Overview

**Part I: Who is it for?** Every user should read this part at least once to get a general sense of the structure and use of the SPM software.

# Part Contents

<b>1</b>	<b>Software Environment</b>	<b>4</b>
1.1	Basic software components	4
<b>2</b>	<b>Typical Use</b>	<b>9</b>
2.1	Launching the AR SPM Software	9
2.2	Collecting and Saving Data	10
2.3	Keeping Things Organized	11
2.4	Shortcuts and Tricks	12
<b>3</b>	<b>All Things Video</b>	<b>14</b>
3.1	AR Video Panel	14
3.2	Video Calibration Panel	15
3.3	Guesstimating Scan Area	16
<b>4</b>	<b>Navigating your Data</b>	<b>17</b>
4.1	Igor Graph Primer	17
4.2	Viewing, Analyzing, and Resaving Stored Data	18
4.3	Printing, Exporting, and Presenting Data	19
<b>5</b>	<b>Software Customization</b>	<b>22</b>
5.1	Creating User Panels	22
<b>6</b>	<b>Miscellaneous Procedures</b>	<b>23</b>
6.1	Lifting the Tip Off the Surface	23
6.2	Engaging on Reflective Surfaces	24
6.3	Adjusting Mis-Aligned Mirror to View Cantilever	24
6.4	Tip NudgerXYZ	24
6.5	Smart Start™	25
6.6	Determining the MFP-3D Software Version of a Saved Expt	25
6.7	Deactivating the Closed Loop Scanner	25
6.8	Sum and Deflection Meter Panel Freeze	26
6.9	Rescan Bus	26
6.10	User Callbacks	26

# 1. The AR SPM Software Environment

CHAPTER REV. 2428, DATED 08/23/2021, 18:52.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

## Chapter Contents

1.1	Basic software components . . . . .	4
1.1.1	Igor Pro . . . . .	4
1.1.2	The Asylum Research XOP . . . . .	5
1.1.3	The AR SPM Template . . . . .	6

Before we get into describing the specific functionality of the AR SPM software, it helps to understand a little about how it is built and how it differs from other instrument software you may have used in the past. The AR SPM software was not written completely from scratch but builds upon a program called Igor Pro, written by Wavemetrics. This quote from Wikipedia nicely sums up its history:

“IGOR Pro is a scientific data analysis software, numerical computing environment and programming language that runs on Windows or Mac operating systems. It is developed by WaveMetrics Inc., and was originally aimed at time series analysis, but has since then evolved and covers other applications such as curve fitting and image processing. It comes with a comprehensive programming language with compiler, but many functions are also accessible through menus. IGOR Pro’s strengths include its graphics capabilities, and the possibility of extending the built-in functions with external operations XOP allowing data acquisition, manipulation and analysis features, communication with external devices and in principle any other task that can be programmed in C or C++.”


In a nutshell, Asylum Research writes C++ code which communicates with and controls the SPM controller on a very basic level. This code cannot be altered by the instrument user. Asylum Research also writes code in the Igor Pro programming language which executes commands in this Xop to send parameters to the controller and retrieve data from the controller. More Igor code creates a complete user interface to give you, the user, control over the instrument. Display of the data, curve fitting, filtering, much image analysis, and many more features take advantage of a huge number of built in Igor Pro functions.

Let’s open up the software and identify some of these basic components.

## 1.1. Basic software components

### 1.1.1. Igor Pro

You start up Igor Pro from the Start Menu on your PC computer: *Start > All Programs > Igor Pro > Igor*. You will see something like [Figure 1.1 on page 5](#).



**Figure 1.1.:** Basic Igor Pro starting screen.

Notice the menu bar has items *File*, *Edit*, *Data*, *Analysis*, *Macros*, *Windows*, *Table*, *Misc*, and *Help*. These menu items belong to Igor Pro and give you access to much of the Igor Pro data graphing and analysis functions. Note that the *Table* item will change depending on what type of window you click on (graph, panel, notebook, table, etc.)

If you are planning to do serious analysis on your AFM data, this is a good time to take a tutorial on Igor itself. If you are only after a few AFM images to look at, this step can be skipped. However, if you plan to use the AFM much, the tutorial is a very worthwhile investment of your time.

**Just do it!**


The Igor Pro *Volume I - Getting Started* manual (especially the Igor Guided Tour portion) can be found from the Igor menu bar: *Help > Manual*. Then click the 'Open Online Manual' button, and a pdf file of all of the Igor manual volumes opens. It takes two to three hours to complete and is an excellent way to learn about the basic graphing and analysis functionality of Igor Pro. Although it is not necessary to complete the Igor Pro portion of the "Getting Started" manual at this time, it is a highly recommended part of all new user training.

### 1.1.2. The Asylum Research XOP

You don't really need to know much about the XOP. It is an Igor Pro term which stands for eXternal Operations. They are typically computer files based on raw C++ code which allow Igor to cooperate with external instruments and software. Asylum Research writes this software, and you can be largely oblivious.

You will interact with the Asylum Research XOPs on some level. Notice the strip of icons at the bottom of [Figure 1.1 on page 5](#). You will not find that in the Igor Pro manual. These are some small controls which are a result of an Asylum Research XOP and tell you about the status of the AFM controller. We'll refer to this strip of icons as the "software tray" or "system tray".

While we don't want to get bogged down in the detailed purposes of all of these controls, it is good to have a basic grasp.



① **Igor Pro status** If Igor is ready to accept a software request from you, it will say Ready. If not, it's thinking or calculating indicative of the Abort Button with a rotating quartered circle. This one is an Igor Feature and not part of the AR XOPs.

② **Rescan Smartstart Bus** Click this when the

SPM software gets confused. This may be necessary if Igor freezes up or when adding new components to the system (e.g., heaters, different cantilever holders that have identity resistors in them) that are not recognized.

③ **Smart Start Status** Click this to see what MFP-3D components are communicating with the Controller and PC computer. Furthermore, individual information (e.g., temperature at, serial number, etc.) on each component can be accessed by clicking the triangular button to the right of each component icon.

④ **MFP-3D status** This displays if the AFM instrument is ready or if the controller or other hardware isn't found.

⑤ **Quicklink to AR online support** AR engineers can control your MFP-3D over the web during online help sessions.

⑥ **CCD camera** Click this to get camera screen to come up.

⑦ **MFP-3D XOP version** May be asked of you during a support call.

### 1.1.3. The AR SPM Template

**Note** This may still be called the MFP-3D template for some time to come. For many years the MFP-3D was sole instrument which communicated with the AR SPM software. Since it now also drives the Cypher AFM, the software name was changed. However, references to the MFP-3D AFM will continue to remain in the software. Don't be alarmed, as your software will work with any Asylum Research AFM.

The third arm of the SPM software is written in a special programming language defined by Igor Pro, and it looks a bit like C. Nearly all the user interface controls are written in this language.

To load up this Igor Pro computer code, from the Igor menu bar select *File* ▸ *New MFP3D Template*. As explained earlier, this will probably change to New AR SPM Template, or some other name not explicitly associated with the MFP-3D. The Igor Window will now look more like [Figure 1.2 on page 7](#).


Here are some of the basic features, clockwise from the upper left. (Note that if you are viewing this file on a computer, you can zoom into the screenshot for closer scrutiny.)

**Master Panel** Upper left window. It has five tabs with controls and data displays for:

**Main** AFM imaging controls and control parameters

**Thermal** Cantilever thermal spectroscopy

**Force** Cantilever force vs distance curves



**Figure 1.2.:** Typical start up screen for the Asylum Research SPM Software, after the Mode Master panel has been closed. A few image panels have been left off the right of the screen, which usually extends across the second monitor of your system. Note that not all the panels shown here are in their original locations. It helps to manually tidy things up a bit so it looks like this figure.

**Tune** Cantilever vibrational tuning

**Fmap** Maps of force vs distance curves

**Master Channel Panel** During imaging, multiple data streams, such as cantilever amplitude, height, and phase, return from the SPM to the computer. This panel contains information about those streams and allows for some real-time scaling and processing.

**History Window** Also known as the Igor Command line window (Which was the lower window in Figure 1.1 on page 5). It will always have the name of your current experiment, which you can also see at the top of the border for the entire software window. On occasion items executed by clicking software buttons will generate some output here. Power users will type commands at the command line to accomplish a variety of advanced tasks.

**Sum and Deflection Meter** A real time display of various data such as cantilever deflection, amplitude, piezo voltage, and various other user definable channels. Also contains buttons for engaging and withdrawing the AFM tip.

**Image Windows** For each active channel on the Master Channel Panel, one image will appear on the screen. Before the scan starts, they are empty and appear small. They balloon to proper size as soon as scanning starts. The windows display real time, line by line, the sample topography, phase, amplitude, voltage, or any other measured quantity, as the cantilever raster scans the sample. There is usually one such window per active tab in the Master Channel Panel (Lower left window). While these windows are primarily a data displays, right clicking with the mouse can activate various commands such as zoom and translate. A white area at the bottom of this window shows you a real time “oscilloscope view” of the most recent line of image data, which are very useful when tuning imaging feedback parameters.


A few other things of note are:

**System Tray** Along the bottom of the screen. Icon controls relate to the status of connected instrument components. Low level software version is also displayed. This was discussed on 6.

**Menu Bar** Along the top of the screen. Note that there are a number of new items in addition to the ones that “belong” to Igor (See Section 1.1.1 on page 4). The new items, specific to the SPM functions, are *MFP Controls*, *MFP IP*, *Programming*, *Data Folders*, and *User Settings*.

**Q** Oops! I accidentally closed one of the control panel windows. How do I get it back?

**A** You can re-activate the panels via *AFM Controls* (*MFP Controls* in earlier versions of the software) in the menu bar.

When you open the  software you are opening an Igor Pro “Experiment” in which extra software specific to the operation of the Cypher has been loaded. Beside graphing data and curve fitting, Igor Pro provides a programming environment with full graphical user interface capabilities including control panels with buttons, sliders, and fields for data input and output. These panels can be mixed freely with graphs and tables of data. Hence the controls for the Asylum SPMs are a mixture of graphs and control panels.

**Tip**

Nearly each individual item in the software control panels has a small question mark button next to it. Click the button to read the relevant parts of the software help file.

## 2. Typical Uses of the AR SPM Software

CHAPTER REV. 2428, DATED 08/23/2021, 18:52.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

### Chapter Contents

2.1	Launching the AR SPM Software . . . . .	9
2.2	Collecting and Saving Data . . . . .	10
2.2.1	Collecting Data . . . . .	10
2.2.2	Saving Data . . . . .	11
2.3	Keeping Things Organized . . . . .	11
2.4	Shortcuts and Tricks . . . . .	12
2.4.1	Igor Hot keys . . . . .	12

The Asylum Research Scanning Probe Microscope Software (AR SPM Software) is capable of many, many things—more than can reasonably be written in one single manual. You will probably only use a small fraction of its capabilities, and this chapter can help you find that portion of the software that is likely to be of use of you.


### 2.1. Launching the AR SPM Software

1. Boot up your PC (likely the one that came with installed AFM System).
2. On a PC computer select from the windows star menu: *Start > All Programs > Igor Pro > Asylum Research > MFP3D*. (Please see the note on page 6.) In case you have copies of the software installed on your computer, this will launch the original copy installed on the computer

–OR–

Look for the same shortcut in the Windows Task Bar (usually next to the Start Menu). Again, this launches the original installation.

–OR–


Look for the  logo on your desktop. The accompanying number with the icon denotes the major software release version. If there are multiple copies of the software on your computer, there should be one shortcut icon for each copy. Double-click the icon of choice to launch the software.

With the software launched, you'll see something like [Figure 2.1 on page 10](#).

**Note**

The windows in the ARM SPM Software do not have permanently assigned positions. Get in the habit of arranging things neatly, as shown in [Figure 2.1 on page 10](#), so you have an intuitive idea of where to look for certain controls.





**Figure 2.1.:** Typical startup screen after closing the Mode Master window, a few other image windows, and a little re-arranging

A VERY brief guide to what's on the screen, from top to bottom:

**Top window border** This displays your version of Igor and AR SPM software.

**Menu bar** The menu bar links to all of the panels in the AFM software. It rarely directly controls actions of the SPM hardware but rather allows you to launch control panels that appear in the main part of the screen.

**Main Screen** The main screen is filled with control panels for controlling the SPM hardware and software functions (here arranged around the perimeter of the screen) and data display windows (one shown at the center of the screen) for viewing and interacting with data. Note the white window near the bottom.

**System Tray** Also known as software tray. Shows you the status of Igor Pro and of all the connected SPM equipment. It also shows the version of low-level software known as the “XOP”.

Now that the software is up and running (very lower left bit of the screen should say “ready”), you will probably want to either collect new data (see [Section 2.2 on page 10](#)) or view and analyze existing data.

## 2.2. Collecting and Saving Data

### 2.2.1. Collecting Data

1. Run through a tutorial. No amount of discussion will illustrate the basics of the software and how it controls the hardware.

For the Cypher AFM see<sup>1</sup>

For the MFP 3D AFM see<sup>2</sup>

2. Since a basic command of AC mode imaging is necessary before more exotic imaging modes are attempted, the general business of data collection is covered in the AC Mode Imaging chapter.
3. For non-imaging data collection, consider the following chapters:
  - Force Curves:** <sup>3</sup>
  - Lithography:** <sup>4</sup>
  - Spring Constant Calibration:** <sup>5</sup>

In general, data collection is managed through panels with buttons and fields containing numerical parameters which control the data-taking process. Each imaging mode usually has its own panel. Some of those panels, such as the Master Panel, which controls standard tasks such as AC and Contact Mode Imaging, appear automatically when the software boots up. Panels for more advanced imaging modes can be opened via items in the AFM Imaging Menu in the main menu bar at the top of the software window.

### 2.2.2. Saving Data

Images are saved to disk automatically at the end of every image IF you select the appropriate checkbox, such as the one called “Save Images” near the lower left corner of the Master Panel in [Figure 2.1 on page 10](#). DO THIS or you will lose data!

Force curves (see the Save tab under the Force Curve tab on the Master Panel) take up less memory and can be automatically saved to disk or saved in memory.

If, at the end of a data collection session, you are not certain that all data in memory are saved to disk, from the main menu bar select *File > Save Experiment As*. Igor Pro will save every last bit of information currently in memory as a .pxp file (packed experiment file). When you open such a file again, your entire data collection session is brought back up, down to the positioning of every window.

## 2.3. Keeping Things Organized

As mentioned before, the control panels and data display panels in the AR SPM Software do not have fixed positions. Here are a few tips on how to keep things organized and how to find things that you may have “lost”:

- Get in the habit of arranging windows and frequently used control panels in a consistent fashion (such as [Figure 2.1 on page 10](#)).

<sup>1</sup> *Cypher User Guide, Chapter: Tutorial: AC Mode in Air, Std. Scanner..*

<sup>2</sup> *MFP-3D User Guide, Chapter: Tutorial: AC Mode Imaging in Air..*

<sup>3</sup> *Applications Guide, Chapter: Force Curves Acquisition..*

<sup>4</sup> ?..

<sup>5</sup> *Applications Guide, Chapter: Spring Constant Determination..*

- Close windows and panels as soon as you don't need them anymore. A good example is the graph which appears after tuning a cantilever in AC mode.
- Use the Window Manager Panel (Ctrl+1 or select *AFM Analysis Menu > Window Manager*) to see organized lists of all open windows and panels. Click on what you lost to bring it to the front. Similar controls can be found in the Windows Menu.
- It can happen that a window or panel in the AR SPM software is not to be found because it has landed "off the screen". From the main menu bar, select *Windows > Control > Retrieve All Windows* to bring all the windows and panels back into view.
- There are various tools for simplifying the SPM control panels and highlighting items on them, including building your own custom control panels by picking and choosing controls from various existing panels.

## 2.4. Shortcuts and Tricks

### 2.4.1. Igor Hot keys

**ctrl+I** Show Info, Igor A&B cursors

**ctrl+A** Auto scale graph

**ctrl+F** Find

**ctrl+T** Tool mode

**ctrl+W** Closes forwardmost window

**ctrl+D** Duplicates forwardmost window

**ctrl+C** Copy

**ctrl+V** Paste

**ctrl+Y** Window control, can rename windows

**ctrl+O** Open experiment

**ctrl+P** Print

**ctrl+L** Load waves

**ctrl+J** Bring History window forwardmost

**ctrl+1** Toggle tool tips on or off

**ctrl+shift+1** Window manager

**ctrl+2** Do thermal

**ctrl+3** Single force pull

**ctrl+4** Cantilever auto tune

**ctrl+6** Open S&D meter

**ctrl+7** Open Master Channel Panel

**ctrl+8** Open Engage Panel (Cypher or Jupiter only)

**ctrl+9** Open NAP panel

**ctrl+F1** Igor search panel

**ctrl+F2** Open Mode Master

**ctrl+M** Open Igor Procedure window

## 3. All Things Video

### Chapter Contents

3.1	AR Video Panel . . . . .	14
3.1.1	Capturing movies: . . . . .	14
3.2	Video Calibration Panel . . . . .	15
3.3	Guesstimating Scan Area . . . . .	16
3.3.1	Centering tip in center of scan area in CCD image . . . . .	16

### 3.1. AR Video Panel

To capture images (or acquire movies) from the CCD camera, the AR video panel is needed.

1. Go to *MFP-Controls > Video Panel*; the AR Video Panel appears.
2. Click Top view or Bottom view if you have a Dual-View Stand-Alone base.
3. Have the 'Scale Image' checkbox checked if you have calibrated the image ; if the image hasn't been calibrated, the scale on the axis isn't accurate.
4. Confirm the image is in the desired focus.
5. To save images, click the 'Save to Disk' checkbox or 'Save to Memory' checkbox. It is best to save images as .tiff or .ibw file.
6. Click the 'Capture' button; the new image should be displayed immediately.
7. To get an assessment of the scan area relative to the cantilever position, click the 'Show Scan Box' checkbox.
8. On the displayed image, place cursor towards the end of the cantilever, right-click, and select *Set Tip Location*; this displays a red circle indicating the tip location (just as in the 'Point and Click' functions in 'Go There' subtab).
9. Activate the 'Show Scan Box' checkbox. Two squares will appear on the displayed image; the outer square shows the entire XY range (90µm), while the inner square shows the Scan Area relative to the XY offsets.
10. As the Scan Area and XY offsets are changed, the display window is immediately updated. .

#### 3.1.1. Capturing movies:

Movies of CCD images can also be captured with the AR Video panel as follows:

1. Select a capture rate of no more than a couple Hz. Any more than that, Igor will have trouble with collecting a lot of frames per second.
2. Choose the file type to save as. It's recommended to choose .ibw or .tiff as the file type.
3. Give it a filename.
4. Click the 'Start Background Capture' button. This will start collecting at the specified rate.
5. When finished, click the 'Stop Background Capture' button.
6. Once complete, click the 'Make Movie' button to designate a frame rate for the movie. The movie will be strung together into a video file.

### 3.2. Video Calibration Panel

The CCD image can be calibrated such that the scale bar along the X&Y axis of captured / displayed CCD images is accurate. This process is very easy to do, and it's described in full detailed in a User Panel. It's best to use a diving board cantilever for this.

1. Go to *AFM Control > User Panels and Funcs > Video Calibration Panel*.
2. This brings up the AR video calibration panel along with a Help menu to walk you through it. This protocol is very easy to follow along with.
3. Confirm that CCD camera is on and cantilever in focus; the tip doesn't have to be engaged with sample.
4. Confirm that the probe is oriented in its pocket such that the cantilever length is parallel to the X axis of the CCD window.
5. Enter nominal dimensions of cantilever for 'Lever Length' and 'Lever Width'. The default cantilever in the panel is for a Olympus AC 240 (i.e., 240  $\mu\text{m}$  long, about 30  $\mu\text{m}$  wide)
6. Turn off 'Scale Image' checkbox.
7. Click 'Capture' button
8. Click 'Display' button
9. On the new display image, call up the Igor cursors (ctrl + i)
10. Place the cursor at the end of the cantilever; place the other cursor at the base.
11. Click the 'Calibrate Length (X)' button.
12. Place the two cursors at each side of the cantilever to measure its width.
13. Click the 'Calibrate Length (Y)' button
14. Click the 'Video Params' button

The Pixel Size should now be updated. It's not a bad idea to jot that down somewhere in case you'll be using multiple CCDs (i.e., top and bottom views) during the experiment.

Now when the 'Capture' button is clicked, the displayed images should be set to proper scale.

If more precision is needed for the calibration, reticle calibration grids can be used to calibrate the X and Y lengths.

### 3.3. Guesstimating Scan Area

Here's how to roughly determine if the selected scan area will contain a feature of interest as seen in the Top View CCD camera. This method is a quick and easy way to figure out if the scan area will reach the sample area of interest in the CCD image without calibrating the Video Panel (see [Section 14.8 on page 187](#)).

1. Increase the Set Point voltage (if in AC mode) past the 'Free Air' amplitude voltage so the tip comes off the surface ; if in Contact mode, decrease the Set Point voltage past the free air deflection voltage so the tip comes off the surface .
2. Click the 'Frame Up' button and watch where the tip scans over in the CCD camera.
3. Click the 'Frame Down' button and watch where this scans over.

Make adjustments to coarse XY position micrometers or XY scan area if necessary

#### 3.3.1. Centering tip in center of scan area in CCD image

Go to Force tab in the Master Panel. With no spots selected, have suffix on zero, and click the 'Go There' button. This will put the tip in the center of the scan area. The tip will now be in the center of its scan range.

## 4. Navigating your Data

CHAPTER REV. 1984, DATED 11/10/2017, 17:24.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

### Chapter Contents

4.1	Igor Graph Primer . . . . .	17
4.1.1	Installing a grid on image screen: . . . . .	17
4.2	Viewing, Analyzing, and Resaving Stored Data . . . . .	18
4.2.1	Image files . . . . .	18
4.2.1.1	Opening and viewing images . . . . .	18
4.2.1.2	Modifying and analyzing images . . . . .	18
4.2.1.3	Saving changes . . . . .	19
4.3	Printing, Exporting, and Presenting Data . . . . .	19
4.3.1	Quick and Dirty Printing . . . . .	19
4.3.2	Pretty Printing with Layouts . . . . .	19
4.3.3	Exporting Graphs or Layouts . . . . .	20
4.3.3.1	Save as a bitmap to disk . . . . .	20
4.3.3.2	Copy as a bitmap to the clipboard . . . . .	20
4.3.4	Exporting Raw Data . . . . .	20
4.3.5	Igor Electronic Notebooks . . . . .	21

### 4.1. Igor Graph Primer

#### 4.1.1. Installing a grid on image screen:

If it would help you to have a grid on the image window (or any Igor plot) in real time or offline; this is easy:


1. Double-click on one of the axis of the plot/graph. This shows the Igor Modify axis.
2. Click on the Ticks and Grids tab.
3. From the Grid menu, select *On*.
4. Select color and line thickness or any other aesthetic from this tab. For more on the Modify Axis panel, see the Getting Started manual of the Igor Pro Manual.
5. Click the 'Do It' button.



## 4.2. Viewing, Analyzing, and Resaving Stored Data

### 4.2.1. Image files

The entire topic of opening, viewing, analyzing, and saving files is outlined in great detail in [Chapter 7 on page 31](#). Please read at least the first few sections and work through the examples to get familiar with this process. What follows here is a very brief synopsis.



**Figure 4.1.:** Image Analysis Display

#### 4.2.1.1. Opening and viewing images

1. From the main menu bar select *AFM Analysis > Browse Saved Data*.
2. Choose a directory where the images are stored and show a panel with thumbnails of all the images in the selected directory.
3. Double-click on any image to show an image display window ([Figure 4.1a on page 18](#)).
4. The tabs on the image show the various data channels/layers, such as height, amplitude, phase, etc.
5. To view data in colorful 3D, navigate to *AFM Analysis > 3D surface plots*. See [Chapter 11 on page 132](#) for details.

#### 4.2.1.2. Modifying and analyzing images

1. Click the 'M' button to bring up controls to filter or otherwise modify the image channels/layers.

2. Click the 'A' button to bring up controls to analyze attributes (e.g., roughness, histograms, etc.) of the image. For additional details, see [Chapter 7 on page 31](#).

#### 4.2.1.3. Saving changes

1. Choose *Save Image* from the pulldown menu shown in [Figure 4.1a on page 18](#). Note that any modifications you make to your data channels are stored in additional new image layers. The original data channels (layers) will not be altered.
2. If you are concerned about altering original data, first make a copy of the file from within Windows and then work on that copy.

**Note** While there are other ways to save data in Igor, the method above is the only reasonable way to save SPM Image files. Don't be tempted by alternatives from the File and Waves menus.

### 4.3. Printing, Exporting, and Presenting Data

#### 4.3.1. Quick and Dirty Printing

1. Select an image window ([Figure 4.1a on page 18](#)) and tab so that the layer of interest is showing.
2. From the main menu bar, select *File > Print Graph*.
3. Choose the 'Same Aspect' button in the following print dialog box to maximize the image size on the paper.
4. Select your printer and click 'OK'.

**Note** This process also works for any other graphs on the screen, such as cantilever tunes and force curves.

**Note** This process will not work for 3D images displayed via ARGyle. The next section describes how to put the image in a layout first.

#### 4.3.2. Pretty Printing with Layouts

Igor Layouts allow you to arrange images, graphs, text, and other graphical elements on a single page for printing or bitmap export. For more information, refer to the Igor Pro manuals from Wavemetrics. For a barebones rundown, read on:

1. Select an image window ([Figure 4.1a on page 18](#)) and tab so the layer of interest is showing.
2. Choose *Export to Layout* from the pulldown menu shown in [Figure 4.1b on page 18](#).
3. A new layout window will appear with a copy of the image and some text with image parameters that are fine-tuneable by choosing *Set Layout Note Fields* from the Commands pulldown menu.
4. Rearrange or scale the items as needed or add other graphical elements using the toolbar.
5. If need be, repeat the first few steps and add more images to the layout.

6. To add graphs, such as force curves or cross sections, to the same layout, click the 'Layout' button above the graph.
7. To add 3D views, from the Master ArGL panel, select New tab, and click the 'Export to Layout' button.
8. From the main menu bar choose *File > Print Layout* to print the layout.

### 4.3.3. Exporting Graphs or Layouts

#### 4.3.3.1. Save as a bitmap to disk

1. Select the graph, image, 3D view, or layout of your data.
2. From the main menu bar, choose *File > Save Graphics* and use the dialog box to export a bitmap in various resolutions and file formats.

#### 4.3.3.2. Copy as a bitmap to the clipboard

1. Select the graph, image, 3D view, or layout of your data.
2. From the main menu bar, choose *Edit > Export Graphics* and use the dialog box to copy a bitmap in various resolutions and file formats to the Windows clipboard.
3. Choose any other application, such as Paint, Photoshop, Word, e-mail, etc., and choose *Edit > Paste*.

### 4.3.4. Exporting Raw Data

Quick and dirty recipe for exporting line graph data:

1. Right-click on the trace you want to export.
2. Select *Edit* (which is followed by the name of the trace).
3. Right-click on the Table which appears.
4. Select *Digits*, then select the highest number there (16 at the moment).
5. Press Ctrl + A, to select the entire table
6. Press Ctrl + C to copy the contents of the table to the Windows Clipboard.
7. Paste your ascii data into any editor program (suggestion: *Windows Start Button > All Programs > Accessories > Notepad*).
8. Save the data as a text file (.txt).
9. Import the data into the analysis program of your choice.

For a more formal process, or an automated process for many graphs or images, please see on page 173.

### 4.3.5. Igor Electronic Notebooks

The Igor Electronics notebook editor is something of a word processor. You can copy any graph, table, or layout and paste it directly into the notebook editor. Then you can add text between the images and layouts, keeping a rough laboratory notebook of events. The notebook can be saved as html or rich text (.rtf) and edited more from within your favorite word processor.

1. From the main menu bar, choose *Windows > New > Notebook*.
2. In the dialog box, choose 'formatted text'.
3. Click the 'New' button.
4. Select any Igor graphic display window, layout window, or table window and press Ctrl+C (or from the main menu bar *Edit > Copy*).
5. Select the notebook window and paste (Ctrl + v or *Edit > Paste*)
6. Type text between the images if desired.
7. To save, go to *File > Save Notebook* or *Save Notebook As*. Choose .rtf file format for further editing in a word processor.

For more information on Electronic Notebooks see the Igor Manuals.

**Tip**

With a notebook selected (in the foreground) choose from the main menu bar *Notebook > Special > Enable Updating*. In the next dialog box, check the box and click the 'OK' button. Now, if you update any graph that is also in the notebook, you can right-click on the time in the notebook and select *Update Selection*. The graph will refresh to look like the original. With nothing selected in the notebook, you can also select from the main menu bar: *Notebook > Social > Update All Now* to refresh every copied item.

## 5. Customizing the AR SPM Software User Interface

### Chapter Contents

5.1	Creating User Panels . . . . .	22
-----	--------------------------------	----

### 5.1. Creating User Panels

In the main menu of the AFM software, select *Programming > User Panels*. This opens the ARUP Manager Panel (Asylum Research User Manager). The software engineers have included a number of panels, some of which can be used independently and others which must be opened within a Mode Master template.

- In the New Panel Name window at the top, enter the name for the new window.
- Click the 'Start' button. A new panel will open.

At this point, parameters, buttons, and radio buttons can be copied and pasted from other AFM software panels.

It is also possible to include graphs and other objects.

## 6. Miscellaneous Procedures

### Chapter Contents

6.1	Lifting the Tip Off the Surface . . . . .	23
6.2	Engaging on Reflective Surfaces . . . . .	24
6.3	Adjusting Mis-Aligned Mirror to View Cantilever . . . . .	24
6.4	Tip NudgerXYZ . . . . .	24
6.5	Smart Start™ . . . . .	25
6.6	Determining the MFP-3D Software Version of a Saved Expt . . . . .	25
6.7	Deactivating the Closed Loop Scanner . . . . .	25
6.8	Sum and Deflection Meter Panel Freeze . . . . .	26
6.9	Rescan Bus . . . . .	26
6.10	User Callbacks . . . . .	26

This section includes procedures the author finds valuable to perform while using the MFP-3D. Some are commonly used, some are transparent, and others are rarely used; but, hopefully, this collection may be helpful.

### 6.1. Lifting the Tip Off the Surface

Anytime you need to adjust something where it would be a good idea to lift the tip off the surface, but you don't want to disengage (i.e., withdraw the tip), you can adjust the Set Point voltage to the false-engaged side of the free air voltage to lift the tip off the surface. This can be done with the Hamster wheel or manually typing in the setvar parameter window. Specifically, you might want to do this when:

- Shutting the door on the acoustic enclosure
- Changing image resolution
- Adjusting PD deflection: Be careful! This will likely move tip position up to a couple  $\mu\text{m}$ 's just by touching the head
- Adjusting thumbwheel position

If in AC mode, increase the Set Point voltage to above the Free Air amplitude. Note that the Z Piezo voltage will decrease.

If in Contact mode, lower the Set Point voltage to a value more negative than the Free Air deflection.

## 6.2. Engaging on Reflective Surfaces

Aside from monitoring the Free Air amplitude in the Sum and Deflection Meter while engaging, the CCD camera can also be used when engaging on reflective surfaces to monitor relative tip substrate separation based on the cantilever shadow distance from actual cantilever. Focus CCD onto cantilever, then focus just slightly beyond that. This will allow surface to come into focus when the tip nears it. When the tip is engaged, the cantilever reflection should not be seen.

Another tip that Scott MacLaren taught the author: If using a Stand Alone base, you can change the aperture to slightly closed to get a bit of shadow around the edges of the video window. As the tip approaches the surface, the image of the surface will come into sharp focus. This process allows the tip to be brought towards the surface quickly with the thumbwheel for a coarse approach, then slowed for engagement.

## 6.3. Adjusting Mis-Aligned Mirror to View Cantilever

Sometimes the previous user can mis-align the X&Y knobs on the top view optics mirror so much that the cantilever can't readily be found in the top view CCD- this problem becomes compounded if the focus is also screwed up.

Instead of blindly manually rastering the mirror looking for the tip, here's a quick way to fix this problem: make sure a probe is in the tip holder to give a correct target for you.

1. Look straight down into the mirror housing (i.e., plumb), just as the light path does (see Figure 14.10A)
2. Look for crisp clear trapezoid shape of quartz window. To assist this, wave your finger in front of the window and look for the fleshy color of your finger in the mirror (Figure 14.10B); also notice the opaque area of cantilever holder. The author looks for the crisp line between the opaque area and the clear quartz window.
3. Move the XY mirror positioning knobs to bring cantilever into middle of trapezoidal window.
4. Once the mirror is positioned correctly, the cantilever should be easily viewed (Figure 14.10C)
5. Place the head back over the sample, and focus to suit needs.

## 6.4. Tip NudgerXYZ

This is very useful if the tip must be moved to a surface feature, in closed or open loop, that is resolvable in the CCD camera. It's also very intuitive to use, as follows:

1. In the main menu, go to *Programming > Start User Panels*.
2. Double-click on 'NudgerXYZ'.
3. Enter the 'Step Size' that you want to step the tip with.
4. Activate the X, Y, or Z axis checkboxes if you want the nudge to occur under *Closed Loop* control or *Open Loop*.
5. Click the desired axis direction. Depending on the step size, the tip's movement can be seen visually, and there is typically a concomitant scanner noise.

## 6.5. Smart Start™

Sometimes when switching between software versions (e.g., old vs. newer), or if there has been a repair performed on some hardware, there are calibration files that change and are highlighted in red in the far left. Smart Start™ will catch this and display the ParmPanel window asking you to change the calibration in the software so that the measurements are precise.

This panel can seem a bit of a sensory overload initially, but it's relatively easy to deal with, as follows:

1. Look at the Use This column and read the value for any row that is highlighted in red.
2. In that same row, compare the Software and Hardware column values and 'move->' the value that is different from the Use This column.
3. Click the 'Use the 'Use This' Values' button at the bottom of the panel. This will activate the button below it.
4. Click the 'Save the 'Use This' Values to Software' button.
5. The process is now complete, and the window can be removed.

## 6.6. Determining the MFP-3D Software Version of a Saved Expt

Say you have a saved experiment and need to know what version of software it was originally created with. The impetus of this is from trying to reload an experiment offline in which errors are experienced because the version of software you are loading with is not the same as the version the data was acquired with. To resolve this:

1. Go to *Data > Data Browser*. The Igor data browser window appears.
2. Confirm that only 'Variables' and 'Info' checkboxes are selected.
3. Click the 'Browse Expts' button. A browse dialogue appears. Find the saved .pxp experiment needed. A second column appears on the right of the data browser.
4. Go into that column and drill down into root/MFP-3D/Main/variables/VerDate. Once selected, the version will appear at the bottom of the data browser.
5. Now you can match that experiment for opening the experiment offline. If you need to know the software version, contact support@asylumresearch.com .

Another easy way to determine the version is to view the Asylum Research tab in the Properties information using Windows Explorer or My Computer.

## 6.7. Deactivating the Closed Loop Scanner

To deactivate the Closed Loop Scanner:

1. Go to Setup in the main tab of Master Panel.
2. Towards the bottom, click on the 'Show?' check box next to 'Scan Mode' pull-down. It includes the option to turn off the Closed Loop feature.



## 6.8. Sum and Deflection Meter Panel Freeze

This can occur if the cable connection isn't very good. With time, and all that repeated cable twisting from putting the head on its back while attending to the sample, the plug can wiggle its way out enough to cause the meter to freeze.

To fix it:

1. Power the controller off for 10 to 20 seconds, then turn back on.
2. While you're waiting, wiggle the plug back in a bit more.

We have noticed that anytime one of the levels in the S&D meter aren't responding properly, it's a good idea to do some variation of the above procedure.


## 6.9. Rescan Bus

Sometimes if the system experiences a freeze, or the system is hanging because Igor has been confused, click the 'Rescan' button (lower left). This should "ping" all the hardware systems to make sure they are all communicating properly with each other.

## 6.10. User Callbacks

<https://support.asylumresearch.com/forum?t=97&highlight=survey>





## Part II

# Data and Image Analysis

**Part II: Who is it for?** Once the data are collected, this volume describes many of the powerful analysis capabilities of the AR software. It may also be useful with a second copy of the software on a computer not attached to the AFM itself.

# Part Contents

<b>7</b>	<b>Introductory Image Analysis</b>	<b>31</b>
7.1	Opening Stored Images	32
7.2	The Modify Panel	39
7.3	The Analyze Panel	69
7.4	Miscellaneous Operations	85
<b>8</b>	<b>Force Curve Analysis</b>	<b>92</b>
8.1	Opening Stored Force Plots	93
8.2	Display Tab	94
8.3	Pref Tab	99
8.4	Parm Tab	102
8.5	Cursor Tab	103
8.6	Modify Tab	105
8.7	WLC Tab	107
8.8	Analyze Tab	110
8.9	Spot Tab	114
8.10	Elastic Tab	115
8.11	Online Video	129
<b>9</b>	<b>Particle Analysis</b>	<b>130</b>
<b>10</b>	<b>Video Overlay</b>	<b>131</b>
<b>11</b>	<b>ARgyle 3D Data Visualization</b>	<b>132</b>
11.1	Viewing Saved Images in 3D	132
11.2	Overlaying Data on a Surface	133
11.3	Display Settings	140
11.4	Exporting Images	148
11.5	Activating Real Time ARgyle™	149
11.6	Advanced Command Line Control	150
11.7	Argyle Lite Windows App	150
<b>12</b>	<b>Managing Image Data</b>	<b>152</b>
12.1	Introduction	152
12.2	Data Channels	153
12.3	Selecting Image Pixel Density	153

12.4	Master Channel Panel . . . . .	154
12.5	Capturing Images . . . . .	158
12.6	Offset Scan Area . . . . .	160
12.7	Zooming Images . . . . .	161
12.8	Multiple Offset/Zoom Parameters . . . . .	161
12.9	Real Time Section Analysis: . . . . .	164
12.10	Other Features . . . . .	169
<b>13</b>	<b>Exporting Raw Data . . . . .</b>	<b>173</b>
13.1	Exporting Your Data in ASCII Format . . . . .	173

## 7. Introductory Image Analysis

CHAPTER REV. 2436, DATED 09/04/2021, 14:34.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.


### Chapter Contents


7.1	Opening Stored Images . . . . .	32
7.1.1	The Browse Graph . . . . .	34
7.1.2	The List Panel . . . . .	35
7.1.3	The Display Window . . . . .	36
7.1.3.1	Reprogramming Shortcut buttons . . . . .	38
7.1.4	AR Thumbnail Viewer . . . . .	39
7.2	The Modify Panel . . . . .	39
7.2.1	Flatten tab . . . . .	40
7.2.1.1	Flatten Example, features on flat substrate . . . . .	41
7.2.1.2	Flatten and Planefit Example, Lattice steps . . . . .	43
7.2.1.3	Long involved flatten explanation . . . . .	47
7.2.1.4	Making a manual mask . . . . .	48
7.2.2	Planefit tab . . . . .	52
7.2.2.1	Planefit Example . . . . .	53
7.2.2.2	Planefitting a step . . . . .	55
7.2.3	Erase tab . . . . .	57
7.2.4	Mask tab . . . . .	59
7.2.4.1	Mask calculation methods . . . . .	60
7.2.4.2	Hand drawing masks . . . . .	61
7.2.4.3	Copying masks . . . . .	62
7.2.5	Filter tab . . . . .	63
7.2.6	FFT tab . . . . .	64
7.2.7	History tab . . . . .	67
7.2.7.1	How to apply one process to another image: . . . . .	67
7.2.7.2	How to apply one process to many images: . . . . .	68
7.2.7.3	How to remove one step from a multi-step process: . . . . .	68
7.2.8	Prefs tab . . . . .	68
7.2.8.1	Modify Insert Layers . . . . .	68
7.2.8.2	Undo Only Last Planefit . . . . .	69
7.3	The Analyze Panel . . . . .	69
7.3.1	Roughness Tab . . . . .	70
7.3.2	Section Tab . . . . .	72
7.3.2.1	Straight Line Sections: . . . . .	72
7.3.2.2	Free Hand Lines . . . . .	74
7.3.3	Histogram tab . . . . .	76
7.3.3.1	Fitting . . . . .	77

	7.3.3.2	Data Modification . . . . .	77
	7.3.3.3	Other Options . . . . .	78
7.3.4	Particle Tab . . . . .		78
	7.3.4.1	Particle Analysis Example 1: Finding the Length of a Piece of DNA . . . . .	78
	7.3.4.2	Particle Analysis Example 2: Finding Average Height . . . . .	82
7.4	Miscellaneous Operations . . . . .		85
7.4.1	Saving modifications to images . . . . .		85
7.4.2	Custom work on image channels . . . . .		85
	7.4.2.1	Extracting layers: . . . . .	85
	7.4.2.2	Subtracting Image Layers: . . . . .	86
	7.4.2.3	Rotating Images: . . . . .	88
	7.4.2.4	Reinserting layers: . . . . .	88
7.4.3	Cropping Images: . . . . .		89
7.4.4	Igor layouts . . . . .		91

This section discusses how to perform basic image analysis on stored image files (file extension .ibw). Keep in mind that there are not necessarily hard and fast rules or sequences regarding image processing. It can be a trial-and-error process that depends heavily on the data. For this reason, the processing techniques are broken into sections below with some examples provided.

## 7.1. Opening Stored Images

1. **Prepare the software:**
  - Launch the AR SPM software and you will see the Mode Master window.
  - Click the 'Mode Master' button at the bottom of the screen any time you want to load the Mode Master: .
  - Click on the 'Standard...' tile.
  - Then click *Offline Image*.



2. **Open the image browser:**
  - From the menu bar select *AFM Analysis > Browse Saved Data*.

3.

**Select a directory of images:**

- If image data is not already loaded in the experiment, you are asked where you would like to load data from. If there is data already in the experiment, the Browse Graph is opened, showing those images. To open a new folder, you can click on the 'Change Directory' button in the upper-left corner of the Browse Graph (Figure 7.1 on page 34).
- From the AR Load Path dialog, you can click the 'Browse' button to select a new directory.
- OR-
- Use the pulldown menus to visit standard (Default Paths) or recently visited (Recent Paths) directories.
- Click 'That's It' to load the directory.
- The Browse graph will open (Figure 7.1 on page 34) as will the List Panel (Figure 7.2 on page 35) .




4.

**To open an image:**

- Double-click on the thumbnail to open the full-sized image.



## 7.1.1. The Browse Graph



**Figure 7.1.:** The Browse Graph displays the data as thumbnails. This window can be resized by clicking and dragging the corner.

The number labels in [Figure 7.1 on page 34](#) refer to the following:

1. **Image Data Type** Specific data channels (i.e., Height, Phase, etc.) can be selected in the *Image Data Type*. This is convenient if a specific image in a folder containing many images includes a data type that most of the other images do not have. For example, choosing NapPhase could help locate an MFM image among many non-NAP images.
2. **Which Image** Determines whether Trace or Retrace data is displayed. Selecting *Either* allows both types of data to be displayed.
3. **Vertical/Horizontal Number** Thumbnail display array sizes can be altered by changing the *Vertical* and *Horizontal* values. *Vertical* determines the number of rows in a window while *Horizontal* determines the number of columns.
4. **NextPage** This button advances the screen to the next page of thumbnails, much like using the Igor scroll bar on the right side of Browse Graph. The 'Page Up' and 'Page Down' commands are keyboard shortcuts for this function.
5. **Copy to Layout** The 'Layout' button takes a snapshot of the Browse Graph and dumps it into an MFP-3D layout window. See [Section 7.4.4 on page 91](#).
6. **Show Force Data** This is useful for showing the different data types within a directory of force-distance curves.

## 7.1.2. The List Panel




Figure 7.2.: The List Panel

Images can also be opened in the List Panel (Figure 7.2 on page 35). It consists of two columns: *Images on Disk*, which are the files stored in the folder, and *Images in Memory*, which displays the images that have already been opened in the experiment.

- To open an image file, double-click on the file name.
- To open multiple images, hold down the 'shift' or 'ctrl' key while selecting the images you want to open, then click the 'Open Files' button.
- Most of the function buttons are self-explanatory, so only a few are discussed, including:

**Refresh List** This updates the list of images in the currently selected directory. This is used if images are being acquired real time while stored images are subsequently being processed and analyzed.


**Change Directory** This is the same as selecting 'Change Directory' from the Browse Graph. See Section 7.1 on page 32. It allows you to change the data file directory which is displayed in the list.

**Save to Disk** Select one or more images to be saved over the original.

**Note:** Processed image layers are typically added to the original as new layers, and your raw data tabs are preserved.

## 7.1.3. The Display Window

When an image is loaded from the hard drive, it is displayed in the *Display Window*.



**Figure 7.3.:** The Main Image Display Panel

Figure 7.3 on page 36 shows an example of a Display Window image that will be used throughout this section. Some of the controls at the top of this window include:

1. **Range** This adjusts the color scale, thereby changing the contrast.
2. **Offset** This adjusts brightness by changing where the center of the color range is.
3. **Color Bar** This allows the color scale to be viewed or removed from the image.
4. **Auto button** This auto-scales the color scale and offset.
5. **ColorMap** This offers many color tables for displaying the image data.
6. **Stored data channel tabs** The measured data channels, such as amplitude, phase, height, and Zsensor, collected during a sample scan are saved as layers of the image. The tabs bring an image layer to the front for viewing. Tab Naming Conventions include:
  - **Raw Data** Raw Data tabs typically have three-letter names. The first two letters refer to the data channel, such as Ht (Height), Am (Amplitude), Ph (Phase), and ZS (Z Sensor). These two letters are then followed by the letter T for trace or R for retrace. HtT would be the height trace signal, HtR the height retrace signal. The number of tabs present in an image depends on which channels were saved when the data was collected.
  - **Modified Data** Modified or processed image layers include an asterisk in their name and are located to the right of the raw data tabs. Since it is possible to have multiple modified layers of one original raw layer there is a number included in the layer name. For example, the first modified layer of HtT would be HtT\*0. Please see [Section 7.2.8](#) on page 68 for more on this subject.

- 7. Command Functions** The *Commands* menu (Figure 7.3 on page 36) offers a large variety of functions, most of which are described below. Some of the more frequently used functions are repeated as programmable buttons to the right of the *Commands* menu. The default buttons are shown in the margin next to the relevant command.

**Display Panel** This opens a window called the Display Manager, which includes a list of all the image display windows open at the moment. It allows you to bring any window to the front and to close, tile, or stack selected windows.



**Modify Panel** Functions such as flattening, plane fitting, masking, and filtering are performed via the Modify panel. See Section 7.2 on page 39.



**Analyze Panel** This is where roughness measurements, line sections, histograms of image data, and particles are created, displayed, and generated. See Section 7.3 on page 69 for a more in-depth discussion.



**List Panel** This panel allows you to open images, change directories, save images, and rename extracted images. See Section 7.1.2 on page 35.



**Extract Layer** Can export a copy of an individual image channel to facilitate custom processing. See Section 7.4.2.1 on page 85.

**Insert Layer** Reinserts an extracted layer into the display window. See Section 7.4.2.4 on page 88 for examples.

**Delete Layer** Deletes the current image layer.

**Show Note** This is where many parameters associated with the image can be viewed.



**Image Name—Show** This button, which toggles between 'Show' and 'Hide', shows or hides the image name within the display window. It can be useful for graphic exporting.

**Change User Names** Allows you to change the names of images stored as data from the BNC connectors on the front of the controller. "UserIn0" can become something more descriptive, such as "Photocurrent".

**Save Image** Saves any modifications to the display window, such as modified data or inserted layers in new tabs. See Section 7.4.1 on page 85 for more details.

**Save As...** After requesting a new name for the selected image, this function saves the image to a location specified by the user.

**Save then Kill Image** Saves the image, including all associated layers and tabs, then closes the window.

**Kill Image** Removes the image from memory.

**Export to Layout** Appends the current image to a notebook style page. This process is described in Section 7.4.4 on page 91.



**Tiff exports:** Tiff exports the current layer to a directory of choice as a separate image. It includes numbers indicating the resolution of the new file it will have (1, 2, or 4 times the number of pixels as the original image). These files can get very large but are great for graphics that need to be displayed at larger scales, such as on posters. Note that you can also choose *Edit > Export Graphics* from the main menu. This gives you other sizes and file formats as export options.

**Set Layout Fields** Shows a list of display parameters for the 'Export to Layout' function.

**NanoIndenter Tip Area Function** Shows a graph that will calculate Tip area from Z sensor images of either nanoindenter tips or indentation. These tip area graphs can be used in the analysis of indentation force plots on the Elastic Tab of the force review software.

**Crop Image** Once zoomed in on an image, this command discards all data outside the current viewing area. This is done to all layers of the image. Please follow the steps in [Section 7.4.3 on page 89](#), as they include a protocol to prevent the permanent loss of raw data.

**Help** Links to more detailed descriptions of these items.

**New Appearance** Toggles between the new and classical appearances of the 2D graphs.

**Argyle 3D** This button launches a 3D view of the data. Detailed in [Chapter 11 on page 132](#).



For Display Windows with many image channel tabs, the small arrow buttons to the left of the tabs allow you to move left or right through the image layers. If there are only a few tabs, these are grayed out and inactive.

### 7.1.3.1. Reprogramming Shortcut buttons


**Tip** The shortcut buttons in the Display Panel ([Section 7.1.3 on page 36](#)) are programmable.

Any function from the *Commands* menu can be assigned to any of the buttons. For instance, the 'Note' button, marked with the letter N, can be replaced with a 'Crop' button, labeled with the letter C. You can reassign one of the shortcut buttons as follows:

- **Ctrl+click** on any of the six buttons.



- This will produce the dialog box to the right.
- Select the button you want to reprogram.
- Enter a new name for it; up to two letters can be displayed on the buttons. Usually, the name is based on the initials of the function.
- Choose the desired command from the drop-down list.



## 7.1.4. AR Thumbnail Viewer




Figure 7.4.: AR thumbnail viewer extension for Windows

On the Windows operating system side of the software, the AR thumbnail viewer is useful for finding a specific image within Windows File Manager. This feature is installed as a Windows extension along with the AR software.

1. Right-click on an image file, and choose *Properties*.
2. A panel similar to Figure 7.4 on page 39 will appear. The Asylum Research tab will have a larger thumbnail of the stored image, and the lower portion of the panel has the pertinent scan parameters.
3. On the dropdown menu, you can view any of the data channels, including amplitude, height, and phase.

## 7.2. The Modify Panel



M stands for the Modification of raw image data. Usually this entails image processing of the individual image layers to enhance features. The Modify panel can be opened from:

- The menu: *AFM Analysis > Modify Panel*
- The 'M' button on any offline image.
- *Commands > Modify Panel* on any offline image.



## 7.2.1. Flatten tab




Figure 7.5.: Image Flatten controls.

**When to use:** When an AFM image is captured, there are two main things that affect the scan from line to line. One is thermal drift, which smoothly affects the recorded height of the image. The bad part about thermal drift is that it usually moves up and down in a seemingly random fashion. The other thing that can affect the scanning is having the probe change its behavior suddenly. It can break the tip off or pick up a piece of debris, and suddenly the measured height is changed by more than the actual range of the image. Both of these can be dealt with by flattening. Since you are scanning at 1 Hz, the time difference between two adjacent points horizontally is around 2 mS, but the difference between two points vertically is 1 S; flattening is designed specifically to deal with this problem.

**How to use:** Flatten Order 0 removes an offset from each line so that all of the lines have the same average value.

Flatten Order 1 removes an offset from each line and also removes the slope of the line.

Flatten Order 2 removes offset, slope, and also a second order polynomial fit to the line.

Flatten Order 3 removes offset, slope, and a third order polynomial fit to the line.

Use Flatten on flat samples that show height varying in the slow direction, often caused by tip changes or thermal drift. Longer capture times often require flattening of the image. Start with a zeroth order flatten and look for artifacts or improvements. If the scan was slow, a Flatten Order 1 might make a noticeable difference also, as thermal drift can cause individual lines to have slope.

In general, the flatten must be done on a region of the sample that is on a continuous plane. Raised areas or pits must be isolated from the flatten regardless of the image.

Note also that the images are saved with a first order flatten on them. To restore an image to its raw form, on the Modify Panel (in either Flatten or Planefit tab), click the 'Ultra Restore Layer' button.

## 7.2.1.1. Flatten Example, features on flat substrate

**Overview** For a successful application of flattening with minimal data transformation, flattening should only be performed on portions of the AFM image that represent roughly a continuous plane. In this case, we presume the background to be a relatively flat substrate. Nonbackground features are excluded from the calculation with a mask, thereby preventing actual features from being used to calculate the fit of the polynomial. The use of masks is discussed in [Section 7.2.4 on page 59](#). Here we perform a standard iterative process that works well for a wide variety of images:


- Flatten 0
- Mask
- Flatten 1


In this particular example, we use an image of DNA on mica. You can download this image file from here: <http://www.AsylumResearch.com/Files/Data/FlatteningExample1.zip>

1.

**Starting image layer:**

- With an image open, select the desired image channel.
- From the large offsets between lines, you can tell that this image will require at least a flatten 0.




- Click the  button to open the Modify Panel.
- Set the 'Flatten Order' to 0.

4.

**Perform the Flatten:**

- Click 'Flatten'.
- Notice that there are flattening artifacts (black streaks) to the left and right of taller features.




- Go to the Mask tab and select *Iterative* from the *Calc Method* dropdown menu.



6.

**Create the first mask:**


- Click 'Calc Mask'.
- Make sure *Fill Mask* is still selected to show the masked pixels. The red pixels highlighted in this way will subsequently be masked out and ignored.



7.

**Erode the mask [optional]:**


- Click 'Erode Mask' once to eliminate small, undesirable spots that have been masked.
- Notice the mask covers less of the features.




8.

**Dilate the mask [optional]:**

- Click 'Dilate Mask' twice to enlarge the mask.
- Note that eroding and dilating are not always reversible.
- Redo the mask with an iterative mask to restore the image to its state in step 6. This is necessary to continue the example.





9. **Flatten.**
- Select the Flatten tab.
  - Select a Flatten Order of 1.
  - Click the 'Flatten' button.



**Note** This procedure works because it keeps the order to 0 or 1. Sometimes it is useful to go to a second order flatten in order to calculate the mask. Be aware, in this case, that you must undo the second order flatten between creating the mask and applying a first order flatten to eliminate second order modifications to the data.

10.



**Review: [optional]**

- Create 3D images of the starting point and the end point.

### 7.2.1.2. Flatten and Planefit Example, Lattice steps

#### Overview


Sometimes it is useful to use a combination of planefit and flatten. For example lattice steps can prove to be difficult to correct. It will require an iterative process of masking and correcting, while trying to minimize the distortion of the data. The use of masks is discussed in [Section 7.2.4 on page 59](#), and the use of planefit is discussed in more detail in [7.2.2](#). Here the outline is:


- Manual mask
- Planefit X1
- Improve mask
- Planefit X1
- Assess if Flatten 0 is required

- Assess if Flatten 1 is required


In this particular example we will be using an atomic step of gold, sample courtesy of Nir Kampf of the Weizmann institute. You can download this image file from here: <http://www.AsylumResearch.com/Files/Data/FlatteningExample2.zip>

1. **Starting image layer:**
  - With an image open, select the Z sensor channel.



2. Click the  button to open the Modify Panel.

3. **Hand draw a mask on one side of the image:**
  - Go to the Mask tab, click on 'Include Points' and then drag a box from 6 to 8  $\mu\text{m}$  in X and spanning the entire Y axis.
  - Click 'Make Mask' on the Mask panel.



4. Go to the PlaneFit tab.
  - Set the *Offset Uses* to Mask. This will set one level of the lattice to zero, otherwise zero would be between the two levels.
  - Do a first order planeFit in just X.


5.

**Create a better mask:**

- Go to the Mask tab.
- Set the *Calc Method* to Bimodal.
- Click 'Calc Mask'.

-OR-


- Click in the middle of the histogram to set the level between the two peaks.



6.

**Planefit X 1 again:**


- Optional: Do an Ultra Restore to undo the previous modifications to the image, including the saved planefit done to the image.
- Back to the Planefit tab, do another first order planefit in XY.
- Click 'Auto' scale (on image).



7.

**Do a 0 order flatten [optional]:**


- The rest of the instructions in this example are optional. You need to look at your image to figure out if you need to go further or not.
- This particular image shows some discrete shifts in the slow axis, a couple from 0 to 2  $\mu\text{m}$ , and clear one at about 3.5  $\mu\text{m}$ , these will require at least a zero order flatten.
- Make sure that the *Offset Uses* is set to mask.
- Try adjusting your mask again, use the histogram, and mask off most of the histogram.
- Look at the other side of the lattice (doing a flatten will make the visible side look good, but could make the other side worse).
  - Click on 'Inverse Mask' on the Mask panels.
  - Click on 'Auto' on the image to adjust the scale.



8.

**Mask and flatten 0 again:**


- Use the histogram to make a wider mask and make sure you cover the step. But if you get mask showing up on the far side of the image (e.g., the right side is mostly masked, and the left edge of the image is also masked), then that is going to go badly. Try adjusting your mask more, though you may need to undo flattens to fix the mask.
- Do a 0 order flatten again.
- This should do a good job of getting rid of the offsets between scan lines, though there could still be some slope to the lattice steps.



9.

**Mask and Flatten 1:**

- Use the histogram to make a wider mask. Make it really wide to include nearly all of the histogram but leave much of the image unmasked.
- Do a 1st order flatten.
- Look at the other side of the lattice again.
- Take a section across the image.

**7.2.1.3. Long involved flatten explanation**


This is an explanation of flatten where we look at the effects of flatten on just one line.

- The **black line** is the original,
- The **red line** has had the flatten applied to it.

1.

**Flatten 0:**


- This is a line from an image that has Flatten 0 applied to it.
- You can see that the shape and everything stay the same, just the offset is changed.



2.

**Flatten 0 vs Flatten 1:**


- Next is a comparison between Flatten 0 and Flatten 1.
- Though the amount of bow makes it a little hard to see, the tilt has been removed in the Flatten 1 version. This is more obvious in a flatter scan line.



3.

**Flatten 0 vs Flatten 1, cont'd:**

- Easier to see difference between Flatten 0 and Flatten 1, with different data.






4.

**Flatten 1 vs Flatten 2:**


- This is the difference between Flatten 1 and Flatten 2.
- Our scanners don't usually have too much bow; this scan line had bow artificially added to it.
- You can see some bow in the Height layer if you have a lot of tilt in the fast direction.



5.


**Flatten 2 vs Flatten 3:**

- Next is the comparison between Flatten 2 and Flatten 3.
- It is not unusual if there is a lot of bow that there also more complicated distortions.



6.


The actual unmanipulated data used in the above examples. The bow in our microscopes for the most part come from tilt in the sample in the fast direction, so whether the bow is up or down depends on whether the scan was ascending or descending.

**7.2.1.4. Making a manual mask**

Flattening works best on surfaces that are very flat without big features, which is good because that are the conditions that normally need to be flattened the most. If there are big features, then a mask usually has to be applied to keep the big features from distorting the image.


1.

This scan is the ZSensor layer of a sample that has a raised surface and big particles on the right side of the image. It has been saved with Planefit 1, which means that the vertical offset and X and Y tilt have been removed. You can see from the wavy looking height on the left side of the image that it does need to be flattened.




2.

Also, because of the big particles on the right side, the plane fit actually introduces tilt into the flat left side of the image. This is a line through the image that misses a big particle on the right.




3.

So, we need to either do a Flatten 0 and then a correct X Plane fit 1, or just go ahead and try a Flatten 1. Well that didn't go well.



4.

To see what went wrong, let's look at an individual line. The Flatten 1 algorithm does a linefit to each line of data and then subtracts that fit from the line. In this case, the linefit is skewed by the big particle, so the subtraction of the fit messes things up. What is happening is mathematically correct, but you don't really want to be doing this.






5.


To fix this, we need to apply a mask, which is something that limits what parts of the image we are going to use to calculate the fit. Masks limit what is used for the fit, but the resulting fit is always applied to the whole image disregarding the mask. We could start over, but since we are still planning to use a Flatten 1, there is no need to. As long as you are going to be using the same level of fit or higher, you don't need to go back to the raw image. If you are going to just use PlaneFit, and you have used Flatten, then you need to start over. Or, if you want to use Flatten 0, and you have used Flatten 1, then you need to start over. But if you are using the same or higher level of modification, you can just keep going.

This image has a fairly obvious flat surface that we can use for our flattening. Click on the 'Include Points' button, and then draw a box starting outside the image over the left side of the image. After you click the 'Make Mask' button, it will look something like this, depending on what the transparency on the Mask tab is set to .




6.


Now it's time for Flatten 1. That looks better. But now you can see that a couple of particles on the flat area are affecting the flatten. Click on the 'Exclude Points' button and circle them, then click the 'Flatten' button again. If you are in the middle of drawing things, and you click the 'Flatten' button, it triggers the 'Make Mask' button for you.




7. This looks good.



8. Same image with the mask reset. Usually at this point, it would be a good time to adjust the Range so you can see the image better; but if we did that on this image, you could really see how noisy the image is from being captured with an AFM sitting on a random table. So we won't do that.



9. Here is a look at the line that we looked at before now that the fit is correct. The part that we fit to is now flat. The flat part is a little below 0 because the offset of the image is calculated doing the whole image, and the big particles contribution means that most of the image is below 0. If during the flatten (or afterwards if the mask is still there) the 'Offset Uses Mask' button is clicked, then that flat area would be at 0. On most images with big particles, you want to use a fancy calculated mask; but, in this case, the part of the image where the particles were was also raised, so we didn't want anything from there anyway.



## 7.2.2. Planefit tab




Figure 7.6.: Image planefitting controls.

**When to use:** Planefitting is useful when you have sample that is known to be flat but has a tilt to it. Thermal drift and sample mounting are two common reasons for a straight sample to give a tilted height image. By fitting and subtracting a plane to the entire image, the image can be made flat. Measuring topographical features then becomes easier since there is a relevant zero height from which to measure.

Planefitting an image works on the whole plane of an image at once instead of a line at a time like Flatten. Flatten is a more intense modification, if you have used Flatten on an image you do not need to do a Planefit in Y, and if you have used Flatten 1 then you do not have to do Planefit 1 in X. If you have a channel that was captured with a Save Planefit of Flatten 1, which is the default setting for Height, then you need to do an Ultra Restore Layer on it first.

Planefitting is more tolerant of big particles and other things that would need a mask if you are flattening. The particles average in with everything else better since they aren't being isolated in a single line. A flat surface with big particles evenly distributed will planefit just fine, where it would require a mask to flatten correctly.

**How to use:**

Plane fitting is very similar to flattening. It fits the unmasked data to find and subtracts a plane of specified order.

- 0-order is not very common, as this option simply takes the offset of the unmasked image.
- 1st-order is the most common. It subtracts a flat, sloped plane from the image.
- 2nd-order is much less common. Care needs to be taken with this modification as it has much more potential to alter the topography. In this calculation, a second order polynomial plane is subtracted from the image.

- 3rd-order probably should not be used for any publication. It is included in the software just for tinkering with data. In this calculation, a third order polynomial plane is subtracted from the image.
- Histogram is a much more complicated algorithm; it is more of a recipe that seems to do a good job on a wide variety of images. It:
  - Does some filtering on a temporary copy of the image.
  - Calculates a mask based on the histogram of the filtered data.
  - Does a 1st-order plane fit to the image.

### 7.2.2.1. PlaneFit Example


**Overview** In this example, we go through the various plane fit orders and compare their results on a highly curved sample.

We will be using an image of human hair. You can download this image file from here: <http://www.AsylumResearch.com/Files/Data/PlaneFitExample1.zip>

1.

#### Starting image layer:


- This is the original Z sensor image of a human hair.



2.

#### Do a 0 order plane fit:


- Note that a new layer with an average of 0 was created.



3.

**Do a 1st order plane fit:**


- Since we are increasing the order, we can simply do this on the modified layer.
- Note that the left edge of the hair is lower than before.



4.

**Do a 2nd order plane fit:**


- Note that the image is drastically different.
- This can be helpful in that details on the surface of the hair are clearer, but caution is required as those features could be artifacts from the 2nd order plane fit.



5.

**Do a 3rd order plane fit:**


- The 3rd order plane fit is similar to the 2nd order plane fit, but more extreme.



6.


**Histogram plane fit:**

- You first need to undo the plane fit. Click 'Undo Plane fit' or 'Restore Layer'. Note that the latter option will undo modifications prior to the plane fit.
- Set the order to *Histogram* and click 'XY'.
- Note the similarity to the 1st order plane fit. The hair was masked out, so the side of the hair is more planar as a result of this modification.

**7.2.2.2. Plane fitting a step**


1.

An example of where you need to mask for plane fitting is if you want to image a step correctly. The shown image is of a gold step on a substrate. The sample has been imaged correctly, with the scan crossing the step in the fast direction and enough flat area on both sides of the step for plane fitting. It has been saved with Plane fit 1, which is the default save plane fit for Z sensor.



2.


But if we look at a line in the middle of the image, it doesn't look right. Since the image had a first order plane fit applied, but it has a large step in the middle of it, the plane fit has leveled out the high and low regions on either side of the step.






3.

To correct this, we need to do a plane fit on just one of the flat parts of the image. We'll use the substrate as our plane, as it looks smoother than the top and sometimes deposited films have a slope near their boundary. Use the Include Points button and draw a box on the substrate.




4.

Then a Plane fit 1 in X. We'll change the scaling a little so we can see the top of the step.




5.

Let's look at the same line that we looked at earlier. Compared to the earlier line, you can see that the substrate is now flat, and any slope in the film should be real.




6.

But looking at the image, it looks like there is a horizontal step in the image where the probe changed properties. It is hard to see with image scaling that we have now, so let's change the scaling to see the substrate clearly. We could just change the Offset manually, or we can put our mask back on and check the Offset Uses Mask button on the Flatten tab. This sets the offset by what is not masked instead of using the whole image. Once we do that, we can see there is a step and some pretty serious thermal drift.




7.

So, let's fix that by doing a Flatten 0. We have already done things in the X direction, so we do not need to do a Flatten 1. If we did, that would actually work fine. That gets the substrate flat.



8.

Now, let's look at the whole image again. Looks good.



### 7.2.3. Erase tab




Figure 7.7.: Erase controls.

**When to use:** Abnormal scan lines can make adjusting your color scale difficult, as well as be problematic in masking and skew some analysis (such as sections). Note that erasing too many lines in publication grade data can work against you; experienced SPM users expect aberrant noise lines in publication images. Please use this feature with discretion. Appropriate use may include the



removal of lines caused by disturbances in a lengthy image scan (such as your lab mates slamming doors while the best image of your life is being collected).

**How to use:** Select the line or lines that you wish to have blurred out, adjust the width to make sure the lines are covered, and click 'Do It'. The marked lines will be replaced with averages of the lines remaining on either side of the block to be erased.


1. With a Display Window image open and the desired image channel as the forwardmost tab (for example, Height Trace (HtT) in Figure 7.3 on page 36), click the 'M' button to open the Modify Panel.



2.

**Locate the bad line:**


- Identify the offending line.
- In this case, it shows up particularly well in a 3D view of the data.



3.

**Zoom in:**

- Expand a small area around the lines to be erased.
- Do this by selecting the area of interest in the image and then right-clicking on the area.
  - Select *Expand* or *Vert Expand*.




4. Go to the Erase tab of the Modify Panel.

5. Click the 'Draw Lines' button.


6.

**Select the offending line:**

- While holding the left mouse button down, place the pointing finger cursor over the scan line you wish to remove, then release.
- Notice that the line goes from red while positioning and to blue when set.



7. If necessary, increase the *Line Width* value (Figure 7.7 on page 57) to cover more lines.
8. If the placement of the line is unsatisfactory, click the 'Clear Lines' button to remove them. Repeat the attempt.
9. Click the 'Erase Lines' button. The line will disappear from the image, indicating that the erasure has taken place.
10. With the cursor over the image, use **Ctrl + A** to restore the original zoom and assess the results. It is also possible to zoom out by right-clicking and selecting *Autoscale Axes*.


**7.2.4. Mask tab****Figure 7.8.:** Image masking controls.

**When to use:** Flattening (Section 7.2.1 on page 40) and Plane Fitting (Section 7.2.2 on page 52) should be applied exclusively to the part of the image that is known to be planar. Image masks allow the non-planar areas to be omitted during the fitting. The mask can also be used for image analysis exclusively inside or outside the mask (see Section 7.3 on page 69).

**How to use:** Masking is used to exclude some portion of the image, typically by setting a threshold value to include the pixels above or below that threshold. There are a couple of ways to make finding the threshold value easier, as well as a number of other methods to fine tune the resulting mask. The general steps for a simple manual mask:

1. Drag the blue vertical line through the histogram graph to set the threshold for the mask.
2. Optionally, invert the mask using the 'Invert' checkbox.
3. Optionally, use 'Range' to mask a range of values in the middle of the histogram; it is also possible to unmask this range by checking 'Invert'.

Next, a simple automated method:

1. With a Display Window image open and the desired image channel as the forward most tab (for example, Height Trace (HtT) in Figure 7.3 on page 36), click the 'M' button to open the Modify Panel. 
2. Go to the Mask tab of the Modify Panel (Figure 7.8 on page 59).
3. Confirm that "Iterative" is chosen as the *Calc Method*. This method is a good mask to start with to find the proper threshold value.
4. For typical AFM samples with a flat background and protruding features, check the 'Inverse' checkbox. For a sample with pits, do not choose 'Inverse'.
5. Click the 'Calc Mask' button. A threshold will be determined that may work for your image. Notice that the *Calc Method* has changed to "Manual", and the determined threshold has been entered.
6. If the mask seems close to where you want it but should encompass a bit more adjacent area, click the 'Dilate' button. This will add pixels around the existing perimeter of the mask. Click 'Erode' to shrink the masked area. You can also eliminate lots of small masked areas by eroding one or two steps, then dilating back again. Only large blobs will survive that process. For examples of this, see Figure 7.9 on page 61.
7. If the mask seems completely off, you can always edit the threshold and range values and manually recalculate the mask. Alternately, use the histogram to set the threshold manually, as described in the simple manual mask above.

Once you have your mask, flattening and plane fitting will only use the unmasked areas for the fit and will apply the subtracted fit to the entire image; the unmasked areas are affected, but they are not used for the fitting.

Filters such as smoothing or blurring will only be applied to the unmasked area. Masked areas will remain exactly as before the filter.

#### 7.2.4.1. Mask calculation methods

There are many different mask algorithms to choose from:

**Manual** The user may enter a value into the Threshold range. This will put the mask above or below that value, depending on the state of the 'Inverse' checkbox. There is also a histogram on the Mask panel that represents the image data. You can drag a blue vertical line on this histogram to set the threshold graphically.

**Iterative** This calculation automatically picks a Z threshold range using an iterative method.<sup>1</sup>


**Bimodal** The mask is calculated based on the assumption that the image histogram is a simple bimodal distribution<sup>2</sup>.

**Adaptive** Adaptive thresholding evaluates the threshold based on the last 8 pixels in alternating rows<sup>2</sup>.

**Fuzzy Entropy** This function uses entropy as the measure for fuzziness<sup>2</sup>.

**Fuzzy-Mean Gray:** Fuzzy thresholding uses a method that minimizes a fuzziness measure involving the mean gray level in the object and background<sup>2</sup>.

Typically, one uses manual or iterative. Playing with the more exotic methods is encouraged. They can sometimes surprisingly lead to the result you were looking for.



**Figure 7.9.:** Creating image masks: A) after zero order flatten (i.e., no mask); B) after Iterative mask; C) Panel B after Erode mask button click; D) Panel B after Dilate Mask button click.

#### 7.2.4.2. Hand drawing masks

When all the automated mask algorithms fail, or if they work except for one area, it is possible to add or subtract areas from the mask by hand. This can be done either to an existing mask or to an unmasked image.

1. To add points to the mask, click 'Exclude Points'. This will cause a set of drawing tools to pop up to the left of the Image Display window (Figure 7.3 on page 36). The points will be excluded from subsequent operations.
2. Select a tool and draw a shape around the area of the graph you want to add to the mask. You can draw multiple areas with multiple tools during this step.


**Advanced Note:** Each tool can be right-clicked to reveal more tools.

3. Notice that the 'Exclude Points' button has turned into 'Make Mask'. Click that button and the outlined areas will be added to the mask.

<sup>1</sup> T. W. Ridler and S. Calvard, IEEE Transactions on Systems, Man and Cybernetics, SMC-8, 630-632, 1978.

<sup>2</sup> This function [ImageThreshold] was designed by Wavemetrics Inc. and is built into the Igor Pro software. Please refer to the Igor Pro Software manuals for more information.

**Note:** To remove areas from a mask, follow the same steps using the 'Include Points' button. They will be included in subsequent image processing.



**Figure 7.10.:** Using 'Exclude Points' masking feature: A) Click on exclude points, and select the free hand tool to make shape; B) Draw the shape on the image; C) Change the tool to ellipse and drag around the height feature to also be excluded—blue indicates it can be moved on image with arrow keys; D) Click on 'Make Mask'. (Data courtesy of Keith Jones, Asylum Research; sample courtesy J. Schlenoff, FSU Chemistry.)

#### 7.2.4.3. Copying masks

In some instances, performing two separate iterative mask calculations on two separate image channels in a data file can give two different results. This can be problematic, especially when overlaying two channels in ARgyle (Section 11.2 on page 133). In these instances, it can be a good idea to copy the mask from one channel and paste it into the other channel.

Another reason to copy the mask is to save processing time when trial-and-error approaches are being employed with plane fitting or in finding a flattening order.

1. Make a mask on one of the image layers using the methods described above.
2. Click the 'Copy Mask' button.
  - a) Note that a pull-down menu appears to the right of the 'Copy Mask' button. This allows you to copy to multiple "clipboards".



3. Select another image layer.
4. Click the 'Paste Mask' button. This will paste the mask on top of whatever mask was already present on that layer.
  - a) If you are using multiple clipboards, there will be menu to the right of the 'Paste Mask' button allowing you to select which clipboard to use.

### 7.2.5. Filter tab




Figure 7.11.: Image filtering controls

**When to use:** Flattening and plane fit are very good “filters” for AFM images, they apply a modification to the data that can remove typically problems in AFM images. But sometimes you need a bit more. The filter tab has more general use filters that can be applied to any image data. In fact, many of these algorithms were developed for photography. But if you are having trouble making out edges, and you know not to use that resulting image for quantitative analysis, it is perfectly reasonable to apply a harsh edge finding filter.

**How to use:** This tab applies a variety of standard image processing techniques to the unmasked part of the image (see 7.2.4). The two filter types primarily used are the NxN and 3x3 matrix filter. At each pixel of the image, they look at the nearest neighbors of that pixel and apply an operation on those neighboring pixels to calculate the new pixel. A very simple 3x3 filter would take an average of the 9 pixels and put that value into the center. A description of these matrix filters is defined better in the software Help menus.

- **NxN filter:** This matrix filter has various filter method types including: Median, Average, Gauss, min, max, and NaNZapmedian.

- **Size:** Defines the size of the  $N \times N$  matrix; the larger the number, the more the pixel will be blurred (influenced by more neighbors).
- **Passes:** Defines the number of iterations that filtering process undergoes. More passes mean more blurring but can sometimes work better than increasing the size.
- **3x3 Filter:** Additional filtering options are available here, including: surface edges, find edges, point, sharpen, sharpen more, and gradient filters in each direction (N,S,E,W).

Consult any standard image processing textbook to see which filter is best for the type of noise or feature you are trying to suppress or enhance. Also consult the Igor Pro user guides, as most of the image filters are built-in Igor Pro functions.

**Note**

To remove spiky noise on surfaces, consider the Median filter. Unlike Gaussian filters which tend to blur everything, median filters can remove noise while preserving edges.

An example of how to apply a 1 pass  $3 \times 3$  median filter on an image:

1. With a Display Window image open and the desired image channel to modify as the forward most tab (for example, Height Trace (HrT) in [Figure 7.3 on page 36](#)), click the M button to open the Modify Panel.
2. Select the Filter tab.
3. Select *Median* from the  $N \times N$  Filter dropdown menu.
4. Set *Size* to 3 and *passes* to 1.
5. Click the 'Do It' button, located just below  $N \times N$  Filter.
6. A new tab will appear in your Image Window with the filtered result.

**7.2.6. FFT tab**


**When to use:** FFT filtering can be very handy when trying to reveal an underlying periodic structure in a sea of noise, as well as removing that periodic noise from images. Doing a fast Fourier transform (FFT) on the image gives you a measure of the frequency components of the image. The noise is removed by masking out regions of the FFT and decreasing their strength. The result and the difference are reported before you confirm putting the iFFT (inverse FFT) back into the original image.

**How to use:**

1.

**Select Image, Open Mod:**


- With a Display Window image open and a desired image channel to modify as the forwardmost tab, click the 'M' button to open the Modify Panel.



2.

**Select FFT tab and prepare the panel:**


- Select the FFT tab.
- Make sure the 'Auto' Scale check box is selected. Now change the *Image* dropdown menu to make sure Auto is selected for all three values (FFT, iFFT, Diff).
- Click the 'Make Graph' button.



3.

**Create mask, run the filter:**


- Click the 'Exclude Points' button.
- Have the 'FillMask' check box selected.
- Choose a drawing tool from the pop-up to the left of the window.
- Draw rectangles or ovals around the periodic peaks in the FFT window.
- Click 'Make Mask'.
- Leave the FFT multiplier at "1.00e-9".
- Click 'Do It'. The two square images to the right of the elongated FFT window will refresh with new contents.





4.


**Result:** The window on the upper right was formed by only the spatial frequency components inside the two selected ovals. All the other rejected (masked out) frequency components are shown in the image below that. The lower image looks largely like random noise. There is a tiny bit of structure left in it, perhaps warranting another attempt at drawing the masking ovals a bit larger. But, this will do for the purposes of this demonstration, which preserved the majority of the low frequency components and gave the image its long range structure. Comparing to the unfiltered image shows that some of the atomic steps near the center of the image are not truthfully represented in the filtered result.



5.

**Reinsert the result into the original image:**

- In the FFT tab (seen in Step 2 above), click the 'Insert Layer' button.
- Click 'Do It' in the next window, and the filtered result is added as a new layer to the original image.



## 7.2.7. History tab




Figure 7.12.: History Tab.

**When to use:** This tab shows what modifications have been done to the current image. This is where various steps and filters that were applied can be undone. In addition, once a series of image processing steps have been performed on one image, they can be copied and applied to a batch of images.

**How to use:**

## 7.2.7.1. How to apply one process to another image:

1. With a Display Window image open, click the 'M' button if you have not already done so.
2. With an image layer that has some image processing applied to it, select the History tab on the Modify Panel. Here you will see a list of applied filters and processes.
3. With all or some of these processes selected, click the 'Copy' button. This will copy the selected steps. A new item will appear under the *Set* dropdown menu, indicating that the selected history may be applied elsewhere.
4. Select the image and layer you want to apply this process to. Select the appropriate history choice from the *Set* menu.
5. Click the 'Execute' button, and the series of steps will be applied to the new image.




**Note** If the process is particularly elaborate, save the current experiment history so it can be used to process future images. If the images have been killed and reopened, the history window will not show what modifications were made. FFTs are not listed in the history.

**7.2.7.2. How to apply one process to many images:**

1. Perform steps 1 through 4 above [7.2.7.1](#).
2. Check the 'Global Execute' checkbox at the top.
3. Click on the 'D' button to bring up the display manager.
4. Select the images you want to operate on. Note that only displayed images are listed here; you may need to display more images from the List panel see [Section 7.1.2 on page 35](#).
5. Go back to the History tab and click 'Execute', and those steps will be applied to all of the images selected in the Display Manager panel.

**7.2.7.3. How to remove one step from a multi-step process:**

1. Make sure that only the steps you do want to do are selected in the history list.
2. Make sure the 'Global Execute' check box is NOT selected.
3. Make sure the 'Restore Prior to Execute' check box is selected. This will undo all the modifications to the current image first, and then go through the history, applying the selected steps.
4. Make sure *Set* is set to current layer.
5. Click execute.

**7.2.8. Prefs tab**

**Figure 7.13.:** Modify Panel preferences

This tab has two complicated controls that do not seem to belong in any other section.

**7.2.8.1. Modify Insert Layers**

This is the checkbox shown in [Figure 7.13 on page 68](#). We already mentioned in *Tab Naming Conventions* on [Step 6 on page 36](#) that modifying a layer places the result in a new layer with an asterisk in its tab. This new layer creation is the default behavior when 'Modify Insert Layers' is checked. If a modified layer already exists for a certain channel, more modifications will overwrite that layer. In other words, filtering a raw data layer automatically creates a scratch copy. Subsequent filtering of that scratch copy will keep overwriting the results of the scratch copy. This preserves your original data while preventing the pileup of countless scratch copies.

**Example:** Your data has Height Trace and Deflection Trace. You apply a flatten to the height layer. With *Modify Insert Layers* not checked, the original height layer would be overwritten. If *Modify Insert Layers* is checked, then you will get a new layer called HtT\*0 which has the flattened height data. If you were to then view the HtT\*0 layer and do a filter on it, the data in that layer will be overwritten, regardless of the state of this control, with the flattened and filtered height data. Finally, you go back to the original Height Trace layer and do a plane fit. The HtT\*0 layer will be overwritten with the original height data, followed by the plane fit; you lose the flatten and filter.

#### 7.2.8.2. Undo Only Last Plane fit

This check box changes how 'Undo Plane fit' works.

- When not checked, everything works the way it used to. When the 'Undo Plane fit' button is clicked, it undoes all the consecutive plane fits applied to the selected image. In addition, when this control is not checked, it is impossible to undo a plane fit once another operation, such as a flatten, has been performed. The plane fit can only be removed if the layer is restored.
- When checked, clicking on the 'Undo Plane fit' button will undo only the last plane fit. If a first order plane fit is applied after a second order plane fit has been calculated, only the first order plane fit will be removed. When this control is checked, a plane fit can be undone at any time, and no other modification will be reversed. The only other way to do this would be to use the History tab (see [Section 7.2.7 on page 67](#)). Finally, if this control is checked, then you can undo the RealTime Plane fit; this function is completely backward compatible. You can load up old data and undo the RealTime Plane fit. There is no other way to reverse the RealTime Plane fit.

### 7.3. The Analyze Panel

The Analyze panel is where roughness, line sections, histograms, and particle analysis of image data can be performed.

## 7.3.1. Roughness Tab




Figure 7.14.: The Roughness tab of the Analyze Panel.

**When to use:** This tab provides statistics on images. When used with the [Section 7.2.4](#) on [page 59](#), it can also provide statistics of the unmasked portions of the image. The topmost image will have its statistics displayed. See the gray help buttons (question marks on right side) for descriptions of statistics. If the image is masked, then the Masked Image column will have the statistics excluding the masked pixels.


**How to use:**

1. With a Display Window image open and the desired image channel as the selected tab (for example, Height Trace HtT in [Figure 7.3](#) on [page 36](#)), click the 'A' button to open the Analyze Panel.
2. Select the Roughness tab.
3. View statistics of the currently selected channel.
4. Click the 'Export Roughness' button to create a text file of the statistics and save it to the directory from which the image originated.




An example of image roughness measurements relative to surface features contained by the mask is shown in [Figure 7.15](#) on [page 71](#). You can download this image from here: <http://www.AsylumResearch.com/Files/Data/RoughnessExample.zip>. Panel A shows an image of an elastomeric mold that contains no image mask. Notice the Masked Image column of the Roughness tab shows Nan, “Not a Number”, because there is no mask. Panel B shows the result of an iterative mask with ‘Inverse’ checkbox NOT activated. The values in the Masked Image column represent

the areas outside the mask that include the lighter areas of the image. Finally, in Panel C, the mask was manually set to include the smaller features of the mold above the plane of its main surface; these statistics represent the image with the 'Inverse' checkbox activated.



**Figure 7.15.:** Roughness Panel statistics: A) Unmasked image; B) Iterative mask results; C) Manually adjusted to mask smaller features at top of pattern.

### 7.3.2. Section Tab



**Figure 7.16.:** The Section tab of the Analyze Panel.

**When to use:** The Section tab allows you to make line sectionsto extract image data. Both straight lines and curves can be defined. The lower portion of the panel provides measurements alon the image section.

**How to use:**

#### 7.3.2.1. Straight Line Sections:

1. With a Display Window image open and the desired image channel as the selected tab (for example, Height Trace, HtT in [Figure 7.3 on page 36](#)), click the 'A' button to open the Analyze Panel.
2. Open the Section tab.
3. Make sure *Mode* is set to *Line*.
4. Make sure the 'Angle' checkbox is unselected.
5. Click the 'Draw' button.






6.

**Draw a line:**

- Click the spot on the image where you want the section to start. Hold the mouse button down.
- Drag to the endpoint of the section and release the mouse button.
- Note the two cursors at the ends of the line section. They can be dragged to other locations in order to reorient the line section.

**Optional** Select the 'Full Width' checkbox if you want the line section to extend to the edges of the image.




7.

**View section:**

- Once the line is drawn, a Section plot appears.
- Drag cursors from the bottom of the trace graph onto the trace; corresponding markers appear on the image.
- Positions and distances between the cursors are displayed in the lower half of the Analyze Panel.

**Optional** Click the 'Export Table' button on the Section Panel if you want to store the information.




8.

**Average over multiple lines:**

- Increase the pixel width of the line to average the section over a broader path. The control is located next to the 'Full Width' checkbox.
- The section graph will update automatically.


**Note** The width uses units such as "n" for nanometers or "P" for pixels.






**Leveling Line Sections** If an image happens to possess some tilt, as is usually apparent from section data, using the Igor cursors to measure a meaningful height difference is not an easy process. One possibility is to perform a first order XY plane fit on the image before making sections, as was discussed in [Section 7.2.2 on page 52](#). An alternate solution based on manual subtraction of a linear background can be executed if the section data are already on screen. The latter is described below.

1. Hold down the 'Ctrl' key while moving the mouse pointer over the curve. An up/down arrow cursor will appear.



2. Click and drag vertically to change the underlying slope of the curve. Then, manually level the background.



**Note** Be sure not to confuse these leveled curves with real, unmodified data. Only the section profiles are altered, and this modification will be undone if the section is updated.

**Adjust the XY Angle of the Section Line on the Image** There are two ways to adjust line section angles with the section drawn:

- Move one of the end cursors on the section line by left clicking and dragging it to a new point.
- Select the 'Angle' checkbox, then manually adjust the angle by typing values into the *Angle* field on the Section Panel ([Figure 7.16 on page 72](#)).

**Note** Once the 'Angle' checkbox is selected, the section cannot be rotated with the cursors.

### 7.3.2.2. Free Hand Lines

Free hand lines are curved lines that are drawn by hand.


1. With a Display Window image open, and the desired image channel as the forwardmost tab (for example, Height Trace (HtT) in [Figure 7.3 on page 36](#)), click the 'A' button to open the Analyze Panel.
2. Open the Section tab.



3.

**Set up for drawing:**


- Select *Free Hand* from the *Mode* menu.
- Click the 'Draw' button.



4.

**Draw the curve:**


- Click and drag to draw the cursor path through the surface features of interest.
- Once the mouse button is released, the line becomes a series of square points.
- [Optional] By click-dragging these points, you can modify the curve.



5.

**Finalize the section:**

- Click anywhere on the image away from the curve.
- The square points disappear, and a section graph will appear on screen.

**Plotting Multiple Line Sections on One Section Plot:**

1. Make a section of an image layer using one of the methods described above.

2.

**Create a copy:**

- In the section graph window click the 'Take a "Snap Shot" of active trace and append to upper graph' button.
- The current trace is copied to a second graph, while the bottom graph is empty and awaits a new section.




**Optional** Double click the upper graph and change its color.

3. Make another section of an image layer as described above. A new section will appear in the lower graph.

#### Overlay multiple sections:

- Click the 'Snap Shot' button again to join the lower graph with the other curve on the upper graph.
- Repeat this process as many times as needed.
- 'Ctrl'+click the 'Snap Shot' button to append the section to an invisible left axis. This is good for comparing sections of different data types, such as phase and amplitude, as the section will not need to share the same axis range.
- 'Shift'+click the 'Snap Shot' button to simultaneously append the marker distances (located in the bottom half of the Sections panel) to a notebook.



**Exporting Sections as ASCII:** Click the 'Edit' button on the Section Graph. This will create an X,Y delimited table from which the columns can be copied and pasted as text into any other program.

### 7.3.3. Histogram tab

**When to use:** This tab allows you to view histograms of image data. Histograms are the result of taking data arrays and counting how many times a value falls within various ranges of values, called bins. This can give an estimation of the probability distribution of the data.

#### How to use:


1. With a Display Window image open and the desired image channel as the selected tab (for example, Height Trace HtT in [Figure 7.3 on page 36](#)), click the 'A' button to open the Analyze Panel.



2.

**Histogram Analysis Panel:**


- Click the Histogram tab.
- Click the 'Make New Histogram' button, and a histogram graph automatically appears.
- If you have a histogram graph already up, you can append to that graph with options to the right of the 'Make New Histogram' button.



3.

**Histogram Result Graph:**

- You can put multiple histograms on a single graph from the Histogram tab. There will be a dropdown menu to the right of the Make New Histogram, listing the existing histograms. Select one of those to append the current data set to that histogram.
- The Select Histogram control on the histogram graph will set which distribution the rest of the controls work on.

**7.3.3.1. Fitting**

- There are three distributions built into the software: Gaussian, Poisson, and Lorentzian.
- Select the model in the *Fit Type* field, then click the 'Fit' button. There are more options available in the dropdown menu to the right of the this button.
- The parameters from the fit are displayed under the 'Fit' button.
- To append the fit parameters to the graph, use the dropdown menu to the right of the 'Fit' button and select *Label Fit*.
- To fit a subregion of the distribution, mark the start and stop with the cursors. Press 'Ctrl'+i' to show the info window, then drag the circle and square to mark the start and stop of the fit region.

**7.3.3.2. Data Modification**

To do offsets of the image data, place the cursors at the positions you want to be zero. Press 'Ctrl'+i' to show the info window. Then drag the A (circle) or B (square) cursor where you want zero to be. Then go to the *Commands* dropdown menu and select "Set cursor A to 0" or "Set Cursor B to 0".

### 7.3.3.3. Other Options

**Edit** This button displays a table of all the values of the image that created the current histogram. You can press 'Ctrl'+ 'A' to select all the values, then copy and paste them into any program.

**Bins** You can adjust the slider to change the number of bins in the histogram. If you have multiple histograms, the 'Link Bin Sizes' checkbox can ensure that all the histograms have the same bin size.

**Mask** Place Igor cursors around the region you want to have masked. From the *Commands* drop-down menu, select "Make Mask from Cursors". Now that there is built-in histogram on the Mask panel, this is not as useful as it once was.

**Crop** You can crop data outside the cursors from the *Commands* dropdown menu. This will remove data from the histogram outside of the cursor range set. This is useful if the wide distribution of a few data points makes it difficult to see details in the distribution.

### 7.3.4. Particle Tab

**When to use:** This tab allows you to analyze particles or features of an image defined by a mask.

**How to use:**

1. Create a mask that defines the particles on an image. To bring up the Mask panel, click the 'M' button on the image window and select the Mask tab.
2. Open the Analyze Panel by clicking the 'A' button on the image window, then select the Particle Analysis tab.
3. Click the 'Analyze Particles' button. When the software is finished analyzing, all of the particles will be selected. To select only the particles you want, hold down the left mouse button and drag the cursor over them.
4. Click the 'Detailed Stats' button. This brings up the Particle Analysis Stats Panel that will display the attributes of the currently selected particles.
5. From this panel, you can graph the attributes by clicking on the buttons to the right. The leftmost button is for a histogram, the next button is for a distribution, and the last two are for plotting Y attributes vs 1 or none X attributes. See [Step 4 on page 83](#).


#### 7.3.4.1. Particle Analysis Example 1: Finding the Length of a Piece of DNA

In this particular example, we will be using an image of DNA on mica. You can download this image file here: <http://www.AsylumResearch.com/Files/Data/FlatteningExample1.zip>. We will use skeleton erosion to find the length of a strand of DNA. This method works best on thin particles with curves, where a straight line measurement is not appropriate, making DNA a good candidate.

1.

**Flatten image:**


- This is the same image used in the flatten example. For details on flattening the image, see [Section 7.2.1.1 on page 41](#).
- For this experiment, the image should already be flattened and it should appear as it does to the right.




2.

**Apply mask:**


- From the Mask tab on the Modify Panel, apply an appropriate mask that covers the DNA completely, as shown in the figure at right.
- For details on applying a mask, see [Section 7.2.4 on page 59](#).




3. From the Particle Analysis tab on the Analyze Panel (Figure 7.17a), click the 'Analyze Particles' button. A progress bar will appear (Figure 7.17b), and a blue overlay will appear on top of the mask (Figure 7.17c).



(a) Analyze Panel



(b) Progress Bar




(c) Blue Overlay

**Figure 7.17.: Analyzing Particles**

4.

**Selection of particles:**


- When analyzing of the particles has reached completion, all of the particles will be selected in a green outline.



5.

**Reset the mask:**


- On the Mask Panel, you can reset the mask to make it easier to see the particles.



6.

**Un-select the particles:**

- Left-click anywhere on the surface to un-select all of the particles.
- Then, move your cursor over to the upper-left strand of DNA.
- When the cursor is over it, the particle will turn green.






7.


**Detailed view:**

- Right-click on the particle and choose *Detailed View of Particle*.
- This opens a detailed view of the particle, as shown in the figure to the right.

**Start Skeleton:**

- Click the 'Start Skeleton' button.
- The particle will change to a lighter orange, as shown in the figure to the right.
- Click on the ends of the particle with the mouse cursor. These two points will be shown in blue.

8. **Note** The skeletonizing process works by eroding the particle until it is a single line thick. It is possible and often occurs where the resultant skeleton will have multiple branches. The two points choose are considered safe and can never be eroded, so they preserve their respective branches, while any other possible branches are completely eroded away.




9. When you like the end points selected, click the green checkmark button.



10.

**Completed skeleton:**

- Another progress bar will appear. When it is finished, you should see something like the figure at right.
- In our case, we see the skeleton consists of 4 segments: a light blue segment at top, a green segment at the bottom, and two tiny segments in the middle (yellow and violet). This happened because the mask itself had a tiny hole in it that caused a loop.
- If you mouse over the different segments, they will highlight with their individual distances displayed at the mouse cursor. Clicking on the segment will show it and other clicked-on segments combined distance in the upper right.



**Note** The software automatically calculates the shortest path and displays it in white underneath the segments. In this case, it is a combination of the blue, violet, and green. If you zoom in, you will see the yellow segment is longer than the violet and not needed to make the shortest path. The combined distance is shown as a tag. In this case, we got 696.81nm long. Your results will vary depending on the end points chosen.

**7.3.4.2. Particle Analysis Example 2: Finding Average Height**


In this example, we will be finding the height of the raised features (particles) of a SEB sample. We will use the modified height layer, which has already been flattened.

1.

**Create Mask:**

- To start, we first must create a mask to define the particles just like in [Section 7.3.4.1 on page 78](#).
- In the figure on the right, the mask was calculated as an iterative mask, then eroded twice to make sure to exclude the boundary portion.

**Note** Since we are only looking at the height of the particles we do not have to check the "Ignore Edge Particles" checkbox. If we were interested in the geometry of the particles we would want to ignore the edge particles.




- Click the 'Analyze Particles' button. After it has finished, you may click the 'Reset Mask' button on the mask panel to remove the mask from the image. By default, all of the found particles will be selected with a green outline.

3. Click the 'Detailed Stats' button. A new panel will come out with a list of attributes. The value displayed next to each attribute is always the average of the selected particles. Since we currently have all of the particles selected, we are looking at the average of them all.


#### Find the Height of Raised Features

- 4.
- We want to know the height of the raised features, and have 3 different values to use:
    - Maximum Z: The largest height found on the particle
    - Center Z: The value of the exact center of the particle
    - Average Z: The height value averaged over every pixel of the particle
  - The first two options are more appropriate for spherical particles. In this case, we will use the Average Z. So at this point we have the Average of all the particles' Average Z. The value is 21.242nm as shown in the figure to the right.



#### Distribution Plot:


- 5.
- If we click the second button to the right of the Average Z value, we can bring up a distribution plot.
  - Notice that over 90% of the particles are below 30nm, with a few that are much higher. These high values are due either to impurities on the sample or to the tip having trouble tracking.



6.

**Change the Range:**


- If we look at the image with a large range, the particles we are interested in become grayed out, but the offending areas remain white.



7.


**Remove Offending Particles:**

- We could simply shift-click on the three offending particles to remove them, but doing it by hand is not always ideal.
- Instead, from the Particle Analysis Panel, click the 'Filter' button. A filter panel (PA Search Panel) appears.
- Click the checkbox next to Average Z.
- Change the filter type to *Less than <*.
- In the high value, enter 40nm.
- Click 'Find Particles'. The software will then select only those particles that have an Average Z of less than 40nm.



8.


- The selected particles should look like the figure to the right, with the offending particles not selected.



9.

**Histogram of the Data:**

- Recheck the PA Stats Panel (see step 4) and note that the *Average Z* value has changed to 20.170nm.
- Click the button just to the right of the *Average Z* to bring up a histogram of the data.
- Change the bin size and then apply a fit of the histogram data to get a Mean and Standard deviation value.



## 7.4. Miscellaneous Operations

Below are some useful operations for image processing.

### 7.4.1. Saving modifications to images

Any image modification can be saved by going into the *Commands* dropdown menu and selecting one of the “Save...” options. *Save as* will ask you to rename the file, then to select the path you want to save the image to. *Save* and *Save then Kill* will both save over the original file, but *Save then Kill* will close the image after saving it.

You can also save multiple images from the List panel. To do so, select the images from the Memory column on the right, and then click on *Save* or *Save as*. *Save as* will ask for a new location in which to save the data.

### 7.4.2. Custom work on image channels


#### 7.4.2.1. Extracting layers:

A copy of an image layer can be extracted into its own image containing just that layer. Math can then be done on that layer, or the layer can be packaged into individual files.

1.

**Select layer to extract:**


- With the desired channel as the selected tab, select *Extract Layer* in the *Commands* menu.
- [Optional] ‘Ctrl’+click one of the six letter buttons to set it to *Extract Layer*, as described in [Section 7.1.3.1](#) on page 38.



2.

**View the extracted layer:**

- Open the List Panel (see Section 7.1.2 on page 35).
- To view this file, double-click *LayerData*.

**Note**

If you intend to extract several layers in succession, you must rename the LayerData name manually between extractions. The LayerData name gets reused and overwritten with every extraction. To rename the layer, select the LayerData file, and choose 'Rename Images' from the List Panel. There is an advanced command line technique for extracting multiple layers without renaming, described below in Section 7.4.2.2 on page 86.

**7.4.2.2. Subtracting Image Layers:**

Subtracting image layers entails a pixel-by-pixel subtraction of one image layer from another. This can be useful for comparing lateral trace and retrace data.


You can download the file for this example here: <http://www.AsylumResearch.com/Files/Data/SubtractExample.zip>

1.

**Extract the first image layer:**

- In the Display Window of the image layer you wish to perform the subtraction on, select the lateral retrace image, and then select *Extract Layer* from the *Commands* menu.
- The images in memory column of the List Panel will now have *LayerData* as a choice.

**Optional** To see the image, select LayerData and click the 'Open Images' button.



2.


- Back in the List Panel, select *LayerData* in the right column. Click the 'Rename Images' button. The New Name dialog will appear; enter any name to finish the renaming process. In this example, the new name is "LatRetrace".
- Alternatively, in the Command line, rename that extracted layer to whatever you want. For this example, the command would be:

```
Duplicate/0 layerdata LatRetrace
```

3.

**Extract the second image:**

- Go to the lateral trace image and select *Extract Layer* from the *Commands* menu.



4.

- In the Command line, subtract one channel from another:

```
Layerdata = Layerdata -LatRetrace
```


-OR-

```
layerdata -= LatRetrace
```

5.

**Reinsert the result:**


- Reinsert the layer (see [Section 7.4.2.4](#) on page 88) into the original image.





### 7.4.2.3. Rotating Images:

In some rare instances, an image must be rotated, typically by increments of 90°.



**Figure 7.18.:** Image Rotation and the Command window

To rotate an image:

1. Extract the desired layer (see [Section 7.4.2.1 on page 85](#)).
2. In the Command line, type:

```
RotateImage(LayerData,Degrees)
```

where Degrees is the angle you want to rotate the image.

### 7.4.2.4. Reinserting layers:


Once an extracted image layer has had some changes applied to it, it can be reinserted into the original image.

1. Only the last extracted layer can be reinserted.

2.

**Open the original image:**


- From any Image Display Window select *Insert Layer* from the *Commands* menu. The software automatically sends the LayerData back to the original image.



3.

**Insert Layer dialog box:**

- The Insert Layer window will pop up.
- Either choose to overwrite the original layer from which the layer was extracted, or to direct the data to a new modified layer tab.
- Notes can be added to the new layer in the large text box.
- Click 'Do It'.

**7.4.3. Cropping Images:**

Cropping images is relatively straightforward. This section will discuss basic image cropping by eye first, and precision cropping second.

- Either copy the .ibw file in Windows or extract a layer (Section 7.4.2.1 on page 85) before cropping an image; cropping will delete the excess data, so it is generally safe practice to crop a copy of the data lest in order to protect raw data.




2.

**Zoom in on the area of interest:**

- Click+drag an area of interest.
- Right-click and select *Expand* or *Square Expand*.
- You can then 'Alt'+drag to pan the view.

**Note** Square Eexpand only works for offline data, not real-time imaging.

**Notice** The scales do not have their origin at zero.




3.

**Crop and Save:**

- Go to the Image Window *Commands* menu and select *Crop Image*.
- [Optional] Go to the Image Window *Commands* menu and select *Save* or *Save As*.

**Notice** The scales have their origin at zero after the cropping. The image here was taken just prior to cropping.

**Precision Cropping:**

There is a way to crop more precisely:

1. Double-click on the number labels of an axis of the image or select *Graph > Modify Axis* from the menu bar. This opens the Modify Axis Panel.
2. Manually enter the cropping area under the Axis Range tab. Note that it is necessary to enter values in meters, so 1e-6 represents 1  $\mu\text{m}$ .
3. Go to the Image Window *Commands* menu and select *Crop Image*.
4. Optional: Go to the Image Window *Commands* menu and select *Save Image*.

### 7.4.4. Igor layouts

Many graphs have a useful item referring to exporting or appending data to a layout (see [Section 7.1.3 on page 36](#)).

- Any data you have worked on in the offline analysis can be put into one of these layouts by clicking the 'Layout' button.
- To save the layout, go to *File > Save Graphics*. In the dialog box, you can determine the name, path and file type of the layout.
- Use the Igor tool box ('Ctrl' + 'T'; upper-left corner) to add text or shapes to the layout.




Figure 7.19.: Igor Layouts

## 8. Force Curve Analysis

CHAPTER REV. 2436, DATED 09/04/2021, 14:34.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

### Chapter Contents


8.1	Opening Stored Force Plots . . . . .	93
8.2	Display Tab . . . . .	94
8.2.1	Loading and Deleting . . . . .	94
8.2.2	Adjusting the Display . . . . .	95
8.2.3	Navigating the data . . . . .	97
8.3	Pref Tab . . . . .	99
8.4	Parm Tab . . . . .	102
8.5	Cursor Tab . . . . .	103
8.6	Modify Tab . . . . .	105
8.6.1	Modifications . . . . .	105
8.6.2	Undo, Restore, and Save . . . . .	106
8.7	WLC Tab . . . . .	107
8.7.1	Force WLC Models . . . . .	108
8.7.2	Getting Started . . . . .	108
8.8	Analyze Tab . . . . .	110
8.8.1	Single Force Plot Tab . . . . .	110
8.8.2	Entire Force Map Tab . . . . .	112
8.8.3	Examples . . . . .	113
8.9	Spot Tab . . . . .	114
8.10	Elastic Tab . . . . .	115
8.10.0.1	Elastic tab hertz example . . . . .	118
8.10.0.2	Elastic Tab JKR Example . . . . .	123
8.11	Online Video . . . . .	129


This Section discusses how to perform basic force analysis on stored force (.ibw) and force map (.ardf) files. Keep in mind, there are not necessarily hard and fast rules / sequences regarding data processing; it depends on the data, and can be a trial and error process. For this reason, the processing techniques are broken into sections, and some examples given.

## 8.1. Opening Stored Force Plots

1.

### Prepare the software

- Launch the AR software to open the Mode Master window. (You can also click the 'Mode Master' button at the bottom of the screen: )
- Click the 'Standard...' tile.
- Then select *Offline Force*.



2.


### Open the Force Review Panel directly:

- Alternatively, you can open the Force Review Panel from the menu bar. Select *AFM Analysis* > *Master Force Panel*

3.

### Select a directory of force plots:

- If there is no force data already loaded in the experiment, the software asks you where you would like to load data from when it is opened.
- For a new directory, click the 'Load' button (bottom left corner of panel) and select the directory.




4.

### Select a force plot to open:

- Select a force plot from the list near the middle of the panel (see Box 2 in [Figure 8.1 on page 94](#)).
- You can make use of the various data display options outlined in [Section 8.2.2 on page 95](#).

## 8.2. Display Tab



**Figure 8.1.: Box 1: Adjusting the Display, Box 2: Navigating the Data, Box 3: Load and Delete**

### 8.2.1. Loading and Deleting

#### Loading

- When the Force Review Panel is first opened, you are asked where you want to load data from. You can reopen this dialog by clicking the 'Load' (folder) icon, located in the lower-left corner of the panel.
- You can also:
  - Double-click a force plot from Windows (operating system) to show it in the force review. These force plots will be stored loose in the "memory" folder.
  - Double-click force plots from the offline image browse graph to load the force plots into the force review.
- The list in the middle of the panel includes the available force plots, sorted by folder. To select one or more, hold down the 'Shift' key for contiguous section or 'Ctrl' for multiple noncontiguous.

- There is a distinction between listed in the software and actually loaded. Loading all of the force plots at once can present many problems with computer memory, so there is a title that specifies how many are loaded and how many are simply listed. Once you need the data (generally for display), it loads the force plot, and it remains in memory until removed. You can load all of the listed force curves into the experiment with the load listed curves button.
  - This is dangerous with large data sets on older computers, as it can most likely result in running out of memory errors.




Figure 8.2.: Load and Delete

### Deleting

- Deleting refers to removing the force plots from the Igor experiment; it never deletes ibw files from the hard drive. Force plots can be saved two ways: directly to memory and/or to the hard drive as ibw files (or ARDF force maps).
- The default is to save force plots to both disk and memory and force maps only as ARDF to disk.

**Note** If you change the default to save only to memory and then click on the 'Delete All' button, you will be deleting your only copy of the data.

- You can also delete subsets of the listed data using the 'Select & Delete' button and the 'Delete Displayed' button. 'Select & Delete' opens a dialog to delete some of the force plots. 'Delete Displayed' will remove the currently displayed force plots.

### 8.2.2. Adjusting the Display

**Y Data** The controls at the top of the panel are for selecting what is shown on the Y axes of the graph. When you select a force plot for display, the text color of the channels are changed as follows:

- **Black** All of the selected force curves have the channel.
- **Blue** Some of the selected force curves have the channel.
- **Grey** None of the selected force curves have the channel.

All of the selected channels from a given axis can be removed with the 'Select None' button, located at the bottom of each list. You can display additional channels on any axis, but it generally works better to give each channel their own axis. You can show up to five axes with the control at the very top of the Force Display Panel.




Figure 8.3.: Y Axis

**X Data** All of the Y data channels are plotted against one common X data channel, which is set with the X axis popup. Typical X data types are:

- **Zsnsr** Heavily filtered Z sensor
- **Raw** Minimally filtered Z sensor
- **Time** Starting from the beginning of the force plot
- **Sep** Tip-sample separation calculated from Z Sensor and deflection
- **Ind** Tip-sample indentation, calculated from Z sensor and deflection

**Note** Both Sep and Ind require that Involts is correctly calibrated.




Figure 8.4.: X Axis

**Section** The force plot is broken into various sections:

- **Extend (Ext)** Where the tip is approaching the surface.
- **Retract (Ret)** Where the tip is withdrawing from the surface.
- **Dwell towards (Towd)** the surface.
- **Dwell Away (Away)** from the surface.

You can control which sections are shown in the middle-left of the Display Panel.

**Advanced** There are more advanced, lesser-used display options on the Prefs tab. For more information, see [Section 8.3 on page 99](#).

## 8.2.3. Navigating the data




Figure 8.5.: Navigating the data


## Viewing Force Plots

- Use the list to directly select which force plots you want to see. It can be useful to collapse the list when switching between folders (force maps).
  - Click the '+ Expand List' button to list all of the folders. The button will change to '- Collapse List'.
  - Click the '- Collapse List' to collapse all folders in the list. The button will change to '+ Expand List'.
  - Click the '±' sign to the left of the folder to expand/collapse a single folder.
- Use the ribbon of buttons to adjust the list:
  - Use the arrow buttons to jump to the next (or previous) force plot in the list.
  - Use the - and + buttons to add or subtract from the number of displayed force plots.
- You can also use the keyboard to navigate the list:
  - Right or Down arrows move to the next force plot.
  - Left or Up arrows move to the previous force plot.
  - 'Page Up' and 'Page Down' move up or down, respectively, in the list by larger step sizes (12 on most computers).
  - 'Home' and 'End' move to the first or last force plot in the list.
  - '+' or '-' (on a keypad) add or subtract one of the force plots to the graph.

**Note** These can all be done from any tab of the Force Review Panel or from the Force Review Graph.

**Movie Player** Use the 'Movie' button to display the FP Move Player to play through the force plots in a slide show. For more detailed information, refer to Igor software 'Tool Tips'. (To toggle 'Tool Tips' on or off, press Ctrl + 1.)





**Figure 8.6.:** FP Movie Player

### 8.3. Pref Tab



Figure 8.7.: Pref tab on the Master Force Panel

#### General

- This panel is for advanced formatting of the Force Review Graph. If you are familiar with Igor, you know there is very large set of display properties that you can set, including line width, line style, and markers.
- This panel has a system of determining what user actions have been done to the force review graph and records those actions to apply again. This panel refers to that as the *macro*.
  - Note that macro has many other meanings in a different context.
- The top of this panel allows you to turn the Macro on or off. The bottom of the panel allows you to update the macro, as well as open the macro in a table for direct editing.
- The 'Reset Macro' button will put the Macro back to the default state. This is a good safety

measure, because the changes you make to the macro will automatically be saved to disk and reloaded the next time you start Igor.

### Line Style

- The default line style is to display the force plots with solid lines when only a single data type is on a given axis and to switch to data type specific line style when multiple data types are on the same axis.
- You can set this to always be solid, or always be set by data type, or to call a user-specified function to set the line style.
- The data-type specific line styles are also set here.




Figure 8.8.: Line Style

### Line Color

- The default line color is to display each section (Extend, Retract, Dwell) as a different color when there is only one force plot on the graph.
- When there are multiple force plots, each force plot is automatically assigned a different color from a color table. The color table can wrap around; so if enough force plots are shown, it will reuse the same colors.
- You can set this to always be color table, always be color by section, or to call a user-specified function to set the line color.
- If you have multiple force plots on the graph, you can double-click on the trace and then change its color. Then on the force graph, there is a split button labeled 'Review Graph'.
  - In the Review Graph dropdown list, you select *Update Color Table*. This will look at the graph, see what has changed, and then update the color table so that the next time you build a graph, it will use your color. You should also see the lower list on the Force Display Panel update with your new color.
  - There is also an *Update Color Table* button on the Prefs tab.




Figure 8.9.: Line Color

**X Orientation**

- 'Keep surface to the left' When selected, most data types will have the surface on the left.
- 'Keep surface to the right' When selected, most data types will have the surface on the right.
- Indentation is the opposite of other data types (surface will be on the right when set to keep on left).
- Time is increasing when the surface is set to the left.




Figure 8.10.: X Orientation

**Miscellaneous**

- Auto Stack Force Axes is on by default, which means the Y axes will be placed so they do not overlap. When this is turned off, then all the Y axes will be overlapping on top of each other.
- The Macro controls are discussed at the beginning of this section (see [Section 8.3 on page 99](#)).




Figure 8.11.: Miscellaneous

## 8.4. Parm Tab

The Parm tab lists the parameters associated with each force plot. Some parameters are editable and, when changed, update the data. Most parameters are simply viewable.

**Master Force Panel**

Display | Pref | **Parm** | Cursor | Modify | WLC | Analyze | Spot | Elastic

**Editable Parameters**

	Line0000Point0003	Line0000Point0004	Line0000Point0005	Line0000Point0006	Line0000Point0007
Involts	129.82 nm/V	129.82 nm/V	129.82 nm/V	129.82 nm/V	129.82 nm/V
SpringConstant	1.00 nN/nm	1.00 nN/nm	1.00 nN/nm	1.00 nN/nm	1.00 nN/nm
AmplInvolts	141.50 nm/V	141.50 nm/V	141.50 nm/V	141.50 nm/V	141.50 nm/V
Amp2Involts	20.00 nm/V	20.00 nm/V	20.00 nm/V	20.00 nm/V	20.00 nm/V
ForceNote					
Locks	NaN	NaN	NaN	NaN	NaN
TipSerialNumber					
LateralGain	1.00 V/V	1.00 V/V	1.00 V/V	1.00 V/V	1.00 V/V
LateralUnit	V	V	V	V	V
UserCalcForce					
UserCalcBForce					

**User Gains** UserIn9

	Line0000Point0003	Line0000Point0004	Line0000Point0005	Line0000Point0006	Line0000Point0007
Gain	1 V/V	1 V/V	1 V/V	1 V/V	1 V/V
Units	V	V	V	V	V
Name					

FP: Line0000Point0003, Parm: SpringCon Set All Set Selected

Search  Copy

	Line0000Point0003	Line0000Point0004	Line0000Point0005	Line0000Point0006	Line0000Point0007
ARDFForceMap	1	1	1	1	1
NumOfSegments	3	3	3	3	3
Direction	Nan,1,-1,0	Nan,1,-1,0	Nan,1,-1,0	Nan,1,-1,0	Nan,1,-1,0
VerDate	120419	120419	120419	120419	120419
Version	120419B	120419B	120419B	120419B	120419B
XOPVersion	20121.31.4.5	20121.31.4.5	20121.31.4.5	20121.31.4.5	20121.31.4.5
OSVersion	Professional S	Professional S	Professional S	Professional S	Professional S
IgorFileVersion	6.2.2.2	6.2.2.2	6.2.2.2	6.2.2.2	6.2.2.2
XLVDT	-0.30227	-0.28045	-0.25862	-0.2368	-0.21497
YLVDT	-1.4867	-1.4867	-1.4867	-1.4867	-1.4867
ARDoIVCurve	0	0	0	0	0
ForceNote					
TipSerialNumber					
MostNegZvoltage	-10	-10	-10	-10	-10
MostPosZVoltage	150	150	150	150	150
ExtendZ	1.0837e-05	1.0837e-05	1.0837e-05	1.0837e-05	1.0837e-05
RetractZ	-7.2248e-07	-7.2248e-07	-7.2248e-07	-7.2248e-07	-7.2248e-07
StartDist	3.006e-06	3.006e-06	3.006e-06	3.006e-06	3.006e-06
ForceDist	1e-06	1e-06	1e-06	1e-06	1e-06

Setup

Figure 8.12.: Parm Tab

### Editable Parameters

- The top half of the panel deals with parameters that can be changed, mostly calibration values, such as spring constant and involts, or scaling factors for user channels. Only the first five displayed force plots have their parameters shown.
- Once you have made a parameter change, you can apply that change either of following ways to:
  - selected force plots with the 'Set Selected' button

- all listed force plots with the ‘Set All’ button.

**Note** This is very memory intensive, as it tries to load all of the force plots in order to apply the modification. Older computers typically cannot deal with 1000’s of force plots.

### Parameter Search

- The bottom half of the tab reports the rest of the parameters. It works well if you can guess part of the name of the parameter you are interested in, type it into the ‘Search’ box, and the resulting list will be reduced to only those that match the search string.
- The ‘Copy’ button takes the resulting list of matches (or selected lines) and put it in the clipboard. You can still use the keyboard shortcuts when this panel is topmost (left/right arrows go to the previous/next force plot). Only the first five force plots displayed have their parameters shown.

**Note** The ‘Search’ box can be turned on and off in the Setup.

## 8.5. Cursor Tab




Figure 8.13.: Cursor Tab

The Cursor tab is used primarily for marking regions of interest (ROIs) for the Worm Like Chain (WLC) tab (see [Section 8.7 on page 107](#)).

You can set the ROI at any time on the Force Review Graph.


- **To set a new point:** Shift + left-click on the trace at a new point.
- **To select an existing point:** Shift + left-click on an existing point.
- **To deselect an existing point:** Shift + left-click on the selected point.
- **To change the selected point:** Drag the A cursor after selecting the point of interest.
- **To remove a point:** Ctrl + Shift + left-click on an existing point.

The Cursor tab allows you to modify the behavior of how the above Force Review Graph shortcuts work, as described below.

**Limiting** You can limit where you put new markers by segment and by force plot. This can be useful if the data is overlapping and you want to make sure you put the markers on one specific trace. To limit by section, use the list box; the checked segments are segments where the markers are allowed. Turning off a section removes any markers in that segment. You can also limit by trace, but this is only useful if you are plotting multiple force curves.

**Clearing** You can clear markers from traces by:

- Selecting the trace you want to remove from the 'Clear Locks From Trace' menu.
- Turning off sections. The markers will then be removed from the excluded sections.
- Ctrl+clicking on the markers you want removed.



**Figure 8.14.:** Clearing and Highlighting

**Highlighting** You can make the markers on a trace larger by using the 'Highlight Locks From Trace' control.

### Editing

- You can modify a marker by Shift+clicking on it (you may need to zoom in if there is another marker near the one of interest). You should see a circle around the marker. You can then drag that to another location or use the arrow keys (left or right) to move the position of the marker.
- For more advanced editing, the *Locks* field in the editable parms can be used to update the markers. See [Section 8.4 on page 102](#).

## 8.6. Modify Tab




Figure 8.15.: Box 1: Modify Box 2: Undo/Redo Options Box 3: Save Options

### 8.6.1. Modifications

You can remove offsets and lines from the data. These operations work on the currently displayed force curves or all of the force plots listed. Operating on all the force plots listed can be extremely memory intensive. It will need to load the force curves into memory, do the modification, and keep the data in memory. Older machines running Win XP cannot process very large data sets. It is better to break the data into smaller data sets (move groups into sub folders in windows) or only modify the data as needed.

#### Automatic Offsets

- These operations are in terms of X and Y axes, whatever the force plot is displayed as. The offsets find the zero point in the Retract or the extend data (based on the state of the radio buttons to the left). see also: Offset algorithm below.

#### Manual Offsets

- You can right-click on the Force Review Graph to select X, Y, or XY offsets. The position where you click will be set to zero for the axes selected.

#### Line



- You can subtract a line from the data. Press Ctrl+i (i for info) when the graph is topmost and drag the A (circle) and B (square) cursors to define the region you want to fit to a line. The line is extrapolated over the entire force curve and the line subtracted.

### Offset Algorithm

- The Y Offset takes the average of the last 10 points for retract data or the average of the first 10 points for extend data as the zero point. If your force plot was taken using the Indenter panel (or Indentation Master panel), the offset is taken as the min of the dwell towards the surface.
- The Raw (X) offset uses the function GetContactSlope to find the rough estimate of the surface contact point. It does this by fitting the slope of the deflection vs. ZSnsr in small (~7 points) chunks until the slope is  $> 0$  (see function FindSurfaceIndex). It then fits the deflection from 20% force to the contact point and extrapolates that line to where it intersects the zero force.
  - If you do not have deflection data, then this function will not modify to the data.

**Note** This is the simple algorithm. The elastic panel has a much more involved iterative method of better finding the contact point. See [Section 8.10 on page 115](#).

### 8.6.2. Undo, Restore, and Save

These operations all work on:

- The last modified force plot
- A selected subset of the force plot
- All of the force plots listed

Save FMap works on:

- A single selected force map
- All of the force maps listed

**Attention** Typically, operating on all force plots listed is memory intensive, but the **Save As All** operation will remove force plots that it loaded as it goes, and so the total amount of memory used for this operation should be minimal.

### Undo

- These buttons toggle between undo and redo. They undo all of the last type of change.
- For example, you change the InVOLS on a force plot 20 times. Clicking 'Undo' last will undo all the InVOLS changes and revert the force plot to its state before any InVOLS changes have been made.



Figure 8.16.

**Restore**

- These buttons restore the specified force plots to their original state, which is typically the hard drive state. For force plots only saved to memory, there is an extra copy of the data before modifications were done to the data that is used.



Figure 8.17.

**Save**

- 'Save' and 'Save As' save IBW files. If the source is an ARDF, it will convert to the old loose style of ibw files. Be careful with this since you may end up with a large number of files, making things slow; however, it can be useful for exporting the data into other packages.
- 'Save As' asks you once for the root folder you want to save to. You cannot rename files, as it is possible to save many, many files with this operation.
- 'Save FMap' will save ARDF files, which are packages of force plots. If the selected force map is really from ibw files, it will still save as an ARDF, so this is a good way to convert old style force maps to newer ARDFs. Note that Save FMap only works on a single selected force map or on all of the force maps listed in memory.




Figure 8.18.

**8.7. WLC Tab**

This tab is used to fit polymer stretching events to the Worm Like Chain (WLC) model.




Figure 8.19.: WLC Tab

### 8.7.1. Force WLC Models

The Multi Chain checkbox switches between the Single Chain (multiple attachments) and the Multiple chains (single attachments) models.

The difference between these two extremes is in how the stretching events (Hitches) are added together:

- Single Chain, then it adds the WLC response only in the error region.
  - Applicable if each response is NOT influenced by the other stretches.
- Multi Chain, then it adds the WLC response to the entire stretch.
  - Applicable where the observed response is the sum of the individual responses.

The difference between these two cases can be made more apparent by turning on the 'Show Individual Stretches' check box.

If the Single Pers checkbox is selected, all the stretching responses are fit to one persistence length. If you believe that all the stretching events are from the same type of polymer (in particular the single chain extreme), selecting this will force each response to have the same persistence length.

### 8.7.2. Getting Started

How to fit force plots with the WLC model:

1. Find the force plot you want to fit.
2. Plot it as Force Vs. Sep (Sep is short for tip-sample separation). The X axis is controlled with the Force X axis popup menu beneath the Y axis Data Type list, located on the Display tab.
3. You need to define where 0,0 is. To set the offsets, you can either:
  - a) Right-click on the graph to offset X and Y so that the force is zeroed and the tip hits the surface at zero separation, or
  - b) Go to the Mod tab and click on 'Y offset' and 'Raw (X) offset'.
4. Go to the Cursor tab and select that you only want cursors on Ret (short for retract) Force Cursor Sections.
5. Shift + left-click on the graph to define the stretching regions. You want to define the error region of the fitting function; so pick one point where the polymer ruptures its attachment to the tip and pick the other point at the lowest extension as the data looks like it will be fit. You can fit multiple stretching events in a single force plot. If you want to do that, pick two "locks" (cursors) to define each stretching event (the software calls polymer stretching events "hitches").
6. Open the WLC tab. Make sure that the initial state of the parameters is close to the data; the fitting function often throws up its hands in despair if the initial guess is too far from the data:
  - a) The **pers** is the persistence length; it is basically a measure of the polymer's stiffness. Lower persistence lengths have much more non-linear responses.
  - b) The **extension ratio** is the fraction of the contour length that the polymer chain is extended. The **fit parameter** is the extension ratio where the polymer chain ruptures its attachment to the tip (max extension ratio). Therefore, you can get the contour length of the chain from the rupture length divided by the max extension ratio.
7. Once you have your fit parameters fairly close to the data, you can then click 'Fit'.
8. If you have multiple hitches in a force plot that you want to fit, then there are two questions you have to answer:
  - a) Is the observed response due to one chain with domains or loops or is each hitch the result of stretching a separate chain? If you think you have a single chain, you want to unselect 'Multi Chain' (upper-left corner) Force WLC Models. If you think that each hitch is from a separate chain, then select 'Multi Chain'.
  - b) Do you think that each hitch should have the same persistence length? Generally, this is the case, but it is not required by the fitting function. If you want to fit each hitch to a single persistence length, select 'Single Pers'; if you want each hitch to have its own persistence length parameter, then unselect 'Single Pers'.
9. You can also hold and constrain each of the fit parameters using the hold and constrain checkboxes:
  - a) **Hold** means that the fitting function will hold that parameter constant and not fit it.

- b) **Constrain** means that it will keep the fit parameter between the upper and lower limits. The limit controls show up after you constrain one of the parameters. The constrain option does not always work well; often when the fit fails, it will run the parameter into the limit.
10. Once you like the fit, you can show the individual responses by clicking on the 'Show Individual Stretches' checkbox. You can then label the contour length of each response with the 'Tag Stretches' checkbox.
  11. You can also tweak your fit with the *Max Iterations* and *Fit Tolerance*. You can specify the temperature at which the stretching was done with the temperature control.

## 8.8. Analyze Tab

### 8.8.1. Single Force Plot Tab

This tab runs various analysis (built-in or user-specified) on collections of force plots or force maps to create histograms, plots, and images of calculated parameters from the force data.




Figure 8.20.: Single Force Plot Tab

**Histogram** This check box creates a histogram of calculated values. On the Single Force tab, you can set which force plots will be used to generate the histogram in the upper-left corner. This includes either all of the force plots listed (memory intensive) or a selected subset of the data. The various parameters that can be calculated are:

- **Adhesion** Simple minimum of the force curve, appropriate for clean adhesive ruptures with no polymer tethers.

- **DC Invol** Calculates the slope of deflection vs. Zsensor in the contact region and can use extend or retract data.
- **AC Invol** Calculates the slope of Amplitude vs. Zsensor in the contact region and can use extend or retract data.
- **Indentation** This takes the difference in the indentation at two marked points in the force plot. You need to mark the start and end points of the indentation with the cursor panel controls.
- Pers
- Contour Length
- Extension Ratio
- Rupture Length
- Rupture Force
- Young's Modulus E1
- OliverPharrFit
- OliverPharrStiffness
- OliverPharrReducedDepth
- **OliverPharrArea** These are all report values that were calculated from other methods (WLC or elastic tab). When these calculations are done, they store values in the wave notes, the histogram or scatter plots are simply reporting those stored values.

### Scatter Plots

- You can create scatter plots of a calculate value vs. another parameter. This calculation shares most of its controls with the histogram; the force plot selection and the calculation type are the same as the histogram.
- The X axis of the scatter plot is set by the X parm popup, just under the calculate popup. These are generally parameters that are extracted from the force plot note. For example, you can plot the adhesion vs. seconds to see how the adhesion changed over time.


### Force Average

- This takes the currently displayed force plots and puts them on a second graph that will be averaged along the X axis. It divides the X axis into evenly spaced bins and averages the Y data that falls into each bin.
- You can adjust the number of bins from a slider on the graph.
- It is critical for each of the force plots to be overlapping when you do this. If there is an offset in the deflection from one force plot to the next, that is what the averaging will be showing you.

### PFM Hysteresis Loops

- PFM Force plots taken in the DART spectroscopy controls (driving the DC tip bias with a triangle square wave while driving the AC tip bias with a sine wave) can be further analyzed with this button.
- It converts the raw data in the currently displayed force plot to a hysteresis plot of the average of the data while the bias is off, plotted vs. the DC bias of the previous step pulse.
- More detail on this can be found in ARApplicationsGuide in the PFM chapter.

## 8.8.2. Entire Force Map Tab



**Figure 8.21.:** Entire Force Map tab. **Box 1:** Input and Calculation **Box 2:** Optional Advanced Controls **Box 3:** Starts the Calculation

On this tab, there are controls that create images from force map data. Select the 'Input' and the force map you want to use to calculate the image. Select the function. Below are descriptions for each of the functions:

**FMapCalcAdhesion** Calculates the adhesion by returning the difference of the minimum of the data and the average of the last 10 points.

**FMapCalcJKR2Point** Calculates Young's modulus of the sample using the 2 point JKR method Sun Y. et. al. Langmuir 2004, 20, 5837-5845. The First element of the Parm Wave is the radius of the spherical apex of the paraboloid tip, in meters. The second (optional) element of the parm wave is the Poisson Ratio of the sample [Default = 0, meaning that the result will be in Reduced elastic modulus, not sample modulus].

**FMapCalcHeight** Calculates the height by returning the negative max of the data.

**FMapCalcMax** This function simply returns the max value of the input channel.

**FMapCalcMin** This function simply returns the min value of the input channel.

**FMapCalcXMaxLoc** This function returns the value of the Data where the DataTypeB is maximum. EX: Data is deflection and DataTypeB is ZSnsr, then it would return the trigger point of the deflection.

**FMapCalcXMinLoc** This function returns the value of the Data where the DataTypeB is minimum. EX: Data is ZSnsr and DataTypeB is deflection, segment is Ret, then it would return the Z position of the adhesive rupture.

**FMapTransitionMap** Heavily filters the data, then returns the value of Data where the derivative of DataB is at a minimum.

**FMapCalcPFM(On/Off)** These functions extract the Hysteresis loops of the PFM measurements. The On version takes the data during the On cycles of the triangle square wave, the Off version takes the data when the DC bias on the tip is zero. These functions are specialized in that they build four layers: Coercive Imprint Negative Positive Bias is rising, is the loading cycle; Bias is Falling is the unloading cycle. Positive (Vp) is the bias at which the input channel is at a minimum during the Loading cycle. The negative (Vn) is the bias where the input channel is a minimum during the unloading cycle. The Coercive is  $(\text{abs}(Vn) + \text{abs}(Vp))/2$  Imprint is  $(\text{abs}(Vp) - \text{abs}(Vn))/2$ .

**FMapCalcInvol** This function calculates the inverse slope of the contact region, between 10 and 90% of the Max force for that segment.

**FMapWork** Integrates DataType vs DataTypeB. The default section is ALL, so that it is returning the difference between the trace and retrace segments.

**FMapContactWork** Integrates DataType vs DataType B. Finds where DataTypeB is zero, and integrates the positive values of DataType vs DataTypeB from zero to the trigger point.

**FMapCalcPlasticity** Uses FMapContactWork to calculate the work for the extend data and retract data. Returns the ratio:  $(\text{AreaExt} - \text{AreaRet}) / \text{AreaExt}$ .

- These functions live in XCalculated.ipf. You can add your own functions to UserCalculated.ipf. XCalculated.ipf will be overwritten by the installer the next time you update your software. UserCalculated will not.
- For the most part, when using existing functions, you don't need to worry about Data, DataB, Section and Edit Params, those are for doing custom calculations, selecting the function will set the dataType, dataTypeB, and section to the default state. The auto name and output controls set where the image is going once it is calculated. If you point it to an existing image with a different number of points and lines, the calculated image will be interpolated to fit in the existing image. This way you can calculate images from force maps, and put them into typical AFM images.
- The 'Do It' button will start the calculation. The force map calculation will unload any force plots that it loaded to do the calculation, so that the memory footprint should not grow significantly while doing the calculation.

### 8.8.3. Examples

**Adhesion histogram** In this example, we analyze a batch for force curves to obtain a distribution of adhesion. The software measures the adhesion by taking the difference between the minimum and the zero point. The zero point is defined here to be the average of the last 10 points in the noncontact portion of the curve. The minimum is the lowest point in the retract portion of the curve. Keep this in mind when using this application for polymer stretching force curves, where there may be a larger adhesion away from the surface, this analysis may not be all that useful.

1. In the Analyze tab, select *Adhesion* from the 'Calculate' dropdown menu.
  - a) If ALL the curves are to be analyzed, activate the 'Use all FPs' checkbox

-OR-



- b) Click the 'Select FPs...' button to bring up a panel for selecting the force plots.
2. Click the 'Histogram' checkbox.
3. Click the 'Do It' button. Igor will process and produce a histogram.
4. Within the Histogram window, a variety of information can be acquired or exported. Fits can also be acquired from the 'Fit Type' menu.

**Invol Statistics** Small improvements can be made in the precision of the Invol by taking a statistical average of many invols values. In this example, we calibrate the invols by taking statistics on a batch of force curves.

1. Collect a reasonable number of force curves (say 10-100) on a hard surface. Make sure they all have the same trigger point.
2. Load just those force plots into the force review software.
3. In the Modify tab, click the 'All Y' offset.
4. Plot the force curves as DeflV and note what voltage range you want to do the fit over.
5. On the Analyze tab, select *Calculate DC Invol Ret* or *DC Invol Ext*, depending on which segment you want to use.
6. Two controls will appear, allowing you to set the deflection voltage range to do the line fit over.
7. Click the 'Histogram' check box.
8. Click the 'Do It' button. Igor will process and produce a histogram.
9. Within the Histogram window, a variety of information can be acquired or exported. Fits can also be acquired from the 'Fit Type' menu.
10. Back on the Analyze tab, there are two buttons:
  - a) 'Set to fit' sets the invols of the analyzed force plots to the fitted mean.
  - b) 'Set to each' sets the invols of the analyzed force plots to each have their own Invol value determined when creating the histogram.

## 8.9. Spot Tab

This tab allows you to mark where a force plot was performed relative to the image coordinate system.





Figure 8.22.: Spot Tab

1. Turn the spots on and off by using the 'Activate Spots' check box in the upper-left corner.
2. Select the force plots from the left column. This list is the same list as the one on the force Display tab, and the offline force graph will update to display the selected force plots.
3. The locations of the force plots will be marked on the images selected on the right. Note that only the displayed images are listed. You may need to go to the list panel to open the images (note the List Panel... button which will open the list panel). To have the spot updated on the images, you need to have the 'Activate Spots' checkbox turned on in the upper left corner.
4. The Display Prefs (bottom part of tab, click the '+' to expand) section allows you to set the text, color, and marker properties of the force plot locations.

**Note** Moving the sample or tip in the holder (changing tips) will make this analysis useless. Additionally, when this analysis is appropriate, there may also be small errors arising from thermal drift that are not accounted for.

## 8.10. Elastic Tab

This tab fits single force plots or entire force maps with various elastic models to provide elastic modulus.




**Figure 8.23.:** Elastic tab. Major interface areas are numbered and outlined below.

1. **Fit Input** The upper-left corner of the tab is where you can set it to work on a single force plot or calculate an entire elastic map from a force map. You can also set which segment of the force plots you are working on: extend, retract, or dwell.
2. **Calibration** Here you can override the invols and spring constant stored in the file. Clicking the red 'X' button in the corner of this area will undo the override values and revert the data back to the stored value.
3. **Warnings** This area posts any problems that the software finds with your data. Clicking the buttons here will take you to the Help section that describes what the problems are, their implications, and how to fix and avoid the problem in the future.
  - a) If your input is set to force map, then initially this will be a button with an estimate on how long it will take to load up a handful of the force plots. Once this has been done, any problems found with the examined force plots will be listed.
  - b) Also, if you have analyzed the force map, clicking the 'Warnings' button will bring up a mask that shows which parts of the image had this problem.

**Model Selection Guide:**

- 4.
- This area posts suggestions on which models may describe the data better.
  - It calculates the plasticity factor, the Force ratio, and the Tabor coef; and from those parameters, selects one of the four available models that best works with those parameters, as seen in the flow chart at right.
  - When the Fit Input is set to Entire Force Map, the model selection guide changes to report statistics on the force map that has been examined so far. The colors to the left of these controls reflect the colors of the mask applied to the image.
  - There is a control to check the map that loads and examines force plots randomly to get more statistics. This operation can be canceled at any time from the progress bar that appears.



5. **Models** The tab selects which elastic model you are using. The **Model Selection Guide** (above) can give some guidance as to which models may describe your particular data. Models include:

- Hertz** The classical theory of contact mechanics. From 1886 to 1889 Heinrich Hertz studied how lenses deform under load, and his equations led to the foundation of contact mechanics. See also 8.10.0.1. This model assumes:
  - Negligible adhesion between the tip and sample
  - Strains are within the elastic limit
  - The contact area is much smaller than the radius of the tip
- DMT** Derjaguin-Muller-Toporov; takes into account adhesion outside the contact area. Applies to samples with a small Tabor coefficient.
- JKR** Johnson-Kendall-Roberts; takes into account adhesion inside the contact area. Applies to samples with a large Tabor coefficient. See also 8.10.0.2.
- Oliver Pharr** Relatively new elastic model that fits the unloading curve. More appropriate for samples that deform plastically.

6. **Model Assumptions** This area describes the tip and the sample. The elastic models calculate the reduced Young's modulus, which is comprised of the tip and sample's Young's modulus and Poisson ratios. You need to describe the tip and sample in order to determine the sample modulus from the reduced modulus.

$$\frac{1}{E_c} = \frac{1 - \nu_i^2}{E_i} + \frac{1 - \nu_s^2}{E_s}$$


$\nu$  is the poisson ratio,  $E$  is Young's modulus,  $E_c$  is the reduced Young's modulus  $i$  denotes indenter,  $s$  denotes sample.

- a) The elastic models also deal with pressure. This means you need to know area and tip geometry, which is set in the geometry area. The Oliver Pharr model has more complicated and flexible controls to specify the area.
  - b) The Tip Geometry controls are above the Model Assumptions. There is also an estimation of which geometry will describe your data the best; it is useful as a double-check.
- 7. Initial Guess:** The specifics of this area depend on the Elastic Model. In general terms, you can type in values for the fit parameters and see how that alters the calculated line. Then fitting will take that initial point and try to converge on a better description of the data. If the initial guess is very far from describing the data, the fit will not be able to converge.
  - 8. Results:** Here the sample modulus is calculated from the tip and sample properties and the reduced modulus from the fit.

#### 8.10.0.1. Elastic tab hertz example

This example's instruction set uses a force map collected on a homogenous polyacrylamide gel with a modulus of ~700 Pa. The instruction set will refer to specific force curves from this sample. You can download this sample data set from here:

<http://www.asylumresearch.com/Files/Data/PolyacrylamideGel.zip>.



(a) Master Force Panel, Display Tab




Figure 8.24.: Master Force Panel, Elastic Tab


1. Open the AFM software. In this example, we are using Igor PRO 6.32A and the MFP3D 120804+0806. Earlier versions of the AFM software (especially anything before 101010) may not work or have the features we use in this document. Please check the Asylum support site (<https://support.asylumresearch.com/forum/content.php?4-Software>) for the latest version.
2. Load your Force Map data. From the menu bar, select *AFM Analysis > Master Force Panel*. A new panel called the Master Force Panel should appear, along with the AR Load Path window. Using the browse button on the AR Load Path window, select the folder that contains your force map data, see 8.2.1.

### Open a Force Curve:

3.
  - You should see a list containing all of your force curves towards the bottom of the Display tab of the Master Force Panel (Figure 8.24a on page 118).
  - Simply click the name of any one of your curves, and a graph of the force curve will appear. In this example, we've picked the force curve located at Line 1, Point 1.

-OR-

- If an image came up with the force map, you can click the image to display the corresponding force plot.




4. Select the Elastic tab on the Master Force Panel. The panel should now look like Figure 8.24b on page 119.

**Note** By default, the Elastic tab starts in Single Force Plot mode, as evidenced by the blue Single Force Plot subtab towards the top of the panel. This is a good place to start, even if you are analyzing an entire force map.

**Inspect the Force Curve:**

- 5.
- The force curve will be displayed with Force on the Y-axis and Indentation on the X-axis.
  - Note the crosshair on the graph showing the location of the software's current estimation of the contact point.
  - Also note the software's first fit to your chosen elastic model, indicated by a brown line. In the example, the first fit does not match the data very well, indicating that we have more work to do before we get a better fit on our data.



6. You can change the selected curve either by going back to the Display tab or by pressing the arrow buttons under the Fit Input grouping on the Elastic Panel. You can also use the keyboard arrow keys to navigate across adjacent curves. See 8.2.3. Also, we will use the extend portion (not the retract portion) of the force curve for this model, which should automatically be selected in the drop-down box in this grouping.
7. Check your warnings: When the fit is made, the software checks for a series of common problems that are found in indentation data. The top-right corner of the panel will display any problems that have been detected with the data in yellow. Click the name of the warning for a description of the problem. If there are no issues detected, a green button labeled No Problems Detected will be displayed.
8. Input your model assumptions. Next, you will want to assume a tip shape, size Poisson ratio for your sample. These need to be entered into the model assumption grouping. The tip material is not required for soft samples, when the sample is in the kPa to a few MPa range, then the difference in a 60 GPa indenter and a 360 GPa indenter is negligible.
- a) In this case, we used an unmodified silicon nitride AFM cantilever, which we assume is a cone with a half-angle of  $36^\circ$ . For a PA gel, the literature reports various Poisson ratios- here we will use 0.4. Note that assumptions are very difficult to make, and in real-life situations it might not be feasible to expect that your assumptions are constant throughout a single data set.
9. Adjust your initial guess, if needed. The Initial Guess grouping allows you to assist the model in making first-order guesses. The Force Offset and Depth Offset are the X- and Y-locations of the contact point, relative to the 0,0 point in the force curve data. They describe the center of the crosshairs displayed on your force curve graph. When the boxes next to these parameters are checked, the program will try to solve for these offsets by iterative fitting.
10. Adjust the fit region, if needed. In many cases, only a subset of your curve will be amenable to a single model. By default, the fit is made between 10%-90% of maximum indentation. To select a smaller or larger region, press the Show Fit Region Controls Checkbox. Using




the parameters shown, you can select a lower and upper bound by inputting a specific value or just a percentage of either the maximum force or indentation that will be fit by the model. If the low parameter is grayed out, that means that you are trying to fit the force offset, and the non-contact data is required to fit for that. Note that as you change the ROI parameters, your data graph may auto-zoom the graph window.

- a) When looking at these example force curves, you can see a small kink in the extend portion of the curve at about 50% of the maximum force application. This might be due to various factors that are not described by the model. In order to exclude this section of the curve from the model analysis, we can use the ROI controls to select a sub-set of the data for analysis.

11.

### Fit your Data

- You will see the model fit represented as a brown dashed line on your force curve, and the numbers in the Results grouping will be updated. To the right we show a good fit on our data, which was obtained by rejecting all data above the kink in the force curve.



12. Look at the Results part of the panel. The modulus of the sample is highlighted in yellow. The reduced chi-square ( $\chi^2$ ) value is an estimate of the fit quality, where values closer to 1 indicate a better fit. The reduced modulus is the actual measured modulus, which is a convolution of the sample modulus and the indenter modulus.
13. 'Start Over', if needed. This button removes the region of fit restrictions and allows you to start again.
14. Select 'Entire Force Map': Selecting this tab should result in the border of the panel turning green (Figure 8.24c on page 119).
15. Fit an Entire Force Map: We started by analyzing one force curve in order to adjust our fitting parameters and assumptions. We can now analyze an entire force map.

**Note** When you do this, you apply the same assumptions and fit parameters to every curve in the force map. This might not be appropriate as the variability of the sample and the tip during the experiment might require different assumptions for different areas.

16. Check for warnings: The warnings area of the panel should now give you the option to check a small subset of your data for fit problems. It will also list the estimated length of time that such a check will take. The warnings that appear here are identical to those described in the previous section.
17. Select your data output: Using the top-left section of the panel, you can select the name of the image output. The 'Image Output: AutoName' check box is selected by default and will produce an image with the same name as the folder that contains your force map data. This name will be listed just under the check box. Be careful if you do not check the box; adding the calculated data as a layer into a new image can be tricky, especially if you are adding it to an image file that has a different pixel count or scan size.




18. Fit your data: Click the 'Fit' button. A progress bar with estimated time to completion will appear as the program fits the individual curves.

#### Explore your Data:

- A figure with results should appear. Each pixel represents a modulus value.
- To see a point in more detail, click the point on the image. The force curve display will show you that force curve and the Results area of the Elastic tab will show you the value.

19.

**Note** Very poorly fit points are displayed as dark red pixels. If you choose one of the problems noted in the Warnings section, a bright red mask indicates which areas represent the reported warnings.




20. Refit a specific point, if needed: If a single point has been fit poorly, or requires different assumptions or parameters, you can select that point and adjust the fit parameters in the Elastic tab, just as you would when analyzing single force curves. If your new parameters and result are satisfactory, you can then update the force map with the new value by clicking the 'Update Value in Map with Current Fit' button.

#### Further Analysis (Optional)

- In this example, the Force Map was taken on a homogenous gel. This was done so that many data points could be analyzed to achieve an average value for the sample. If you had a large data set (> 100 force plots), you could do a Gaussian fit to the distribution from the histogram panel. This example data set is only 16 force plot, so reading the average and standard deviation from the roughness panel is more appropriate (see [Section 7.3.1 on page 70](#)). The histogram example is described here as an example.
- To show the Histogram tab, open the Analyze Panel by clicking 'A' at the upper-right of your result map (see [Section 7.3.3 on page 76](#)).

21.




22.

**Further Analysis (Continued)**


- Click the 'Make New Histogram' button.
- The resulting histogram can be adjusted using the slider to control the number of samples per bin.
- A Gaussian fit can be applied to the data by selecting 'Do it' on the top right of the Histogram window. You can also fit a subset of the data by bracketing the region of interest with the Igor cursors (press Ctrl + i to reveal these) and then doing the fit.

**Note** In this example, the Gaussian mean value was 664.18 Pa with a standard deviation (Width) of 7.56 Pa. Compared to the actual average of 682.9 Pa and standard deviation of 8.183 Pa from the Roughness tab.


**8.10.0.2. Elastic Tab JKR Example**

This examples's instruction set uses a force map that was collected on Dow Corning Sylgard 186 Silicone Elastomer <http://www.dowcorning.com/applications/search/default.aspx?R=118EN>. The instruction set will refer to specific force curves collected on this sample. You can download this sample data set from here:


<http://www.asylumresearch.com/Files/Data/JKRData.zip>.



(a) Master Force Panel, Display Tab



(b) Single Force Plot



(c) Entire Force Map


Figure 8.25.: Master Force Panel, Elastic Tab

1. Open the AFM software. In this example, we are using Igor PRO 6.32A and the MFP3D 120804+1118. Earlier versions of the AFM software will not have the features we use in this document. Please check the Asylum support site (<https://support.asylumresearch.com/forum/content.php?4-Software>) for the latest version.
2. Load your Force Map data. From the menu bar select *AFM Analysis > kMaster Force Panel*. A new panel called the Master Force Panel should appear, along with the AR Load Path window. Using the browse button on the AR Load Path window, select the folder that contains your force map data.

3.

**Open a force curve:**

- You should now see a list containing all of your force curves towards the bottom of the Display tab of the Master Force Panel (Figure 8.25a on page 124).
- Click on the name of any one of your curves, and a graph of the force curve will appear. In this example, we've picked the force curve located at Line 1, Point 3.



-OR-

- If an image came up with the force map, you can click on the image to display the corresponding force plot.


4. Select the Elastic tab on the Master Force Panel. Then, click on the JKR tab. The panel should now look like Figure 8.25b on page 124.

**Note** By default, the Elastic tab starts in single force plot mode, as evidenced by the blue Single Force Plot subtab towards the top of this panel. This is a good place to start, even if you are analyzing an entire force map.

5.

**Inspect the force curve:**


- The force curve will be displayed with Force on the Y-axis and Indentation on the X-axis.
- Note the crosshair on the graph showing the location of the software's current estimation of the contact point.
- Also note the software's first fit to your chosen elastic model, indicated by a brown line. In the example, the first fit does not match the data very well, indicating that we have more work to do before we get a better fit on our data.



6. You can change the selected curve either by going back to the Display tab or by pressing the arrow buttons under the Fit Input grouping on the Elastic Panel. You can also use the keyboard arrow keys to navigate across adjacent curves, see Section 8.2.3 on page 97.

7.

**Select the region:** The JKR model can fit either Extend, Retract, or both (Ext+Ret). The JKR model does not include visco elastic deformation, so when there is significant separation between extend and retract (as shown in these force plots), you will not end up with a good fit to both extend and retract. The fit curve is between extend and retract, describing neither. This should be an indicator that your model may be incorrect, and you may want to do a force rate dependency study to see what effects that loading rate has on the shape of the force curves. These curves do, however, have a very clear decrease in restoring force just before the rupture on retract, this is a signature of the JKR model that is absent in Hertz and DMT. So for the sake of this example, we **set the fit input to retract** to see how well the JKR model could describe that portion of the curve. You can see it immediately does a much better job of describing that data.



8. Check your warnings: When the fit is made, the software checks for a series of common problems found in indentation data. The top-right corner of the panel will display any problems that have been detected with the data in yellow. Click the name of the warning for a description of the problem. If there are no issues detected, a green button labeled No Problems Detected will be displayed.
9. Input your model assumptions: Next, you will want to assume a tip shape, size, and a Poisson ratio for your sample. These need to be entered into the model assumption grouping. The JKR model is derived for a paraboloid shaped tip. The tip material is not required for soft samples, when the sample is in the kPa to a few MPa range, then the difference in a 60 GPa indenter and a 360 GPa indenter is negligible.
  - a) In this case, we used an unmodified silicon nitride TR800 PSA (Short)<http://www.asylumresearch.com/Probe/TR800PSA>, Olympus that we assume is a spherical apex with a radius of 20 nm. We also assume a Poisson ratio of 0.5. Note that assumptions are very difficult to make, and, in real-life situations, it might not be feasible to expect that your assumptions are constant throughout a single data set.
10. Select the fit parameters and adjust your initial guess, if needed: The more fit parameters, the more freedom the fit function has, making it easier for it to get lost. If you have a difficult time converging (getting a good fit), try turning off some parameters and adjusting them by hand. For example, adjust the Depth Offset, so that it looks correct to you and then turn off that Fit parameter. Adjust the Pull Off Force to be the minimum of the curve and turn off that parm. Adjust the Force Offset so that the base line looks correct and then turn that parameter off. Then do a Fit and see where that gets you. Once you are close, you can turn on more fit


parameters, and they should improve when you fit.

11. Adjust the fit region, if needed: In many cases, only a subset of your curve will be amenable to fitting. By default, the fit is made between 10%-90% of maximum indentation. To select a smaller or larger region, select the 'Show Fit Region Controls' check box. Using the parameters shown, you can select a lower and upper bound by inputting a specific value or just a percentage of either the maximum force or indentation that will be fit by the model. If the low parameter is grayed out, that means that you are either trying to fit the contact point (see the previous step) or you are fitting both extend and retract. In the latter case, the upper limit is found in both the extend and retract data, and the force above that level is excluded from the fit.

12.

#### Fit your data:

- You will see the model fit represented as a brown dashed line on your force curve, and the numbers in the results grouping will be updated. To the right we show a good fit on our data.




13. Look at the results part of the panel. The modulus of the sample is highlighted in yellow. The reduced chi-square ( $\chi^2$ ) value is an estimate of the fit quality, where values closer to 1 indicate a better fit. The reduced modulus is the actual measured modulus, which is a convolution of the sample modulus and the indenter modulus.
14. Start Over, if needed: This button removes the region of fit restrictions and allows you to start again.
15. Select Entire Force Map: Selecting this tab will result in the border of the panel turning green (Figure 8.25c on page 124).
16. Fit an Entire Force Map: We started by analyzing one force curve in order to adjust our fitting parameters and assumptions. We can now analyze an entire force map. Note that when you do this, you apply the same assumptions and fit parameters to every curve in the force map. This might not be appropriate as the variability of the sample and the tip during the experiment might require different assumptions for different areas.
17. Check for warnings: The warnings area of the panel should now give you the option to check a small subset of your data for fit problems. It will also list the estimated length of time that such a check will take. The warnings that appear here are identical to those described in the previous section.
18. Select your data output: In the upper-left section of the panel, you can select the name of the image output. The AutoName checkbox is selected by default and will produce an image with the same name as the folder that contains your force map data. This name will be listed just under the check box. Be careful if you do not check the box; adding the calculated data as a layer into a new image can be tricky, especially if you are adding it to an image file that has a different pixel count or scan size.
19. Fit your data: Click the Fit button. A progress bar with an estimated time to completion will appear as the program fits the individual curves.

20.

**Explore your data:**

- A figure with results should appear. Each pixel represents a modulus value.
- To see a point in more detail, click the point on the image. The force curve display will show you that force curve and the Results area of the Elastic tab will show you the value.

**Note** Very poorly fit points are displayed as dark red pixels. If you choose one of the problems noted in the Warnings panel, a bright red mask indicates which areas represent the reported warnings.




21. Refit a specific point, if needed: If a single point has been fit poorly, or requires different assumptions or parameters, you can select that point and adjust the fit parameters in the Elastic tab, just as you would when analyzing single force curves. If your new parameters and result are satisfactory, you can then update the force map with the new value by selecting the 'Update Value in Map with Current Fit' button.

22.

**Further analysis (optional):**

- It may be possible to distinguish between two different samples by looking at the histogram of the obtained Young's modulus.
- To show the Histogram tab, make the Analyze Panel by clicking 'A' in the upper-right section of your result map, see [Section 7.3.3](#) on page 76.




23.

**Further analysis, continued:**

- Click 'Make New Histogram' button.
- The resulting histogram can be adjusted using the slider to control the number of samples per bin.
- A Gaussian fit can be applied to the data by clicking 'Do It' in the upper right section of the Histogram window. You can also fit a subset of the data by bracketing the region on interest with the Igor cursors (click Ctrl + i to reveal these) and then doing the fit.

**Note** In this example, the Gaussian mean value was 1.25 MPa with a standard deviation (Width) of 0.0287 MPa.

**8.11. Online Video**

Consider watching this video: [Hertz Fitting for Nanomechanical Analysis](#) (requires an internet connection).



## 9. Particle Analysis

CHAPTER REV. 2049, DATED 09/27/2018, 11:13.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

Currently the Video Overlay software feature for AR SPM Software Version 14 is described in this online video tutorial: [Particle Analysis](#).

## 10. Video Overlay

CHAPTER REV. 1984, DATED 11/10/2017, 14:24.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

Currently the Video Overlay software feature for AR SPM Software Version 14 is described in this online video tutorial: [Video Overlay - Combining Optical and AFM data](#).

# 11. ARgyle 3D Data Visualization

CHAPTER REV. 2436, DATED 09/04/2021, 14:34.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

## Chapter Contents

11.1	Viewing Saved Images in 3D . . . . .	132
11.2	Overlaying Data on a Surface . . . . .	133
11.2.1	Basic Examples of Color Overlay . . . . .	134
11.2.2	Overlaying Data from Difference Sources . . . . .	136
11.2.2.1	Overlay of Height on Height . . . . .	137
11.2.2.2	Overlay tricks . . . . .	138
11.2.2.3	Flattened Surface Overlay on Topography . . . . .	138
11.2.2.4	Error Signal Overlay on Topography . . . . .	138
11.3	Display Settings . . . . .	140
11.3.1	Axes, Prefs, Lights, and Windows Tabs . . . . .	140
11.3.2	Color Tables, Rotate, Zoom, Pan, and Lighting . . . . .	141
11.3.3	Z Axis Range, Offset, and Clipping . . . . .	144
11.3.4	3D Graphics Aspect Ratio . . . . .	146
11.3.5	Anaglyphs, Auto Rotation, and Fly-Over. . . . .	147
11.4	Exporting Images . . . . .	148
11.5	Activating Real Time ARgyle™ . . . . .	149
11.6	Advanced Command Line Control . . . . .	150
11.7	Argyle Lite Windows App . . . . .	150

ARgyle™ (a play on ArGL or Asylum Research Graphics Library) is the MFP-3D's 3D image rendering software. It uses open GL code and the graphics card on your PC to produce stunning 3D representations of your data. In real time or with offline analysis, custom colors and specular lighting can bring out features in your data that are not easily perceived in the standard 2D view. It also allows multiple data channels to be overlaid as color on another data channel.

There are two ways to get started using ARgyle™, depending on whether you are imaging in real time or doing offline analysis. 3D representations of saved data will be covered first, so this tutorial can be used even if you are not actively imaging a sample. For the same capabilities with reference to data being collected in real time, please see [Section 11.5 on page 149](#).


## 11.1. Viewing Saved Images in 3D

1. Open the Image Browser Panel as outlined in [Section 7.1 on page 32](#).

2.

**Open a 3D Window:**

- Open an image of choice. This will bring up the 2D image display window.
- Select a tab of the image: probably Height or, as shown here, Zsensor.
- Click the '3D' button.



3. A 3D rendered image will appear in a new window. Figure 11.1 on page 133 shows a typical example. The data is scaled in a box of aspect ratio 3 to 1. The color bar corresponds to the total range of the Z axis.
4. Try a few things with your mouse:
  - a) Click-drag on the 3D surface with the left mouse button to rotate the view.
  - b) Use the mouse wheel to change the level of zoom. Clicking and dragging the left mouse button while holding down the 'Ctrl' key has the same effect.
  - c) Right-click-drag to change the angle of the lighting.
  - d) 'Shift' + left-click-drag to pan the image within the XY plane.




Figure 11.1.: Typical ARgyle window when it first opens.


## 11.2. Overlaying Data on a Surface

Figure 11.2 on page 134 shows an example of an overlay. Unlike the images covered so far, the color map overlay has no relation to the data channel it is overlaid on. In this case, the color is proportional to surface potential at the cantilever tip as measured by Kelvin Probe Force Microscopy, but the shading and 3D rendering is proportional to the sample topography. The potential as measured at each pixel is painted over the sample topography. It is not unlike a 3D

rendering of a geographical map in which colors represent vegetation type. This effect can be quite useful for depicting many quantities measured by AFM, such as magnetic field, conductivity, tip sample capacitance, and cantilever phase.

**Note**

AFM images commonly have spikes in the data, which can be very distracting in surface plots. Use of the filter panel can be a quick and easy way to clean up the image for use as a surface plot. See [Section 7.2.5 on page 63](#) for detailed instructions.




**Figure 11.2.:** Surface potential overlay on topography.

### 11.2.1. Basic Examples of Color Overlay

1. Open the Image Browser Panel as outlined in [Section 7.1 on page 32](#).

2.

**Open a 2D Window:**

- Open an image of choice, which displays the 2D Image Display window. Here we show a polymer topography image and the phase channel. Both are flattened and slightly filtered for noise.


3.

**Set the Color Channel:**



- On the 2D image graph:
  - Select the channel you want for the color. This will usually be phase, current, or some other measured quantity.
  - ‘Ctrl’ + Click the ‘3D’ button. This will set the currently displayed layer as the color for the next 3D image.

-OR-

- Open the Master ArGL Panel by selecting *AFM Analysis > 3D Surface Plots* from the menu bar:
  - Click the New tab.
  - Click the ‘Color’ checkbox.
  - From the ‘Color’ menu, select the same data file as you selected for the surface.
  - From the ‘Data Type’ menu, select the desired channel in that image.



4.




#### Set the Surface Data:

- On the 2D image graph:
  - Select the channel you want for the surface. This is usually Height, Zsensor, or a modified version thereof.
- OR- Click the '3D' button.
- From the Master ARgyle panel:
  - Click the New tab.
  - Select the topography data file from the 'Surface' menu.
  - Select the desired channel of that image from the 'Data Type' menu.
  - Click the 'Do It' button to create the 3D view.
- The image shows a nice example of how one value of phase associates with ridges and another value with furrows.

**Note** Once the 3D image is rendered, neither the surface nor the color channel can be reset without closing the window.

### 11.2.2. Overlaying Data from Difference Sources



**Figure 11.3.:** Force Map overlay on LFM image.

The previous section overlaid the phase channel from one image onto the height channel from the same image. As you may have guessed from the dropdown menus when you selected the Surface and Color sources, the channels do not have to be from the same image; surface and color sources can be selected from any open images. Figure 11.3 on page 136 shows an example of a lateral force image which is used as the surface source (topography actually represents force) with a color overlay based on a map of force curves.

The sample is a micro contact printed mercapto undecanoic acid (-COOH terminus), backfilled

with dodecanethiol (-CH<sub>3</sub> terminus). The Au coated SiN<sub>x</sub> tip had the acid thiol on it as well, and it was acquired in a pH 4 buffer standard. The example was done in moderate haste to show this type of experiment can be done, and seems to correlate reasonably well.

### 11.2.2.1. Overlay of Height on Height

An overlay of topography as a color on topography in a 3D image can be very useful when trying to create a certain data landscape. In a standard 3D image, the color bar of the height is directly linked to the range of the Z axis. The only way to get a certain color to occupy a certain altitude in the topography is to translate offset the data, usually causing some of it to go outside of the frame box.


By associating color with a secondary channel, the ArGL Master Panel allows the color scale to be offset freely. Please follow this brief demonstration.

1. Open the Image Browser Panel as outlined in [Section 7.1](#) on page 32.

2.

#### Open a 2D Window:


- Open an image of choice which brings up the 2D Image Display Window.
- Select a tab of the image data channel of interest, probably Height or, as shown here, Zsensor.
- 'Ctrl' + click the '3D' button.



3.



#### Open a 3D Window, Topography:

- Go back to the image.
- Click on the '3D' button, but do not hold down 'Ctrl'.





4.



#### Manipulate the Color Bar Scaling and Offset:

- Select the View tab.
- Set 'View Mode' to "Color Scale".
- Use the sliders to change the color scale range and offset. Notice that the total range on the color scale is 2 microns, while the total range on the Z axis is 500 nm. This cannot be achieved in a non-overlay type graph.

#### 11.2.2.2. Overlay tricks

These are not intended for quantitative scientific visualization but can be very helpful when trying to communicate features of an image.

#### 11.2.2.3. Flattened Surface Overlay on Topography

You may want to do a 3D rendering of some material that has smaller features but some larger curvature of radius, for instance, large spheres or hair. The problem is that the finer details of the topography features do not get colored well in ARgyle due to the extreme changes in Z scale of the material.<sup>1</sup>


To solve this problem, flatten a height image channel with a second order flatten to take out the large-scale curvature of the material, then overlay it onto a Z sensor or copied height channel that has nothing greater than a 0-order flatten applied to it.

Using ARgyle, overlay the 2nd order flattened image onto the other, such that the 2nd order image is the colored layer. This will help to bring out the smaller features much better. [Figure 11.4 on page 139](#) shows an example of this technique. Note that this technique is mostly to enhance visualization. The color bar may best be left off since its scale will only reflect, at best, local variations with respect to the surface of the object.

#### 11.2.2.4. Error Signal Overlay on Topography

The amplitude signal in AC mode imaging or the deflection channel in contact mode imaging, also known as the "error signal", is a measure of the inability of the Z feedback mechanism to actually


<sup>1</sup> This idea courtesy of Scott Maclaren, Center for Microanalysis of Materials, Frederick Seitz Materials Research Laboratory, University of Illinois, Urbana Champaign.



**Figure 11.4.:** Overlay of a second order flattened height image on the topography of the same image.

keep the cantilever amplitude fixed. When the cantilever encounters edges in the sample, it will take some time for the feedback to react. During this time, the amplitude or deflection will not be at their respective set points. The resulting image of this signal often looks appealing since it accentuates edges and ignores slower variations in topography. It contains no direct topographical information, only convoluted information which, in a way pleasing to the human eye, brings out small surface detail. It is not unlike some of the traditional edge finding filters that can be found in the Filter subtab of the Modify tab.

When imaging cells, users commonly overlay the amplitude signal onto the overall topography to visually enhance the presence of small surface features. Always clearly state how the image was processed, since this borders on “Photoshopping” the actual data.



**Figure 11.5.:** Fibroblast cells on a glass substrate. Overlay of the amplitude image channel on the topography. Some blue colored specular lighting added.

## 11.3. Display Settings


When the '3D' button is clicked, the ArGL Master Panel is opened along with the 3D representation of the data. This gives additional controls over the 3D image beyond the mouse clicks and drags of the last section. Each of the tabs will be discussed but out of order for the sake of clarity.

### Note


On the first activation of a 3D image, the panel will have certain default parameters (i.e., Zoom, scale, rotation, color, etc.) that may appear off-scale (in Z) if the surface features are large relative to the XY scan area, or close to scale (in Z) if the features are small. If you have already processed a 3D window in the current experiment, the new ARgyle image will have the same rotation, zoom, etc. as the previously opened 3D window. It will have the same Z data scale as the 2D image in the Display Window.

### 11.3.1. Axes, Prefs, Lights, and Windows Tabs


The tabs that are self-explanatory will be discussed first. Take a look at [Figure 11.6 on page 140](#) and peruse the settings. Within Igor, you can always click the '?' buttons next to any item for more explanation.




(a) Axes tab



(b) Lights tab



(c) Prefs tab



(d) Windows tab

**Figure 11.6.:** Master ArGL Panel Tabs

Here are some of the more popular actions on these panels, even though many users will rarely use any of them:

**Remove text or axes** Go to the Axes tab. Uncheck the 'Show Axes' check box or the 'Axes Labels' check box.

**Arrange or close multiple 3D image windows** Go to the *Windows* tab. Hold 'Shift' and click on windows in the list to stack, tile, or close them en masse.

**Make surfaces appear shiny** Go to the Lights tab and check 'Shiny'. The percentage next to the 'Shiny' checkbox will adjust the intensity of specularly reflected light. Apply changes to the last viewed image or check 'All' to apply to all current and future 3D windows.

**Change illumination color** Go to the Lights tab and change the light color using the 'Color' drop-box. This only affects the specularly reflected light, so if 'Shiny' is set to 0%, then the light color is irrelevant.

**Change the background color** Go to the Prefs tab. Use the 'Color' dropdown to change the background.

**Note** Apply changes to the last viewed image or check 'All' to apply the changes to all current and future 3D windows.


### 11.3.2. Color Tables, Rotate, Zoom, Pan, and Lighting

This instruction set will show the relevant panel, usually the View tab, with different modes selected on the left, the resulting image on the right, and instructions below.



1.

#### Starting Point:

- We will apply a series of actions to the image shown at right.



2.





#### Choose a Different Color Table:

- Select the Display tab.
- Select a different color table using the 'Color Map' menu.
- Invert the direction of the color table with the 'Invert' check box.

#### Optional (not shown in the above image):

- Check the 'Negative Color' check box to make the image be a color negative of the initial color scale. Red, for instance, would be displayed as green.
- Remove the color scale bar from the right of the image by unchecking 'Color Bar'.

3.




### Rotate the Image:

- Under the View tab, select “Rotate” mode from the ‘View Mode’ menu.
- Change the view angle using the text fields or the sliders.
- Notice that the colors associated with the pitch and rotation fields are also applied to the sliders.

### Optional

- Or, with the mouse, left click-drag on the 3D image to change the rotation variables.

4.





### Zoom In:

- Select the View tab.
- Set ‘View Mode’ to “Zoom”.
- Enter a number for ‘Zoom Factor’ or slide the left slider back and forth.

### Optional

- Or, with the 3D image window selected, rotate the mouse wheel to zoom in and out.

5.



**Translate the Image:**

- Select the View tab.
- Set 'View Mode' to "Offset".
- Enter XY offset values or move the sliders to slide the image around on the screen.

**Optional**

- Or, left Shift + click to translate the 3D image.

6.

**Change the Lighting Angle:**

- Select the View tab.
- Set 'View Mode' to "Light Rotate".
- Enter values or use the sliders to change the incident angle of the image lighting.


**Optional**

- Or, right-click-drag the mouse pointer over the 3D image.



## 11.3.3. Z Axis Range, Offset, and Clipping

1.

**Change the Z Axis Range:**



- Select the View tab.
- Set 'View Mode' to "Data Scale".
- 'Clip' must be set to "None".
- Enter a new value into the 'Data Scale' field. In our example, we cut the data scale in half, increasing the size of the 3D topography by a factor of 2.
- Alternatively, you can use the yellow left-slider.

**Note** Since the color scale is normally fixed to the vertical range, the colors of the surface that extend above the box stay constant. The only reason the surface still looks like data is due to the shadows cast by the lighting. If you go to the Master ArGL Panel and turn off the lights under the Lights tab, then the out of range parts of the image will become featureless.

**Note** The data expands equally up and down from the vertical midpoint of the Z range.

**Note** The extremes of the data scale can give visually confusing results or may even seem to erase the data. If the data seems irrevocably distorted, the 'Fix' button will restore the default settings.



2.



### Clip Data Extending Outside the Frame:

- Select the View tab.
- Set 'View Mode' to "Data Scale".
- Set the 'Clip' to "Clamp". The data extending beyond the vertical range are capped with solid surfaces.
- To leave an open hole in which you can see the backside of the data, select the "Clip" option from the 'Clip' dropdown menu.

3.




### Change the Z Axis Offset:

- Select the View tab.
- Set 'View Mode' to "Data Scale".
- The 'Clip' can still be set to "Clamp".
- Enter a new value into the 'Data Offset' field. In our example, we offset by 125 nm, which brought the taller features nearly back inside the box and moved lower features outside the box.
- Alternatively, use the blue right-slider to adjust the offset value.





## 11.3.4. 3D Graphics Aspect Ratio

1.

**Change the 3D Graph Aspect Ratio:**

- Select the View tab.
- Set 'View Mode' to "Zoom".
- Change the 'Z Scale' from the default value of "0.33" to "0.66". The result is a stretched graph which resembles the one in [Step 1 on page 144](#). In this case, however, the entire box is taller, and the data still fits inside that box. The data will not be clipped and will not exceed the color scale.
- Alternatively, you may move the rightmost, blue-slider to change the value.

2.



### 3D Graph Aspect Ratio 1:

- Select the View tab.
- Set 'View Mode' to "Zoom".
- Click the 'Aspect 1' checkbox.
- The Z Scale value will assume a value which shows the sample in its true proportions. For many AFM samples that have tens of nanometers of topography over a many micron scan range, it will appear almost flat. In our case, the sample is fairly rough and still gives something to look at.
- Changing the 'Zscale' by entering a number or using the slider will automatically uncheck 'Aspect 1'.

**Note** Aspect 1 only works when the Z axis has units of length. While not common, any quantity, such as phase or amplitude, can be plotted on the Z axis. In this case, the check box has no effect.


**Note** The term *Aspect Ratio 1* refers to the actual sample X, Y, and Z dimensions, not to the way the frame looks on the screen. This is not to be confused with setting the Zscale to 1, which always makes the 3D frame appear to be a cube.

### 11.3.5. Anaglyphs, Auto Rotation, and Fly-Overs.

#### Viewing with 3D Glasses:

- Go to the Display tab.
- Click the 'Anaglyph' button.


- Note** While guaranteed to work on all computers shipped with Asylum AFMs, the Anaglyphs may not always work with unusual graphics processors. This is often the case with laptop computers.



**Note** 3D viewing is most effective with a gray scale color table, as shown in this example image.

2. **Auto Rotation (Spinning 3D Image):**
- Select the View tab.
  - Set 'View Mode' to "Spin".
  - Type in some small values or use the sliders.
  - The 3D image will slowly spin until the spin rates are set to zero again.

**Note** This can be combined with the 3D glasses though may cause motion sickness.



## 11.4. Exporting Images

There are many ways to export images from the ARgyle window, including the following four methods.

1. **Cut and Paste a Screen Shot:**
- Make the 3D image window the forwardmost image.
  - From the menu bar select *Edit > Copy*. This puts a screenshot of that window on the clipboard.
  - Open another program where you can paste this bitmap, such as Paint, Photoshop, or e-mail.
  - From the menu bar of that program, choose *Edit > Paste*.

**Note** Since this is a screenshot, the larger you make the window on screen, the higher your exported image resolution.

2. **Save a Screen Shot to Disk:**
- Make the 3D image window the forwardmost image.
  - From the menu bar select *Edit > Export Graphics*.
  - Select the file name, location, and type.
  - Click 'Save' to export the screenshot to disk.

**Note** Since this is a screenshot, the larger you make the window on screen, the higher your exported image resolution.

3. **Export to Layout:**
- Make the 3D image window the forwardmost image.
  - Select the New tab of the Master ArGL panel.
  - Click the 'Export to Layout' button. To learn more about this, see [Section 7.4.4 on page 91](#).

**Note** This will work only for open 3D windows. All the dropdown menus above the 'Export to Layout' button have no bearing on this feature.

**Save a High Resolution Bitmap to Disk:**

- Make the 3D image window the forwardmost image.
- Go to the Igor Command line (Ctrl+J).
- Type:

```
argl_export2("", "", 0, 1, 3000, 2000)
```

4.

followed by the Return/Enter key for an exported image of 3000 pixels wide by 2000 pixels tall. The pixel values have a 4000 maximum.

- From the dialog box that appears, select the file name, location, and type. In this case, the type should be a bitmap (BMP).
- Click 'Save' to export the graphic.

**Note** If you have the color bar showing and choose a narrow window, the color bar will fall atop your 3D graphic. Choose an aspect ratio similar to the screen. The Command line approach will only show which 3D window is currently on top, not how it is proportioned.

## 11.5. Activating Real Time ARgyle™


The lighting and shading features of ARgyle can be useful when fine tuning imaging parameters. It is possible to update a live 3D image, line by line, as data is being collected. ARgyle does not distinguish between real-time images and saved images; all of the manipulations discussed in this chapter can be applied to either.

1. Start imaging.

**Open a 3D Window, Topography:**

2.

- Open the Master ArGL Panel by selecting from the menu bar *AFM Analysis > 3D Surface Plots*. Select the New tab.
- From the 'Surface' dropdown menu, select 'Realtime'.
- From the *Data Type* pull-down menu to the right, select the desired channel, usually 'Height' or 'Zsensor'.
- Either skip to the next step or, IF you want color overlays, click the 'Color' checkbox and select *Realtime* as the source. For *Data Type*, the usual channel is *Phase*, but any other channel of your choice may be selected.
- Click the 'Do It' button.



You will now have a 3D image on screen that fills up line by line. You can apply any manipulations described in this chapter.

## 11.6. Advanced Command Line Control

ARgyle has advanced command line features that will eventually be incorporated in the Igor GUI interface. If you are curious, please read the following online help file. From within the AR SPM software, select from the menu bar:

*Help > AR Help Files > Xop Help Files > ARgyle help*

## 11.7. ARgyle Lite: Stand Alone Data Visualizer for Windows

Asylum Research has released a free program which incorporates all the functionality of ARgyle 3D image viewing. See [Figure 11.7 on page 151](#) for a screenshot of the application. It can be downloaded from <https://support.asylumresearch.com/forum/content.php?157-ARgyle-Light>

This program is an excellent way to share and view Asylum Research SPM data on computers that do not have the Igor Pro / AR SPM software installed. You can send \*.ibw files to anyone if they install the free viewer.

For true 3D visualization, please contact Asylum Research for a free pair of 3D glasses: [sales@AsylumResearch.com](mailto:sales@AsylumResearch.com).

If you understand the ARgyle functionality described in this chapter, you will find the ARgyle Lite user interface quite familiar. The only limitation is that the software works with only one image (.ibw) file at a time, so it is not possible to overlay data from one image onto another, as in the force map example in [Figure 11.3 on page 136](#).




Figure 11.7.: ARgyle Lite windows application.

# 12. Managing Image Data

CHAPTER REV. 2433, DATED 08/28/2021, 18:39.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

## Chapter Contents

12.1	Introduction . . . . .	152
12.2	Data Channels . . . . .	153
12.3	Selecting Image Pixel Density . . . . .	153
12.3.1	Choosing Enough Pixels . . . . .	153
12.3.2	Aspect Ratio . . . . .	154
12.3.2.1	Pixel Aspect Ratio . . . . .	154
12.3.2.2	Image Aspect Ratio . . . . .	154
12.3.3	Notes on Selecting Pixel Densities . . . . .	154
12.4	Master Channel Panel . . . . .	154
12.4.1	Opening Image Channels . . . . .	154
12.4.2	Adjusting Image Z Scale . . . . .	156
12.4.3	Adjusting Trace/Retrace Scale . . . . .	156
12.4.4	Real Time Flattening . . . . .	157
12.4.5	Save Planefit . . . . .	158
12.4.6	Real-time or Saved Displays . . . . .	158
12.5	Capturing Images . . . . .	158
12.6	Offset Scan Area . . . . .	160
12.7	Zooming Images . . . . .	161
12.8	Multiple Offset/Zoom Parameters . . . . .	161
12.9	Real Time Section Analysis: . . . . .	164
12.9.1	Basic Real-Time Line Section Analysis: . . . . .	164
12.9.2	Adjusting the Angle of the Section Line . . . . .	166
12.9.3	Plotting Multiple Section Lines on One Section Plot . . . . .	168
12.9.4	Export the line section XY data: . . . . .	168
12.10	Other Features . . . . .	169
12.10.1	Windows Manager . . . . .	169
12.10.2	Saving Igor Experiments . . . . .	171
12.10.3	Adding Bars to the Sum and Deflection Meter . . . . .	171
12.10.4	Modifying Igor Graphs . . . . .	171

## 12.1. Introduction

This chapter assumes you are imaging, but sooner or later you'll need to save the data and quickly analyze it a bit while it's still being collected. Hence, this chapter covers topics such as the following:

- How to select the proper incoming data streams for storing in memory.
- How to make sure that you save it to disk.
- How to best visualize the incoming data, including zooming, sectioning, and adjusting colors and contrast.
- How to analyze it on the fly, even before an image is completed and stored.

## 12.2. Data Channels

An SPM image usually has many layers of data. Typically, the sample topography will be recorded, but other concurrent measurements might be cantilever amplitude and phase (in AC mode), cantilever deflection (in contact mode), cantilever current (for conductive AFM), or a host of other signals. These concurrently measured quantities are known as channels, or image channels. For real-time data, each channel can have a real-time image window. For data retrieved from the hard drive, each channel creates a tab in the image display window. We will often refer to the channels as layers also.

## 12.3. Selecting Image Pixel Density

No matter how many layers or data channels an image has, each layer will have the same physical dimension and the same number of pixels. Before we launch into the business of setting up the channels, we'll discuss how to select the pixel density.

The Scan Line and Scan Points (i.e., pixel density) are selected in the Main tab of the Master Panel using the respective variable fields. Note that the tip must be withdrawn to change these values, unless the Delay Update checkbox is activated.



Figure 12.1.: Master Panel Main Scan Lines And Points

### 12.3.1. Choosing Enough Pixels

One important calculation to keep in mind:

$$\text{PixelSize} = (\text{ScanArea} / \# \text{ of pixels})$$

For example, for 512 x 512 pixels in a 5µm scan area, each pixel is 9.76nm<sup>2</sup>; while in a 90µm scan area, each pixel covers 175.78nm<sup>2</sup>.

So make sure the feature of interest has a high enough resolution such that it isn't pixilated, but isn't so high that it doesn't take 2 days to collect the image.

The absolute number of simultaneous data channels at high pixel resolution is difficult to quantify but is directly related to computing power. For example, twelve channels at 2k x 2k can be collected



with the stock PC, or three channels at 4k x 4k. At lower pixel resolution (512 x 512), up to twenty images can be simultaneously collected. Obviously the higher the pixel density, the longer it takes to collect that data.

### 12.3.2. Aspect Ratio

#### 12.3.2.1. Pixel Aspect Ratio

If you are set up to collect an image of 10 $\mu$ m by 10 $\mu$ m and enter in 512 pixels by 128 lines, your image will still be square, but your pixels will be vertically elongated.

#### 12.3.2.2. Image Aspect Ratio

The aspect ratio of the image can also be changed in factors of 2, whole integer numbers. This will automatically update the number of scan lines or points. For example, say a 1:1 image at 512 x 512 was being collected and then changed to 1:4; the scan lines will be changed to 128, while the points will remain at 512.

### 12.3.3. Notes on Selecting Pixel Densities

When the 'Delay Update' check box is activated, you can change the pixel density and have it take place when the scan reaches top or bottom or when the 'Frame Up' or 'Frame Down' button is clicked.

If the 'Delay Update' check box is not activated, the tip must be withdrawn to change these values.

Since Igor is a memory-based program, collecting images with very large pixel densities (i.e., 2 to 4k x 2 to 4k) is most likely to occur with a fresh experiment/template, as opposed to one that has a lot of information already stored in it. To help accommodate this, saving images to memory is not recommended.


You may use a low pixel density to scan quickly to find an area of interest, and then increase the density to acquire a high-quality image.

## 12.4. Master Channel Panel

This panel is where to select data channels to display and capture, adjust Z scales, color tables, and flatten.

### 12.4.1. Opening Image Channels

To open and activate additional image channels before the tip is scanning, go to the Master Channel Panel. The default ARSPM software is set for AC mode with Height, Amplitude, and Phase channels (Figure 12.1 on page 153) to be open in the Master Channel Panel, but it is advised to always also capture the Z sensors data, especially when working with topography features larger than 500 nm (12.1). This is due to the fact that the closed loop sensors are very linear and quiet over their range, whereas the piezo voltage (i.e., height channel) suffers from piezo creep and hysteresis.



(a) Channel Panel Height Tab

(b) Channel Panel 4 Tab

**Figure 12.2.:** Master Channel Panel. Select addition channels from the Input dropdown menu in the Master Panel.

Whatever the first channel tab is selected as, it will be 32 bit; all others are 16 bit.

Additional image channels cannot be activated during scanning.

Dropdown menus in the Master Channel Panel are self-explanatory, and the Help menus (question marks) to the right further explain the features, including Color Map, Live Display (12.4.3), and Live Flatten (12.4.4), which you can select as needed.

It is best to save both Trace and Retrace images, in case there is sample slope or some debris stuck to one side of the tip that will give an imaging artifact. Capturing the trace and retrace data usually increases your chances of obtaining a publication quality image.

**Input** Choose the data channel to display, per respective tab; Note: only activated when tip is not scanning.

**Data Scale** Adjusts the Z scale of the image; or click the 'Fix' button to autoscale (see 12.4.2).

**Auto** Autoscales the Z scale data during imaging (see Section 12.4.3 on page 156).

**Data Offset** Adjusts where the center of the (Z) Data scale of the image, or click the 'Fix' button to autoscale.

**ColorMap** Offers many color tables.

**Live Flatten** (see Section 12.4.4 on page 157)

**Save PlaneFit** (see Section 12.4.5 on page 158)

**Capture What** (see Section 12.4.6 on page 158)

**Live Display** (see Section 12.4.6 on page 158)

**Show Scope** (see Section 12.4.3 on page 156)

### 12.4.2. Adjusting Image Z Scale

There are many ways to adjust the Z scale in the image channel, including:

1. Select the 'Auto' check box in upper-left corner of image channel tabs in the Master Channel Panel.
2. Click the 'Fix' button next to Data Scale.
3. Draw a quadrangle onto the image window, right- (or left-) mouse-click inside the selected area and select 'Fix Scale'. This fixes the scale according to the Z values within the selected area.

**Note** Alternatively, you can select 'Fix All Scales' which rescales the Z data in all imaging data channels. This is a very nice feature.

4. Manually type a 'Data Scale' value, or use arrows to select a value.

### 12.4.3. Adjusting Trace/Retrace Scale

The individual Trace and Retrace images can be viewed and saved by selecting them in the Master Channel Panel. To view Trace or Retrace, confirm that the proper check box is selected.


The red bar/cursor to the left of the image indicates which scan line the tip is located (i.e., which slow scan line it's rastering).

Trace or Retrace oscilloscope scan lines below the image can be turned ON or OFF with the 'Show Scope' check boxes at the bottom of the Master Channel Panel. Trace is Red; Retrace is Blue.

The 'Auto' check box, when selected, autoscales the trace/retrace lines (below the image) with every slow scan line added. Typically, it shows the slope of the sample relative to the tip, because it is auto-scaled. Even if the slope (Z) is sub-nanometer over several tens of microns (XY), it is still going to display slope because its auto-scaled (Figure 12.3).


If the 'Auto' checkbox is unchecked, it scales to the Data Scale value. For example, if the Data Scale value is too low, then trace/ retrace will be clipped (shown below); this also depends on the Data Offset value.

Figure 12.3 shows an example result in the trace/retrace lines from changing the Data Scale value. Notice that the larger peaks are clipped because the Z scale is lower (red ellipse). By increasing the Z scale, the entire traces can be scaled down (Figure 12.3).



**Figure 12.3.:** Effects of adjusting the Data Scale: A) 'Auto' check box activated; B) 'Auto' check box unchecked, smaller defined Data Scale (notice larger features clipped in the Trace/Retrace lines); C) 'Auto' check box deactivated, larger defined Data Scale.


- The Data and Offset Scales of the Trace and Retrace image channels can be separately controlled by clicking the 'Setup' button in the Master Channel Panel (Figure 12.4). The upper set of scale value applies to the Trace image, while the lower applies to the Retrace. From the author's understanding, if just the top is open, it applies to both Trace and Retrace images of the respective channel.
- It is a good habit to always save the data as raw, unflattened data. To do this, select *None* in the 'Save PlaneFit' dropdown menu (see Section 12.4.4 on page 157).



**Figure 12.4.:** Activating Data Scale and Offset setvar windows for Trace and Retrace images. Green arrow shows where to choose live flattening.

#### 12.4.4. Real Time Flattening

The Live Flatten menu allows the real-time scan to be flattened via plane, line, Masked Line, or not at all. This modification is only applied to real-time data, but do not apply to saved data.



**Figure 12.5.:** Live Flatten Choices

**None** Leaves you with your raw data.

**Offset** Takes the average value of the scan line and subtracts that from the data.

**Line** Fits each scan line to a straight line, subtracts this line out of the scan line.

**Masked Line** Fits a line to each scan line but doesn't use any data on the line that is more than 1/4 of the Data Scale away from the line. If you change the Data Scale, the lines after the change are affected but not the lines already done. If you do want to redo the current lines, reselect Masked Line in the popup.

### 12.4.5. Save Planefit




Figure 12.6.: Planefit Choices

**Caution** This will modify the Saved data!

**None** Saves the raw data.

**Offset** Takes the average value of the Image and subtracts that from the data.

**Planefit** Subtracts a first order XY Plane from the Image.

The rest of these options modify your data quite a bit more, but you may return the data to its raw state with the 'Ultra Restore Layer' button on the Modify panel.

**Flatten 0** Removes the offset from each individual line.

**Flatten 1** Removes the offset and slope from each individual line.

**Histo Flatten** Removes the offset from each line by looking at a histogram of each line to determine the offset. The slope is removed in the X direction by a plane fit on the whole image.

**Magic Mask** Performs a first order flatten on each line, then calculates a mask and redoes the flatten. This process is reiterated until it is satisfied. This process is similar to flattening offline images


**Magic Mask (Pits)** Use if the sample happens to be imaging something with pits. This algorithm can't tell if the features are pits or bumps.

### 12.4.6. Real-time or Saved Displays

Use the dropdown menus to display or save the trace, retrace, or both images. The author suggests saving both traces, though rarely displays both, unless the sample is expected to give an interesting result on the trace and retrace.

## 12.5. Capturing Images

The process of capturing image files occurs in the Main tab of the Master Panel.



**Figure 12.7.:** MasterPanelMainSaveOptions

**Base Name** Type a 17-character (or less) filename. The filename can't start with a number.

**Base Suffix** Increases this by one with every additional captured image, as the software is designed to do continuous capture, increasing.

**Note** Allows an unlimited length of text to describe imaging and experimental conditions that can later be recalled in the Show Note option in the Commands dropdown menu of the Display window.

**Save Images** When selected, a dialogue appears to choose where to save the image data. Use the 'Browse' button to change folders or type a new folder name. Alternatively, the path can be typed using colons ( :) or backslashes ( \ ) to separate folders.

Note that, during imaging, this button changes to 'Save Partial', a function that allows you to save a partial image before the entire image has been collected.


**Save Status** Shows if it is saving the current scan or the next.

**Save Prev** Saves the previous image scan if you had not initially saved it to disk. (It is temporarily stored in experiment memory.)

**Browse** Allows you to look at stored images.

### Save Graphics File (standard Igor panel):

If you only want a screenshot of one of the windows, go to *File > Save Graphics*, and the Save Graphics File dialog box appears (Figure 12.8).



**Figure 12.8.:** Saving graphics in Igor allows you to define the size, file type, screen resolution, file-name, and path.

There are many things that you can define yourself in this panel, including:

**Size** The size of the saved graphics image can be adjusted here.


**File Format** Choose the type of file (TIF, JPEG, etc.).

**Resolution** Choose 1, 2 or 4x the resolution- larger values good for larger presentations to reduce appearance of pixilation.

**File Name** Enter a file name.

**Path** The path needs to be defined. You can do this by going to the main menu and choosing:

- *Misc > New Path*. This will bring up the New Symbolic Path dialog box (Figure 12.9).
- Name the path. (The author generally uses one letter, because it's easy.)
- If you need to change the path to the folder, click the 'Path' button. This brings up the Pick a Folder dialogue.
  - Click the 'Browse' button to select the folder you want and click 'OK' when complete. This will update the New Symbolic Path dialog box with the proper path.
  - Click the 'Do It' button, which installs the new path choice into the Save Graphics File window ( Figure 12.8).



**Figure 12.9.:** Choosing a Path in Igor from the main menu: *Misc > New Path*.

## 12.6. Offset Scan Area

Images can be offset in a number of ways, including:

1. Manually typing in a value in each of the X, Y offset parameter fields. For more precision, pull up the Igor cursors (Ctrl + I) to get the coordinates of the center of the offset image.
2. Use the hamster wheel, with activated radio buttons, to move. This method is not recommended since Igor's update time is the rate limiting step in the process.
3. Place the cursor somewhere on the image or the trace/retrace line, left-click the mouse button and choose X, Y, X-Y offset, or Offset Next Scan from the pop-up menu that appears.
4. Use the Igor cursors (Ctrl + I) to determine the exact XY coordinate of the center of the offset; type these values into the X & Y offset fields.
5. During the same imaging session, the tip can be offset by performing a similar operation as described in option 3 on saved data. The tip will go to the area selected, thanks to the superior closed-loop sensors.

**Note** If the head has been touched for some reason between the time the image was captured and type of offset is selected, there is a great possibility that the head had been displaced by a couple microns in the kinematic divots on the baseplate.



## 12.7. Zooming Images

There are multiple ways to zoom on an image feature, including:


1. Manually type a value for the Scan Area. This will zoom with the same image center as the previous scan area, unless you manually offset, which would be a separate step (see 12.6).
2. Use the mouse cursor to draw a box around a feature of interest on the image. Note that the Igor mouse cursor changes shape when a popup menu is available. Right-click and select one of the following menu options:
  - a) *Zoom Zoom* To pick the smallest side of the (imprecise) rectangle that was created. 12.10 shows an example of this.
  - b) *[edit] Zoom Zoom* For more a more precise zoom; manually type in the 'Desired Size' value. Click 'Do It' to zoom to the desired location.



Zoom Zoom Dialog box, opened from an image.

- c) *Nice Zoom Zoom* To round to a “nice” number for the scan area.

An example of zooming can be seen below ( Figure 12.10).



**Figure 12.10.:** This is an example of zooming on a feature. Image: Moth eye in air using AC mode.

## 12.8. Multiple Offset/Zoom Parameters


You can determine multiple offset/zoom parameters two different ways:




1. Use the Igor cursor (Ctrl + I) to determine XY coordinates before offset zoom.
2. From a larger scan area image, zoom on one area, and then use the larger scan area saved data to zoom on second feature.

In some instances, it is desirable to zoom into two separate areas/features of interest from a larger scan area. The only problem is knowing the exact coordinates of the scan area image centers beforehand so they can be typed in as Offset values when the time comes.

This is easy to accomplish in Igor: Press Ctrl + I to bring up the Igor cursors, then place them in the areas of interest to determine the X,Y coordinates within the scan area.



**Figure 12.11.:** Using the Igor (Ctrl + I) cursors to find the XY coordinates within a scan area during real-time imaging.




**Figure 12.12.:** Copy coordinates in the Command line as comments.

It is a good idea to record these numbers so you can use them in the future, especially for the second image zoom/offset. In the Command line, type a comment beginning with two forward slashes (/). You can easily retrieve these values by scrolling up into the History window (Figure 12.12).


Zooms/offsets can also occur from a saved image.

**Note** To be of significance, the saved image must be viewed while the instrument is still scanning within that 90 um XY scan area. Figure 12.13 shows an example of this.




(a) Real time scan to zoom-in on one area.


(b) Using saved data of initial real-time data to zoom on another area of interest.



(c) Large saved image.



(d) Real-time zoom image (panel A).



(e) Zoom from saved image (blue box in large saved image).

**Figure 12.13.:** Zoom/offsets from saved data.

## 12.9. Real Time Section Analysis:

### 12.9.1. Basic Real-Time Line Section Analysis:

To measure a section line analysis during a real-time image acquisition, you can use the Analyze Panel (Figure 12.14).






Figure 12.14.: Section tab of the Analyze Panel.

1. Go to *MFP IP > Analyze Panel*.
2. Select the 'Include RealTime' check box.
3. *Line* will be selected by default for the 'Mode', because it is unable to do free hand lines in real time.
4. Click the 'Draw' button. Pull the cursor over the image feature to draw a line on the area of interest in the image (Figure 12.15).
5. Select the 'Full Width' check box for sections to traverse to the edge of scan areas along the line vector drawn (Figure 12.14).
6. Multiple scan lines can be averaged to eliminate noise or give more statistical relevance to the sectioned line.
7. Once the line is drawn, a Section plot appears, similar to that shown in Figure 12.15. You can place the Igor cursors (Ctrl + I) on the curve to measure points on the line, displayed in the lower half of the Analyze Panel. (Notice the blue cursors on the feature in the image being measured at FWHM.)



(a) Drawn line on image




(b) Corresponding section analysis, with user-positioned Igor cursors. The lower portion of Analyze Panel gives specifics of line relative to cursor locations.

**Figure 12.15.:** Real time section analysis

If the Y axis data values do not match the values in the Sum and Deflection Meter, some type of real-time flattening is probably activated that will offset the data. Switch the real-time flattening ('Live Flatten') to *None* in the Master Channel Panel (Figure 12.2).

If the image happens to possess some tilt, it won't be intuitive to get a proper height from the cursors. To remedy:

1. On the section plot, hold down the Ctrl + Shift keys over the curve to be "leveled" (Figure 12.16).
2. Right or left mouse-click over the line scan; a double arrow appears.
3. Move the mouse up or down to make the section level. From here, you can put the Igor Ctrl + I cursors on to measure the dY.
4. Click the 'Export Table' button in the Section window to get these cursor stats (i.e., Xa, Ya, Za; Xb, Yb, Zb; plus differential values and angle values).





**Figure 12.16.:** Real-time section leveling to obtain more accurate height differentials.

### 12.9.2. Adjusting the Angle of the Section Line


1. Click the 'Clear' button in the Analyze Panel and redraw a new line over the feature of interest.
2. Grab one end of the Igor cursors on the section line in the image by left mouse-clicking and hold while dragging the cursor to new point (Figure 12.17). This will also update the Angle value.
3. Select the 'Angle' check box to manually adjust the Angle value.





**Figure 12.17.:** Adjusting line section angles


To average multiple scan lines, enter the desired amount of scan lines to be averaged for the Width in the Section tab of the Analyze Panel. Lines will be displayed with a solid center indicating center and dotted lines representing the averaging scan line width.



**Figure 12.18.:** Averaging multiple scan lines (a 25 line average in a 512 x 512 image resolution)

### 12.9.3. Plotting Multiple Section Lines on One Section Plot


1. Create a section line as described above.
2. Click the 'Take a "Snap Shot" of active trace and append to upper graph' button.
3. Create an additional section line.
4. Click the 'Take a "Snap Shot" of active trace and append to upper graph' button again.




(a) Create first section plot



(b) Click "Snap Shot" button in Sections panel



(c) Create an additional section plot, and then click "Snap Shot" button again



(d) Resulting multiple section plots (Red: single line, Blue: averaged line, Green: single line) (Sample: collagen in air)

**Figure 12.19.:** Plotting multiple sections on one section plot

### 12.9.4. Export the line section XY data:

Click the 'Edit' button in the Section plot. This will bring up a table with a Y value relative to the points in the Igor Wave (see the Igor Getting Started manual to understand the definition of a wave), and the blue cursor points. Take the following steps:

1. Go to *Windows > New Table*.
2. Select 'SectionWaveX & SectionWaveY'.
3. Click 'Do It'. This will create an X,Y delimited table from which numbers/data can be cut and pasted.


## 12.10. Other Features

### 12.10.1. Windows Manager

The windows manager helps find windows when you have many different panels and tables open.

- Ctrl + Shift + 1 opens the window (Figure 12.20).
- Find the panel, graph, or table you are interested in and double-click it. This will bring that window forward.





**Figure 12.20.:** The Windows Manager Panel; Igor hotkey (Ctrl + 1)

### 12.10.2. Saving Igor Experiments

Saving the Igor Experiment is a great way to begin right where you left off. When you save an experiment, all the panels, graphs, etc. will come up just as it was exited. Other advantages to saving them is that every command executed is stored in the History window. Recall that Igor is memory-based, so many full days of this will make the file large and ultimately slow the program.

One downside to this feature is that .pxp experiments can be fairly large files. If you save pxp's on a regular basis, it is a good idea to back up these files or remove them from the lab PC so it doesn't get bogged down with a full hard drive.


One of the best ways to open a saved experiment is to open the AR software, then go to *File > Open Experiment* to select the saved experiment you want to open.

**Note** The experiment will open if you used the same software version (or slightly later) to open the saved experiment.

### 12.10.3. Adding Bars to the Sum and Deflection Meter

To view additional signals in the Sum and Deflection Meter:

1. Click the 'Setup' button.
2. Change the parameter of interest 'Status' from to "On".
3. Click the 'Done' button.



**Figure 12.21.:** Turning on additional meters in the S&D meter.

Figure 12.21 shows the phase signal being selected, such that the quantitative value of the signal can be seen during imaging. This will show whether the tip is experiencing net attractive or repulsive forces as it interacts with the surface.

### 12.10.4. Modifying Igor Graphs

To define or adjust the thickness, color, or appearance of a line in any Igor Pro window, double-click the mouse on the plotted line. The Modify Trace Appearance dialog will come up (Figure 12.22) where you can make any changes. When finished, click the 'Do It' button.




Figure 12.22.: Modify Trace Appearance dialogue.

## 13. Exporting Raw Data

### Chapter Contents

13.1 Exporting Your Data in ASCII Format . . . . .	173
--	-----

### 13.1. Exporting Your Data in ASCII Format

You can export any wave in Igor. 90% of the battle is finding the wave.

To export a wave:

1. Double click on the wave. You get a Modify Trace dialog.
2. Hold down the left mouse button on the name of the wave that looks like the correct wave. The box below will show the data folder that is holding that wave.
3. A different way to find the folder is to press Shift + F1, then click on the trace. This displays a tooltip type box with the names of the waves and the data folders.
4. Go to the *Data > Data Browser* menu. Locate the data folder that was shown by the dialog. Alt + left-click to the left of that data folder. There should now be a red arrow pointing at that folder letting you know that you are in that folder.
5. Go to the *Data > Save Waves > Save Delimited Text* menu. Select the wave from the list of waves in that data folder.
6. Click 'Do It'. Then, select where you want to save the .dat file.
7. Note that images are layered waves, 3 dimensional. Ascii can only handle 2D waves, so the layers are stacked in the rows.

#### Example:


- Make/N=(3,4,2) LayeredImageWave
- LayeredImageWave = P+Q\*10+R\*100 [P is row index, Q is column Index and R is layer Index, all indexes start at zero]
- Save/J/M="\r\n" LayeredImageWave as "LayeredImageWave.txt" saves the .dat file with this in it:  

```
0 10 20 30 1 11 21 31 2 12 22 32 100 110 120 130 101 111 121 131 102 112 122 132
```
- As you can see, the 4th row and above are > 100; they came from the second layer of the igor wave.
- You can either deal with all the data in one wave, pulling the layers out that you want, or you can save only the layer you care about by first extracting the layer (with the command popup, extract layer) and then saving the wave LayerData as the ascii.

Exporting to matlab:

- Export Data. As images and as other raw formats used with Matlab.
- C:\Program Files\WaveMetrics\Igor Pro Folder\Technical Notes\Igor Tech Notes\TN003 Igor binary format\TN003.ifn a technical note on the Igor file format.





## Part III

# Programming

**Part III: Who is it for?** The AR SPM software can be customized for your needs in a variety of ways. Custom control panels and macros can be set up with no programming needed. But if you want even more control you can write your own Igor code. This section will cover the basics of writing Igor code to control the AFM.



# Part Contents

<b>14 Basic Programming Tutorial</b>	<b>178</b>
14.1 Introduction	178
14.2 Conventions	178
14.3 Help	179
14.4 Low Level Parameter Names	179
14.5 High Level Parameter Names	180
14.6 The Crosspoint Switch	181
14.7 Parameter Control and Collection	183
14.8 Drive the XY Stage in a Circle	187
14.9 Feedback Loops	188
14.10 Setting Up Your Own Z Feedback Loop	189
14.11 Creating Your Own Programs	191
14.12 More Programming Resources	192
<b>15 Using the Command Line</b>	<b>193</b>
15.1 Command Line Favorites	193
<b>16 Help! Bug Reporting</b>	<b>196</b>
16.1 User Guides (Manuals) and Online Help Files	196
16.2 Strategies for Finding the Help You Want	197
16.3 AR SPM Software Help Specifics	199
16.4 Using Igor Pro Help	199
16.5 How to Deal with Error Messages	200

# 14. Basic Programming Tutorial

CHAPTER REV. 2436, DATED 09/04/2021, 14:34.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

## Chapter Contents

14.1	Introduction . . . . .	178
14.2	Conventions . . . . .	178
14.3	Help . . . . .	179
14.4	Low Level Parameter Names . . . . .	179
14.5	High Level Parameter Names . . . . .	180
14.6	The Crosspoint Switch . . . . .	181
14.7	Parameter Control and Collection . . . . .	183
14.8	Drive the XY Stage in a Circle . . . . .	187
14.9	Feedback Loops . . . . .	188
14.10	Setting Up Your Own Z Feedback Loop . . . . .	189
14.11	Creating Your Own Programs . . . . .	191
14.12	More Programming Resources . . . . .	192
14.12.1	Manuals . . . . .	192
14.12.2	Help Files . . . . .	192
14.12.3	User Forums . . . . .	192
14.12.4	Give us a call . . . . .	192

## 14.1. Introduction

The purpose of this tutorial is to give you a start on learning how to control the AFM at a low-level using Igor. You might want to do this to make the AFM do things not available in the graphical user interface (gui) of the software. Before running through this tutorial you should be very comfortable using the AFM and, ideally, have practiced a little bit of Igor programming. At the very least, it is important to understand the difference between a wave's scaling and wave indexes and doing waveform assignments from the Command line. The "Getting Started" tutorial in the Igor manual can teach you this much. This tutorial focuses on commands specific to the Asylum AFMs. At the end of this tutorial, we'll discuss some references for more programming commands.

## 14.2. Conventions

Anytime in this tutorial where you see something indented and in a different font, such as

```
print td_ReadValue("Head.Temperature")
```

it is meant to be executed from the Igor Command line. You can just copy and paste text, but you will learn more quickly if you type them out yourself. The value that the function returns to the Igor history is often shown below the command, for example

```
print td_ReadValue("Head.Temperature")
31.625
```

This is just to show you what you might expect to see. You should not enter the 31.625 into the Command line.

### 14.3. Help

The majority of commands you will be learning about are commands that start with `td_`. The `td` stands for 3D and/or Todd Day, the programmer that created the commands. Help for these commands can be found in the Xop Help file. The easiest way to get to this is to go to the menu *Help > Command Help*. This brings up the Help Browser, and you can just scroll down to get help on the command you are interested in learning about. Make sure the 'Functions', 'Operations', and 'Programming' check boxes are checked.

Another very useful help item is that, if you type the command into the Command line and right-click on it, you can go to the Help for that command or insert a template. Inserting a template is very useful if a command takes several parameters and you don't remember the order.

### 14.4. Low Level Parameter Names

Everything on the instrument that has a value to read or control has a parameter name. These parameter names are constructed as `Device.[Group(s)].Parameter`. For example, to read the value of a thermometer that is built into the MFP-3D head, you would execute the following:

```
print td_ReadValue("Head.Temperature")
31.625
```

The *head* refers to a specific device that is part of the instrument, while the *temperature* is a specific parameter that the head owns. Examples of devices are the head, scanner, and ARC. Additional devices, such as heaters, can also be added to the instrument. You can see all the devices hooked up to the instrument by clicking on the gear icon at the bottom-left side of the screen.

As many devices have a lot of parameter names, there is an additional level of organization called *groups*. For example, the ARC has a lot of analog-to-digital converters (ADC's) and digital-to-analog converters (DAC's) in it. There is a group named "Input" that contains the ADC's, and a group named "Output" that owns the DAC's. Furthermore, the controller device is often omitted from the address, as it is the default device to read from. So, for example, to read the value of the X sensor, you would execute the following:

```
print td_ReadValue("Input.X")
-6.84063
```

Similarly,

```
print td_WriteValue("Output.X", 60)
0
```

would apply 60 volts to the X piezo. Note that when you use commands like `td_WriteValue`, where you don't expect an output, you should still precede them by a `print` command. They will return 0 if everything worked, and a non-zero error code if it didn't. For example, if you misspell something, you will probably see a 3 or 5 returned.

Note that some outputs can also be read as inputs. You might be interested in knowing what the current voltage to the X piezo is before you change it. You execute:

```
print td_ReadValue("Output.X")
60
```

Another way to access these parameters is to use the menu that pops up when you click on the gear icon at the bottom-left side of the screen. A pop-up menu with icons of the devices currently hooked up will appear. For example, to see the changes you just made, float the cursor over the icon of the ARC in the menu, and an additional pop-up menu appears. Float the cursor over the parameters menu, and yet another pop-up menu appears with a list of all the devices owned by the ARC. Select *Output*, and a panel will pop up that contains all the parameters in the output group. Click on the 'read' button, and you should see that the X output has a value of 60 V.

Before you start, the HighVoltage Relays must be checked to see if they are open or closed. Click on the 'ARC' icon in the toolbar and navigate to *ARC > Parameters > Default*. Then, use the 'read' and 'write' buttons to look at and set the state of the *HighVoltageXYRelay* and *HighVoltageZRelay*. You can also change the state of the relays at the command link as follows:

```
print td_WriteValue("HighVoltageXYRelay",1)
0
OR
print td_WriteString("HighVoltageXYRelay","Closed")
0
```

At the top of the Help file, there is a complete list of parameter names as well as the error codes.

## 14.5. High Level Parameter Names

There are also parameters that exist in Igor rather than down at the device level. If you select *Programming > Global Variables > Master*, you will see a table full of parameters and descriptions of what they control. For example, the top item is the *ScanSize* and shows the scan size in meters under the value column.

Though you can read these values from the table, while programming it is useful to access them with the GV (Get Value) routine. Executing

```
print GV("ZPiezoSens")
```

will print the sensitivity of the Z piezo, in meters per volt, into the history. See Table 14.1 for a list of other get functions.

Table 14.1.: Get Functions

Function Name	Description
GDS('Parameter Name')	Get Description String
GFS('Parameter Name')	Get Format String
GUS('Parameter Name')	Get Unit String
GTS('Parameter Name')	Get Title String
GPS('Parameter Name')	Get Panel String
GV('Parameter Name')	Get Value
GVU('Parameter Name')	Get Units
GVL('Parameter Name')	Get Value Low
GVH('Parameter Name')	Get Value High
GVMU('Parameter Name')	Get Value Min Units
GVS('Parameter Name')	Get Step Size

## 14.6. The Crosspoint Switch

To make the instrument much more flexible, the controller has something called a *crosspoint switch* inside of it. The crosspoint switch is a lot like an old telephone patch board. It has 16 inputs and 16 outputs, and you can connect virtual wires between any of the inputs and any of the outputs. You can have multiple wires running from one input to many outputs.

The crosspoint switch can make programming somewhat more complicated because a given ADC isn't always connected to the same things. From the main menu bar, select *Programming > Crosspoint Panel* to bring up a panel that allows you to see and control the wiring of the crosspoint switch. Switch between Contact mode and AC mode to see how some of the signals are rerouted.




Figure 14.1.: Crosspoint Panel

When programming, you may find yourself using the BNCs on the front of the controller. These appear on the crosspoint switch inputs and outputs. To help keep you sane, we named all the extra ADCs and DACs with letters (e.g., InA) and the physical BNCs with numbers. So, for example, you can connect the BNC In0 to the ADC InA using the top dropdown menu in the Crosspoint Panel, and then clicking the *Write Crosspoint* button.

The crosspoint switch can also be set with `td_WriteString` commands. For example,

```
print td_WriteString("Crosspoint.InA","BNCIn0")
0
```

would accomplish the same thing.

Probably the most confusing option on the crosspoint is the 'FilterIn'. By routing a signal to 'FilterIn', you are routing it into a low pass filter of about 36 kHz. The filter signal then reappears as an input to the crosspoint switch named 'FilterOut'. This filter is useful for filtering high bandwidth signals (e.g., Defl) before passing them to a 50 kHz ADC like 'InA'.

#### A simple example with `td_ReadValue` and `td_WriteValue`

You can use the `td` commands to move the X stage around and measure its motion with the X

sensor. The stage is approximately centered at 70 volts. If you execute the following commands in a row, you should see output similar to those listed below:

```
print td_WriteValue("Output.X",60)
0
print td_ReadValue("Input.X")
-.0396402
print td_WriteValue("Output.X",70)
0
print td_ReadValue("Input.X")
-0.698834
print td_WriteValue("Output.X",80)
0
print td_ReadValue("Input.X")
1.51577
```

You just drove the X stage to 60, then 70, then 80 volts. As you can see, the sensor voltage measured this movement and returned three different values. To see what happened in meters instead of volts, we can use the sensitivities to convert the voltages.

```
print GV("XPiezoSens")*20
1.07474e-05
```

This means moving the stage from 60 to 80 volts (20 volt change) should have produced about 9.6 microns of movement. In reality, the sensors measured:

```
print (1.51577 - -.0396402)*GV("XLVDTSens")
1.1625e-05
```

This is about 11.6  $\mu\text{m}$ . The fact that the stage moved more than we expected is the reason we have the sensors. How far the stage movement is off is dependent on scan size. Try the above experiment for a 2 volt change (drive X to 69 volts and then 71 volts), and you should get better agreement.

While you could make the sample move around by stringing together a lot of `td_WriteValue` commands, that is tedious. A better way is to use an Igor wave to control the stage and collect data. This is described in the next section.

## 14.7. Parameter Control and Collection

It is good practice before sending instructions to the ARC to stop current instructions. The `td_stop()` command does this but in general it should not be used as it stops all outwaves, inwaves, and feedback loops.

For now let's just use the `td_stop` command:

```
print td_stop()
0
```

Now, make two Igor waves, one to store voltages to drive the X stage and one to store the sensor values we read:

```

Make/N=1024 PiezoVoltage SensorVoltage
Next make a graph of each
Display PiezoVoltage
Display SensorVoltage

```

Move the graphs so you can see both of them. You'll notice that, right now, the wave scaling of both is just p scaling. The X axis on each wave is from 0 to 1023. To use this wave to drive the piezo voltage, it is much more convenient to have the x scaling in time. While you could set it by hand, the easiest way is to use a td command. To set up an output wave, you use the `td_xSetOutWave` command. You can find detailed help for the command described in [Section 14.3 on page 179](#).

The syntax of the command is:

```

td_xSetOutWave(whichBank, eventString, channelString, wave,
               interpolation)

```

The parameters are:

**whichBank** Determines which memory bank in the DSP the wave will be stored in. There are three banks of memory, labeled 0,1, and 2, and each can hold up to two pairs of waves. These waves can be up to 87,380 points long and must be single-precision (32 bit) floating point.

**eventString** Determines when the data will be output. You can think of them as a trigger. You usually don't want the wave to be output as soon as it reaches the DSP because you might want to synchronize the output with some data collection. The events are group of parameters owned by the controller, so a typical eventString would be "Event.0", but for functions that have a specific event argument, it is simply "0", the "event." can be dropped. When we later change the value of this event parameter, the wave will be output to the DAC that drives the X voltage.

**channelString** A low-level parameter name like we learned about above. In this case, it is "Output.X".

**wave** Name of the Igor wave with the data to be sent. Here it is "PiezoVoltage".

**interpolation** Determines how quickly the wave are output. The DAC's and ADC's on the ARC all run at 50 kHz, so an interpolation of 1 means that each point in the Igor wave corresponds to a 50 kHz sample. For the 1024 point waves you made, this means the entire wave would be output in about 20 milliseconds. This is very fast, so we will choose an interpolation factor of 100 to give us about a 2 second output.

**Note** Older versions of the AR SPM software (prior to version 101010) ran the ADCs at 100kHz. In that case, the output of the above example would take only 1 second.

Before we send the wave to the ARC, we need to populate it with values. Currently, each point is set to zero. In this case, lets make a ramp over the full range of the piezo (-10V to 150V):

```

PiezoVoltage[0,511] = -10 + 160*p/511
PiezoVoltage[512,1023] = 150 - 160*(p-511)/511

```

Putting this all together gives:

```


print td_xSetOutWave(0, "0", "Output.X", PiezoVoltage, 100)
0

```



If you did everything right, you should have gotten a zero (0) in the history.

If you look at the graph of PiezoVoltage that you made, you will notice that the x scale has been updated to reflect the time that the wave will take to output. Shown in [Figure 14.2 on page 185](#).



**Figure 14.2.:** Piezo Voltage vs. Time

Now you need to set up the SensorVoltage wave to record data. The syntax of the command is fairly similar to that of `td_xSetOutWave`:

```
td_xSetInWave(whichBank, eventString, channelString, wave, callback,
              decimation)
```

The first four parameters are conceptually identical except `channelString` will now be an input instead of an output.

We also have a `callback` parameter and `decimation` instead of interpolation. *Decimation* is the opposite of interpolation. Rather than creating points, we want the DSP to reduce the number of points coming back. Sometimes getting data at 50 kHz is a bit like drinking from a fire hose, and you want less information.

The *callback* function is an Igor function that will be called when the input data has been collected. This is useful so that control returns to Igor while the data is being collected (otherwise Igor would freeze during data collection). This function might be something that does some analysis on the data.

To show how the callback works, let's make a very simple example that prints "Done" when the callback hits. To do this, we must write some code not in the Command line. Select the *Windows Menu > New > Procedure*. Leave the name as is and click the 'New' button.

Type the following onto the procedure window:

```
Function InWaveCallback()
  print "Done"
End
```

Now when the callback function hits, it should print "Done" to the Command line.

One major conceptual difference here is that, rather than the data being stored on the DSP like the output data, it is continually streamed back to the PC. This means the banks correspond to "pipes"

on the USB connection rather than memory. There are 3 banks, each bank brings back 1 or 2 channels of 24-bit data.

Putting all this together gives:

```
print td_xSetInWave(0, "0", "Input.X", SensorVoltage, "InWaveCallBack", 100)
0
```

You will see that the x scaling on the SensorVoltage graph has now changed to seconds. You now have set up everything and all you have to do is trigger Event.0. For this, it is easier to use td\_WriteString instead of td\_WriteValue. They are identical except the first is used for passing strings instead of numbers.


```
print td_WriteString("Event.0", "once")
0
Done
```

If all went well, after about a second, you should see the SensorVoltage wave get updated with data, and “Done” printed in the Command window. The curvature shows that the piezo isn’t moving linearly. If the graph shows noise, then the high-voltage relays are in the open state, so the voltage is not getting to the piezos. You can check this by clicking the ‘Controller Icon’ button on the Igor window and navigating to *ARC > Parameters > Default*. Then use the ‘read’ and ‘write’ buttons to look at and set the state of the *HighVoltageXYRelay* and *HighVoltageZRelay*. You can also change this at the Command line as follows:

```
print td_WriteValue("HighVoltageXYRelay",1)
0
print td_WriteValue("HighVoltageZRelay",1)
0
```

If you graph the SensorVoltage vs. the PiezoVoltage, you will see the hysteresis loop traced out by the stage (see [Figure 14.3 on page 186](#)), using:

```
display SensorVoltage vs PiezoVoltage
```



**Figure 14.3.**

## 14.8. Drive the XY Stage in a Circle

The set of commands below will drive the XY stage with voltages corresponding to a circle almost full scale for the stage. The sensors will measure what really happened. The example uses `td_xSetInWavePair` and `td_xSetOutWavePair`, both of which are very similar to the commands you learned above. You can watch the sample movement using the optics of your microscope.

```
make/N=(1024)/0 XVoltage YVoltage XSensor YSensor
display XVoltage; appendtograph/R YVoltage
display XSensor; appendtograph/R YSensor
print td_xSetOutWavePair(0, "0,0", "Output.X", XVoltage,
    "Output.Y", YVoltage, 100)
//0
XVoltage = 70 + 70*cos(2*pi*x/1.024) //X is seconds, 1.024 is the time per cycle
YVoltage = 70 + 70*sin(2*pi*x/1.024)
print td_xSetOutWavePair(0, "0,0", "Output.X", XVoltage,
    "Output.Y", YVoltage, 100)
//0
print td_xSetInWavePair(0, "0,0", "Input.X", XSensor,
    "Input.Y", YSensor, "", 100)
//0
print td_WriteString("Event.0", "once")
//0
print td_WriteString("Event.0", "once")
//0
display YSensor vs XSensor
ModifyGraph width={Plan,1,bottom,left}
```




Figure 14.4.

As you can see by Figure 14.4 on page 187, because of piezo non-linearities and hysteresis, the circle isn't very round.

## 14.9. Feedback Loops

**Note** This section is very advanced, and you run the risk of damaging your microscope sensors if using unstable gains.

For a lot of things that you might want to do, such as make a round circle, you will want to use a feedback loop. The controller is able to run six PIDS loops. These loops have the following gains:

- Proportional gain P
- Integral gain I
- Differential Gain D
- Double-integral gain S

Each parameter belonging to the loop can be accessed;. For example, *PIDSLoop.0.IGain* is the name for the Integral gain of loop 0.

One problem with feedback loops is that, if you pick the wrong gains, the feedback loop can oscillate. Besides being incredibly loud, this can also push the sensors out of alignment.

The easiest thing to do is to look at the values that the software uses. Be careful not to look at the gains in the panels themselves (like integral gain on the Main panel or the X and Y gains on the XY gains panel). These are scaled to make it more convenient for users to think about.

To peek at the value that the software uses for Z in contact mode, click 'Engage' on the meter panel. It doesn't matter if you have a cantilever in or not. This turns on the Z feedback loop. The following commands will then read the three gain values. This result is typical for contact mode.

```
print td_ReadValue("PIDSLoop.2.PGain")
0
print td_ReadValue("PIDSLoop.2.Igain")
1000
print td_ReadValue("PIDSLoop.2.SGain")
0
```

It is easiest to look at their values in the PIDS Loop Panel. Go to *Programming > XOP Tables > PIDS Loop Panel*.

For AC mode, the gain will be the opposite sign and smaller.

Similarly, to see the X and Y gains, click 'Do Scan' and look at the PIDS Loop Panel, or you could execute the following:

```
print td_ReadValue("PIDSLoop.0.PGain")
0
print td_ReadValue("PIDSLoop.0.IGain")
6887.8
print td_ReadValue("PIDSLoop.0.SGain")
1075803
print td_ReadValue("PIDSLoop.1.PGain")
0
print td_ReadValue("PIDSLoop.1.IGain")
6357.8
print td_ReadValue("PIDSLoop.1.SGain")
1007645
```

In almost all software versions, the X loop is on the first PISLoop (PISLoop.0), the Y is on the second (PISLoop.1), and Z is on the third (PISLoop.2). There is a small chance that the X and Y are on PISLoop.3 and PISLoop.4. You would see this if you went to read the gains.

To set up a feedback loop, you use the `ir_SetPISLoop` command. This is because the PIDSLoop panel is automatically updated by the AFM software and does not allow manual overwrite. The syntax of the `ir_SetPISLoop` command is:

```
ir_SetPISLoop(whichLoop, eventString, inChannelString, setpoint, pgain,
              igain, sgain, outChannelString, OutputMin, OutputMax)
```

The parameters are:

**whichLoop** Determines which of the six feedback loops you will use. They all run at 50 kHz.

**eventString** Similar to the events described above. Unlike with input and output waves, with feedback loops you often want them to start running immediately. Passing a string of "Always" for the eventString will make this happen for the start string; there also needs to be a comma and a string for the stop event, in our examples "Never".

**inChannelString** The channel string for the input to the feedback loop. For example, in Contact mode, this is one of the ADCs (InFast usually), which is wired up to Defl.

**setpoint, pgain, igain, sgain** The setpoint and gains for the loop. After the loop is set up, the setpoint can be driven by a wave. This would be used, for example, to move the XY stage around under feedback control. Note the gains are scaled to be mathematically right. A pgain of 1 will give you a 1-volt output for a 1-volt input. Similarly, the same 1-volt input with an igain of 1 will give you a ramp from 0 to 1 in one second.

**outChannelString** The output of the feedback loop. For Contact mode, this would be the Z DAC.

**OutputMin OutputMax** The range of values used for the OutputChannel.

### 14.10. Setting Up Your Own Z Feedback Loop

To set up your own feedback loop in Z, first turn on the meter and make sure you have a deflection close to zero. Also make sure the cantilever is close to the surface by doing an engage the normal way and then withdrawing. Above, you saw that for contact mode, typical gains are P=0, I=1000, and S=0.

First you hook up the Defl channel to the InA ADC by executing:

```
print td_WriteString("Crosspoint.InA", "Defl")
0
```

Next, you set up the feedback loop with:

```
print ir_SetPISLoop(2, "Always, Never", "Input.A", 1, 0, 1000, 0,
                    "Output.Z", -10, 150)
0
```

If all went well, you should see on the meter that the Zvoltage went to a stable value and that the deflection is at 1 volt.

**Fig Closed Loop Circle** We need to figure out what the gains are for X,Y, and Z. This is done with the same read value commands as above. The microscope should be engaged and scanning.

td\_RV is a shorthand function that calls td\_ReadValue.

```
variable/G X_PGain, X_IGain, X_SGain, Y_PGain, Y_IGain, Y_SGain,
      Z_PGain, Z_IGain, Z_SGain
X_PGain = td_RV("PIDSLoop.0.PGain")
X_IGain = td_RV("PIDSLoop.0.IGain")
X_SGain = td_RV("PIDSLoop.0.SGain")
Y_PGain = td_RV("PIDSLoop.1.PGain")
Y_IGain = td_RV("PIDSLoop.1.IGain")
Y_SGain = td_RV("PIDSLoop.1.SGain")
Z_PGain = td_RV("PIDSLoop.2.PGain")
Z_IGain = td_RV("PIDSLoop.2.IGain")
Z_SGain = td_RV("PIDSLoop.2.SGain")
printf "X Gains: P:%.4g I:%.4g S:%.4g\r", X_PGain, X_IGain, X_SGain
printf "Y Gains: P:%.4g I:%.4g S:%.4g\r", Y_PGain, Y_IGain, Y_SGain
printf "Z Gains: P:%.4g I:%.4g S:%.4g\r", Z_PGain, Z_IGain, Z_SGain
```

Now let's make a variable Radius (the size of the Radius circle in volts). The range of the sensor is between -10V and +10V, but they don't use the full range. Try using 3V for a radius:

```
variable Radius
Radius=3
```

We then need to make all the waves as follows:

```
Make/N=(1024)/0 XVoltage YVoltage XSensor YSensor XCommand YCommand
Display/K=1 /W=(5.25,41.75,399.75,250.25) XVoltage
Appendtograph/R YVoltage
Display/K=1 /W=(7.5,275.75,402,484.25) XSensor
Appendtograph/R YSensor
Display/K=1 /W=(409.5,41.75,662.25,250.25) YSensor vs XSensor
ModifyGraph width={Plan,1,bottom,left}
```

Now set up the circle commands:

```
XCommand = Radius*cos(2*pi*p/1024)
YCommand = Radius*sin(2*pi*p/1024)
```

We should do a td\_stop() to make sure nothing else is running:

```
print td_stop()
0
```

We now set up the feedback loops for X and Y. Notice how the initial setpoint is set to the Radius for X and 0 for Y. This is because the XCommand starts at a voltage equal to the Radius and the YCommand starts at 0V.

```
print ir_SetPISLoop(0,"Always,Never","Input.X",Radius,X_PGain, X_IGain,
```

```

    X_SGain,"Output.X",-10,150)
0
print ir_SetPISLoop(1,"Always,Never","Input.Y",0,Y_PGain, Y_IGain,
    Y_SGain,"Output.Y",-10,150)
0

```

We then pass the feedback loop voltages that we would like by using a command signal, and it will vary the XVoltage and YVoltage to try to achieve this. We use the trick to look at outputs (Output.X) as inputs to see what the voltages are and similarly display the input voltages.

```

print td_xSetOutWavePair(0, "0,0", "PIDSLoop.0.Setpoint", XCommand,
    "PIDSLoop.1.Setpoint",YCommand,100)
0
print td_xSetInWavePair(0, "0,0", "Output.X", XVoltage, "Output.Y",
    YVoltage, "", 100)
0
print td_xSetInWavePair(1, "0,0", "Input.X", XSensor, "Input.Y",
    YSensor, "", 100)
0

```

As you can see in [Figure 14.5 on page 191](#), the circle should be perfectly round.




Figure 14.5.

## 14.11. Creating Your Own Programs

While all the examples here have been performed from the command line, except the callback function, it is pretty easy to start copying and pasting commands to make your own functions inside procedure windows.

## 14.12. More Programming Resources

### 14.12.1. Manuals

The AR SPM software is written in a third party software environment called Igor. There is a pdf of the Igor manual in <Igor Pro Folder>\Manual\IgorMan.pdf, where <Igor pro Folder> is the location where the Igor software is installed. The manual starts with the all-important “Getting started” guide. Igor is supplied by Wavemetrics, and they have an excellent website (<http://www.wavemetrics.com>) with additional information like a knowledge base and a user forum.

### 14.12.2. Help Files

All of the AFM-related user interface panels were programmed by Asylum Research programmers using the tools available in the Igor environment. To make that easier, many special functions are subroutines based on Igor commands. An up-to-date list of these Asylum Research functions can be found in the AR SPM software in the main menu bar under *Help > AR Help Files > AR Common Functions*. Note that all these functions are stored in files that you can open and copy from to make your own custom code and user panels. You might want to explore the higher-level macro programming capabilities, but, in the end, you are allowed to delve into the source code of nearly every aspect of the AR SPM Software.

In the same location, there are other useful Help files with respect to programming under *Help > AR Help Files > AR Advanced*.

Finally, there is a list of low level commands that are mainly related to communicating with the Asylum Research hardware. You have already seen some of them in this tutorial (the ones starting with `td_`). These commands supplement the ones supplied by Igor. An up-to-date list can always be found at: *Help > AR Help Files > XOP Help Files > MFP3D XOP Help*. The MFP3D reference is historical, and the files also cover more recent instruments such as the Cypher and Jupiter AFMs.

### 14.12.3. User Forums

Years of Asylum Research SPM Discussion are archived at **[support.asylumresearch.com](http://support.asylumresearch.com)**. Once you join, you can peruse the DIY programming section of the forums or post questions. All Asylum employees are forum members, and, within a day, you will typically get an answer from the experts at Asylum, as well as from other power users around the world. The support site also hosts the latest manuals and software updates for your instrument.

### 14.12.4. Give us a call

You call us if you need additional help. And, always feel free to contact us directly at:

[Support@AsylumResearch.com](mailto:Support@AsylumResearch.com)



## 15. Using the Command Line

CHAPTER REV. -2, DATED 00/00/0000, 00:00.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

### Chapter Contents

15.1	Command Line Favorites . . . . .	193
------	----------------------------------	-----

### 15.1. Command Line Favorites

Here's a list of common commands (source: "[http://wiki.AsylumResearch.com/Common\\_Commands](http://wiki.AsylumResearch.com/Common_Commands)") that may be useful (some may be outdated):

**Subtracting Images**    `Extractlayer()` //extract layer of first image–rename it  
`Extractlayer()` //don't rename it if you are subtracting first layer from this layer  
`Layerdata = renamed layer` //kicks out subtracted image (but its upside down)  
`Newimage/F = layerdata` //makes it right  
//insert this layer back in to do modification of ARgyle on it.

**Rotating Images (saved data)**    (See example in Section [Section 7.4.2.3 on page 88.](#))  
Extract a layer from the image to be rotated, then type in the command line:  
`duplicate layerdata, filename` //what do you want to call the image  
`ilayerdata = filename[q][X-p]` //where X is equal to the number of points in the image per scan line–remember Igor starts counting at zero:  
//insert the layer back into the initial image.  
//perform an overlay in ARgyle.

**Monitor system parameters**    `print td_ReadValue("Temperature%Default@Head")` //reads temp at head  
`print td_ReadValue("Temperature%Default@scanner")` //reads temp at scanner

**Turn ON/OFF the blinking LED at XYlvdT board and head**    `print td_WriteValue("LedBlinking%Default@XYLVDT", 1) //turns ON LED at XYLVDT board`

`print td_WriteValue("LedBlinking%Default@XYLVDT", 0) //turns OFF LED at XYLVDT board`

`print td_WriteValue("LedBlinking%Default@head", 1) //turns ON LED at head`

`print td_WriteValue("LedBlinking%Default@head", 0) //turns OFF LED at head`

**Set/read output voltage to XYZ piezo**    `print td_writevalue("Z%Output@Controller", 70) //sets output to Z piezo`

`print td_readvalue("Z%input@Controller") //reads z piezo voltage`

`print td_wv("y%output", 70) //sets output to Y piezo`

`print td_wv("x%output", 70) //sets output to X piezo`

**Send voltage pulse output (OutA) using DAC (e.g., 2V, -2)**    `print td_WriteValue("A%Output", 2);print td_WriteValue("A%Output", 0)`

`print td_WriteValue("A%Output", -2);print td_WriteValue("A%Output", 0)`

**Print out what 50 pN is in Volts for Spring Constant and Involts**    `//Useful for setting a trigger during force curves.`

Variable `MyForce = 50e-12 //I want 50 pN.`

`Print MyForce/GV("SpringConstant")/GV("Involts")`

**Display captured video image**    `//First set the data folder to :root or create one like root:Video or root:Capture`

`//Then click Capture button`

`NewImage capture`

**Better way of displaying captured video image**    `Display; AppendImage capture; ModifyGraph height={Plan, 1, left, bottom}`

**Rename captured video image so it doesn't get overwritten**    In the Igor Data Browser, Rename capture pollenA.

**Generate a piecewise continuous wave**    `Make/N=65535 wave0, wave1, wave2, wave3`

`wave0 = (p<32768)?(0):(1)`

`wave1 = (p<16384)?(0):((p<32768)?(5):((p<49152)?(-5):(0)))`

**Piecewise linear**    `wave2 = (p<32768)?(-5+p*10/32768):(15-p*10/32768)`

**Output a wave to Out0 and collect from In0 in sync** //write crosspoint to connect Out0 to OutA and In0 to InA

```
print td_Stop() //usually not necessary but safest
print td_xSetOutWave(1, "2,2", "A%Output@Controller", wave2, 2)
print td_xSetInWave(1, "2,2", "A%Input@Controller", wave3, "print 42", 2)
print td_WriteString("2%Event@Controller", "Once")
//wait for it to print 42
Display wave3
```

**Note** Be sure to read the Help for td\_xSetInWave and td\_xSetOutWave on banks, decimation, and interpolation. In particular, you want the decimation and interpolation to be the same if you want to correlate input with output.

**Integrate noise spectra to find biggest contributor** //From noise panel, click Save

//In Test Results:Noise, rename folder to a unique name (X01, not X; FastOffSurf, not Fast)

//Make that the current data folder

Duplicate /O PSD Int

Int \*= Int //omit this line if you are using \*true\* PSD from Mario's new panel

//old Noise panel gives amplitude spectral density as "PSD" and must be squared

Integrate Int

Display PSD; AppendToGraph /R Int

ModifyGraph log(bottom)=1

ModifyGraph RGB(Int)=(0,0,65535)

ShowInfo

//Now look for the biggest step(s) in the blue integrated curve.

# 16. Help! Bug Reporting

CHAPTER REV. -2, DATED 00/00/0000, 00:00.

USER GUIDE REV. 2436, DATED 09/04/2021, 14:34.

## Chapter Contents

16.1	User Guides (Manuals) and Online Help Files . . . . .	196
16.2	Strategies for Finding the Help You Want . . . . .	197
16.2.1	Finding Igor Pro Help . . . . .	197
16.2.1.1	On Your Computer . . . . .	197
16.2.1.2	Web Based Support . . . . .	198
16.2.2	Asylum Research SPM Software Help Finding Strategies . . . . .	198
16.2.2.1	Peruse the User Guide Volumes . . . . .	198
16.2.2.2	On your computer . . . . .	198
16.2.2.3	On the Web or by Telephone . . . . .	198
16.3	AR SPM Software Help Specifics . . . . .	199
16.3.1	XOP Help files . . . . .	199
16.4	Using Igor Pro Help . . . . .	199
16.4.1	Status Line Help, Tool Tips, and Context-Sensitive Help (Windows) . . . . .	199
16.4.2	More Information on Igor Pro Help . . . . .	200
16.5	How to Deal with Error Messages . . . . .	200
16.5.1	Bug Reporting . . . . .	200

## 16.1. User Guides (Manuals) and Online Help Files

As we have discussed, the Asylum Research SPM Software co-exists with the Igor Pro Graphing and Analysis software. Searching for help and information may require looking in one manual or another.

The main sources of information are:

**The Igor Pro user guide** A box set of 5 volumes (books), also installed as a pdf file on any computer with the Igor Pro Software installed.

**Igor Pro “Online” help files** Only accessible from your computer’s hard drive.

**Igor Pro popup tips** Various bits of popup information about menu items, buttons, and dialog boxes.

**Asylum Research SPM Software User Guide** You’re reading it now. Comprised of several volumes and installed on your computer, available at <https://support.asylumresearch.com/forum>. Please inquire at <https://support.asylumresearch.com/> about printed volumes.

**Asylum Research “Online” help files** Blended in with the equivalent Igor Pro help files. Installed on your computer when you installed the AR SPM software extension to Igor Pro.

**Online ≠  
Internet**

The Igor Pro Manual refers to Help files which reside on your local computer as “Online.” For the sake of consistency, Asylum Research will adopt this nomenclature. Do not confuse “Online” as having anything to do with files on the internet.

## 16.2. Strategies for Finding the Help You Want

First you need to determine if it is Igor Pro or AR SPM related. The following table is a good rule of thumb for making this determination:

Igor Pro	AR SPM Software
	Any control panel on the main part of the screen when you first launched the AR SPM Software.
Anything related to menu items to the left of and including the Help menu.	Anything related to menu items to the right of the “help” menu.
Anything related to graphing and curve fitting.	Anything related to Scanning Probe Microscopy

Note that there are exceptions to this table, but this is a decent starting point.

### 16.2.1. Finding Igor Pro Help

#### 16.2.1.1. On Your Computer

To find Igor help on your computer:

1. Open the pdf version of the manual by choosing from the Igor main menu bar *Help > Manual* and click “Open Online Manual”.
2. Look in the index for the topic of your choice. Note that the pdf file allows the index to be accessed easily for each letter in the alphabet. In turn, clicking on any item in the index will jump to the relevant page in the manual volumes.

**Note** The index of the Igor Pro manual is no small feat. In the boxed set of books, it is bound as a separate booklet in two-column text covering nearly 2000 manual pages. With nearly 12,000 index entries, that averages to about six entries per manual page. An impressive resource.

If that proves unsuccessful, try the following:

1. From the main menu bar, choose *Help > Help Topics* and confirm that the ‘Sort Topics’ check box at the end of the list is checked and that the dropdown menu is set to *All Open Help Files*.
2. Look down the list of topics until you find what you are looking for. Note that mixed in will be items starting with AR and MFP3D. These are not Igor Pro related but may be useful, so also read if you are not completely sure if your issue is related to Igor Pro or the AR software.

3. As a last resort, choose from the main menu bar *Help > Search Igor Files* and try to search for your Help topic. The search appears to not sort the results on any relevance, and this method usually results in very many hits. In a pinch, it can be worthwhile to peruse the results.

Be aware of all the varieties of pop-up help described in [Section 16.4.1 on page 199](#). In the print manual, you can look up Help in the index file to learn more about Igor Pro Help features in greater detail.

#### 16.2.1.2. Web Based Support

From the main menu bar in Igor choose:

**Help > Support Web Page** A link to Igor Pro support, including contact information for e-mail or telephone support and access to web-based databases of previously answered questions.

**Help > IgorExchange** A link to Igor Pro customer forum website, including exchanged questions and ideas with a huge number of Igor users world-wide. Also searches all existing posts.

At any point also feel free to call, e-mail, or otherwise contact Asylum Research for free support.

### 16.2.2. Asylum Research SPM Software Help Finding Strategies

#### 16.2.2.1. Peruse the User Guide Volumes

The most general information can be found in the Asylum Research User Guides. Chapter titles like “AC Mode Imaging in Air” are application-driven guides that tell you how to use the software to achieve an imaging goal.

#### 16.2.2.2. On your computer

1. For specific questions on software controls, look for the tiny question mark buttons next to nearly every item on AR SPM control panels. Clicking it will take you directly to the online Help file entry of relevance.
2. Since these help files are mixed in with the Igor Pro help files, you can search all existing posts. Note that Help topics relating to Asylum Research always start with AR or MFP3D.
3. For more general questions on how to use the software, refer to the printed or online user guides. The latest updates are always available for free.

#### 16.2.2.3. On the Web or by Telephone

1. Join the our support website at [support.asylumresearch.com](http://support.asylumresearch.com). It contains FAQs, software and documentation downloads, instructional videos, and a forum where you can pose questions to fellow users and Asylum Personnel alike.
2. Call or e-mail Asylum Support. If needed, we can initiate an online session where we take control of your AR SPM software remotely. See **Getting Help**, [Section on page iii](#).

## 16.3. AR SPM Software Help Specifics

As is also the case with Igor Pro, the online Help files are stored as .ihf (igor help file) files. If your Igor Pro is installed in the standard location, most of those files are found in this directory: **<Igor Pro Folder>\Igor Help Files\ARHelpFiles**, where **<Igor Pro Folder>** is the location where Igor was installed.

- ToolTips, Ctrl+1 will toggle the tool tips on or off, and nearly every control has extra information on how it works in the tool tips.
- Question mark buttons next to nearly every SPM related control in the Software. Click that button, and the relevant topic in the relevant .ihf file will open. Scroll up and down from there to read related material. You need no additional knowledge of the Help file organization and structure to find what you need.
- Larger green Help buttons have additional documentation, typically linking to help pdfs.
- The Help button the bottom of the status bar opens the Asylum Help browser, which will link to videos, tutorials, and manuals.
- Choose *Help > AR Help Files* from the main menu bar. This presents you with a list of topics with reasonable names, organized by subtopic via nested menus. The relevant .ihf help file will open to the relevant section.

### 16.3.1. XOP Help files

You may notice *Help > AR Help Files > XOP Help Files*. These files describe a series Igor programming commands which were written by Asylum Research to allow, among other things, Igor to control and communicate with the Asylum Research hardware. If you are interested this will be discussed at greater length here:<sup>1</sup>

## 16.4. Using Igor Pro Help

### 16.4.1. Status Line Help, Tool Tips, and Context-Sensitive Help (Windows)

(Source: Igor Pro Manuals)

On Windows Igor provides three ways to get help for icons, menu items and dialog items. These are status line help, tool tips, and context-sensitive help.

**Status Line Help** The status line area at the bottom of the main Igor Pro window shows brief descriptions of icons and menu items. This help is shown automatically; you don't have to do anything to make it appear.

**Note** Not currently available for AR SPM software functions.

**Tool Tips** If you point at an icon in an Igor Pro window, a tool tip will appear after a short delay. It contains just a two- or three-word description of the button. You can adjust the delay before the tool tip appears, and the duration of display in the Help page of the Miscellaneous Settings dialog, which you can choose from the Misc menu.

<sup>1</sup> *Software User Guide, Programming and Automation..*

**Context-Sensitive Help** Context-sensitive help (sometimes referred to as *F1 help*) is displayed in a pop-up window and provides more detail than the status line help. It is accessed in several different ways depending on the type of item you need help for.

**Menus** For help on items in a menu, pull down the menu and highlight the item of interest, then press the F1 key to display the context-sensitive help window.

**Icons** For help on icons in windows such as graphs and tables, *hold down Shift and press F1*. This changes the mouse cursor to a question-mark. Click on a button or icon to display the context-sensitive help window.

**Dialog Boxes** Help for individual items in Igor's dialogs can be summoned by *clicking the question-mark button* at the right end of the dialog's title bar, then clicking on the item for which you want help.

#### 16.4.2. More Information on Igor Pro Help

To learn more about Igor Pro Help features, look up "Help" in the index of the print or pdf manual. It describes your options in greater detail.

### 16.5. How to Deal with Error Messages

If you have a problem that you suspect is due to a software bug, it is always a good idea to update your software. [Chapter 16 on page 196](#) describes how to do this in such a way that a copy of your current version of the software remains functioning. We don't require the update, but it is highly likely that in the process of trying to resolve the bug, we will request that an update is made.

#### 16.5.1. Bug Reporting

The preferred method is as follows:

1. From the main menu bar choose *Programming > Bug Report Panel*.
2. Fill in your e-mail address and your name.
3. Click the 'Edit' button and type in a description of your problem.
4. Close the Edit window.
5. Click the 'Send Report' button.
6. You will be prompted to "OK" the saving of your current experiment, which is required so it can be automatically attached to the outgoing e-mail. The e-mail also contains information on the version of Windows you are running and a few other specifics of your computer hardware.

The message will be sent out using Asylum Research's e-mail servers. Note that an experiment file can easily be tens of megabytes, so it may take a little while to upload, depending on your internet connection.



An alternative method is to describe the bug in an email and send it to <https://support.asylumresearch.com/> or call us to discuss the bug. Note that this way we won't have your experiment file containing clues as to why you are experiencing an error.



## **Part IV**

# **Bibliography, Glossary, and Index**

# Bibliography

## Cited Scientific References

## Cited Asylum Research Documents

Applications Guide, Chapter: Force Curves Acquisition., Placeholder

Applications Guide, Chapter: Spring Constant Determination.

Cypher User Guide, Chapter: Tutorial: AC Mode in Air, Std. Scanner.

MFP-3D User Guide, Chapter: Installation.

MFP-3D User Guide, Chapter: Tutorial: AC Mode Imaging in Air.

Software User Guide, Programming and Automation., Should this be Basic Programming Tutorial?

If so, it's a placeholder. Otherwise I can't find it.

